



About DirectMed

DirectMed, Inc, was founded in 1997 as a source for a wide selection of high-quality components for medical devices. Since then, we have expanded our capabilities to offer full customization on any of our standard OEM components. Since 2018, we have been offering a full line of private-label devices for fluid management.

Our excellence in customer service, competitive pricing, and timely delivery are key to our continuing mission to provide the best possible quality on time, every time.

We maintain a wide stock of components and work with our manufacturing partners from around the world to source every possible option, OEM or fully custom-made. Our team will liaise with your product engineers from concept to prototypes to full production. We aim to minimize risk and add efficiency to the supply chain model from idea to market with our customized solutions and logistics support.

Whether you are seeking a source for OEM components or a complete, ready-to-use device, we can help you in the following key areas:

Supply Chain

We work with our contract manufacturing partners to optimize supply chain management. This ensures we can take the device from start to market in one comprehensive and efficient project path.

Fulfillment

Our manufacturing partners' facilities create and deliver millions of units every year. Each unit is tailored exactly to customer specifications including the product itself, packaging, and delivery schedule.

Proactive Troubleshooting

Our team's 25+ years in design and development allow us to identify hurdles associated with your project **in advance**. Predicting upcoming issues with production before they happen prevents delays in getting your product to market.

Partnership

We believe that being your supplier means being a part of your team. This also allows you to focus on what you do well and leave the manufacturing oversight to an experienced partner like us who has relationships with our key vendors developed over 25 years. Our manufacturing partners' facilities are ISO 13485 and ISO 9001 certified. Our products can be sterilized using a number of different methods.

About DirectMed

DirectMed Location

DirectMed, Inc. global headquarters is located at:

98 Pratt Oval, Glen Cove, NY 11542 USA

The facility is home to our extensive warehouse facility and the core of our logistics, sales, and customer service teams. If you need assistance in your native language, you may contact us in English, French, German, Spanish, or Italian.


DirectMed's suppliers consist of manufacturing facilities in Italy, Germany, Mexico, and Taiwan.

US 

+1 (516) 656-3377

sales@directmed.com

**98 Pratt Oval, Glen Cove, NY
11542 USA**



DirectMed

Table of Contents

Category	Page Number	Category	Page Number
Bag Ports	8	Kit Accessories	74
Caps		Luers	
<i>Non-Tethered</i>	11	<i>Male</i>	76
<i>Tethered</i>	18	<i>Barbed Male</i>	83
Connectors		<i>Female</i>	85
<i>Barbed</i>	19	<i>Barbed Female</i>	89
<i>Cardio</i>	23	<i>Specialty Luer Locks</i>	90
<i>Dialysis</i>	25	<i>Luer Collars / Luer Hubs</i>	92
<i>Funnel</i>	28	<i>Luer Slips</i>	93
<i>Large Bore</i>	29	Needle Hubs	94
<i>Suction</i>	30	NRFit® Neuraxial Connectors	95
<i>Tubing Connectors</i>	31	Stopcocks	97
<i>T-Connectors</i>	32	Stopcock Manifolds	100
<i>Y-Connectors</i>	33	Syringes	101
Drip Chambers		Tubing	102
<i>Non-Vented</i>	38	Tubing Clips	104
<i>Vented</i>	43	Tubing Clamps	
<i>Drip Chamber Covers</i>	45	<i>Pinch</i>	105
ENFit® Enteral Feeding System	47	<i>Roller</i>	110
Filters	50	<i>Slide</i>	113
Flow Regulators	58	Valves	
Injection Sites		<i>Check</i>	116
<i>Standard</i>	59	<i>Hemostasis</i>	121
<i>Negative Needleless</i>	63	<i>Drain</i>	122
<i>Neutral Needleless</i>	65	<i>Tuohy Borst</i>	123
IV Spikes		Index	125
<i>Non-Vented</i>	68	Abbreviation Index	128
<i>Vented</i>	71	Unit Conversion	130
<i>Vial Adapters</i>	73		

Upcoming Trade Shows

China
Medtec

Sep. 24-26
Shanghai, China

Booth # 2C308



Medical Technology
Ireland

Sep. 24-25
Galway, Ireland

Booth # 98G



Sep. 30-Oct. 1
Boston, MA, USA

Booth # 820

COMPAMED



MEDICA



Nov. 17-20
Düsseldorf, Germany

Booth # 8AE15

DirectMed

+1 (516) 656-3377 • sales@directmed.com
www.directmed.com

New Products

1



105

APHC0750



8

BA-023



19

BB-113



9

BC-078



25

BC-131



74

BH-011



74

BH-015



16

CA-161



16

CA-162



20

CB-041N



21

CB-041PP



21

CB-042N



21

CB-042PP



21

CB-043



21

CB-044N



21

CB-044PVDF



21

CB-045N



21

CB-046



21

CB-047



21

CB-048



21

CB-049



21

CB-050



21

CB-051



21

CB-052



21

CB-053B



21

CB-054



21

CB-055



21

CB-056



22

CB-057



22

CB-058



22

CB-059



22

CB-060



22

CB-061



22

CB-062



112

CR-024

DirectMed

+1 (516) 656-3377 • sales@directmed.com
www.directmed.com

						
CV-096	CV-097	CV-104	CV-105	CV-106	CV-107	CV-108
						
CV-278	CY-065	DC-161B	DC-193	DC-195	DC-196	DC-198
						
DY-057	DY-013ADHPF	DY-060DEHT	FC-025	FF-077B	FF-108	FF-127
						
FF-128	FF-129	FF-131	FF-141	FF-178	IS-141	IS-146
						
IS-147	IS-148	IS-149	IS-150B	IS-150G	IS-1500R	IS-150R



67



83



83



83



84



89



89

IS-150W

LB-036

LB-037

LB-038

LB-039

LB-040

LB-041



89



88



90



90



91



90



81

LB-042

LF-103

LF-123

LL-006

LL-018

LL-019

LM-139



76



81



77



77



94



94



94

LM-148B

LM-150

LM-151

LM-152

NH-014

NH-015

NH-016



94



94



94



94



94



94



72

NH-017

NH-018

NH-019

NH-020

NH-021

NH-022

SP-059



70



70



99



99



99



101

SP-060

SP-061


ST-083

ST-087

ST-088

SY-960

NEW



**MADE IN
ITALY**

**SIZES FROM
1 to 250 mL**

**HIGH
QUALITY**

MEDICAL GRADE

DirectMed

+1 (516) 656-3377 • sales@directmed.com
www.directmed.com

Features



SYG-010LS

Glass Syringes are made using high-quality raw materials. Production takes place in our facility, starting from neutral glass tubes that are resistant to thermal shock, impact, chemical corrosion, and repeated sterilization.

The use of glass also helps to prevent the risk of samples contamination. The metal nozzles are made from high-quality nickel-plated or chrome-plated brass, ensuring a safe and precise product in every detail. The barrel of each syringe is marked with a volumetric scale and can be customized with the customer's logo upon request.

LOR Syringes



SYG-005LL-LOR

SYG-010LL-LOR

LOR (Loss of Resistance) Syringes are designed to compose kits for epidural anesthesia. In conjunction with the appropriate needle, they allow for easy localization of the epidural space and enable smooth action of the plunger thanks to a low-resistance design that prevents any friction.

Our LOR glass syringes are FDA (Food and Drug Administration) approved, which certifies the high quality of the product. The luer connection complies with ISO 80369-7 standards, and the NRFit connection complies with ISO 80369-6 standards. The standard packaging consists of plastic trays containing 6 pieces, packaged in cardboard boxes of 36. Available in 5 or 10 mL size.

Numerous Options



1-250mL Sizes Available

Our glass syringes are manufactured for use in the medical and veterinary fields. Both the barrel and the piston are made entirely of glass and are completed with different types of nozzles, in glass or metal. Depending on their intended use, our glass syringes are available in various capacities, ranging from 1 to 250 mL.

Our medical syringes comply with the general safety requirements and provisions of Regulation (EU) 2017/745 and its subsequent amendments.

The nozzles of our glass syringes comply with ISO 80369-6 and 80369-7 (formerly 594-1/594-2).



Metal Luer Lock

Material: Brass, Glass, Nickel, PTFE
Color: Natural
Sterilization: Steam

1 mL	30 mL
SYG-001LL	SYG-030LL
2 mL	50 mL
SYG-002LL	SYG-050LL
3 mL	100 mL
SYG-003LL	SYG-100LL
5 mL	150 mL
SYG-005LL	SYG-150LL
10 mL	200 mL
SYG-010LL	SYG-200LL
20 mL	250 mL
SYG-020LL	SYG-250LL



Glass Luer Slip

Material: Glass
Color: Natural
Sterilization: Steam

1 mL	30 mL
SYG-001LS	SYG-030LS
2 mL	50 mL
SYG-002LS	SYG-050LS
3 mL	100 mL
SYG-003LS	SYG-100LS
5 mL	150 mL
SYG-005LS	SYG-150LS
10 mL	200 mL
SYG-010LS	SYG-200LS
20 mL	250 mL
SYG-020LS	SYG-250LS



Metal Luer Slip

Material: Brass, Glass, Nickel, PTFE
Color: Natural
Sterilization: Steam

1 mL	30 mL
SYG-001LSM	SYG-030LSM
2 mL	50 mL
SYG-002LSM	SYG-050LSM
3 mL	100 mL
SYG-003LSM	SYG-100LSM
5 mL	150 mL
SYG-005LSM	SYG-150LSM
10 mL	200 mL
SYG-010LSM	SYG-200LSM
20 mL	250 mL
SYG-020LSM	SYG-250LSM



Metal Luer Lock - LOR

Material: Brass, Glass, Nickel, PTFE
Color: Natural
Sterilization: Steam

5 mL	5 mL (NRFit)
SYG-005LL-LOR	SYG-030NR-LOR
10 mL	30 mL
SYG-010LL-LOR	SYG-030LL-LOR
50 mL	
SYG-050LL-LOR	

Glass Luer Slip - LOR

Material: Glass
Color: Natural
Sterilization: Steam

5 mL
SYG-005LS-LOR
10 mL
SYG-010LS-LOR
30 mL
SYG-030LS-LOR

Metal Luer Slip - LOR

Material: Brass, Glass, Nickel, PTFE
Color: Natural
Sterilization: Steam

5 mL
SYG-005LSM-LOR
10 mL
SYG-010LSM-LOR
30 mL
SYG-030LSM-LOR

BA-009

Bag with Valve and Clamp

Material: PP, PVC
Color: Natural
Sterilization: Eto
Volume: 2 L



BA-023

NEW

Bag

Material: Film, PC, PVC
Color: Clear, Natural
Sterilization: Eto
Volume: 500 mL



BC-001

Bag Port, No Membrane, with Break-Off Cap

Material: PC
Color: Clear
Sterilization: Eto, E-Beam, Gamma, Steam
ID: 6.0 mm



BC-001R

Bag Port, No Membrane, with Break-Off Cap

Material: PC
Color: Red
Sterilization: Eto, E-Beam, Gamma, Steam
ID: 6.0 mm



BC-002

Bag Port, No Membrane, with Break-Off Cap

Material: PC
Color: Clear
Sterilization: Eto, E-Beam, Gamma, Steam
ID: 4.8 mm | OD: 5.2 mm



BC-003

Bag Port, No Membrane, with Break-Off Cap

Material: PC
Color: Clear
Sterilization: Eto, E-Beam, Gamma, Steam
ID: 3.8 mm | OD: 9.7 mm



BC-004

Bag Port, No Membrane, with Cap

Material: PVC
Color: Red
Sterilization: Eto
OD: 41.00 mm



BC-005

Bag Port, No Membrane, with Cap

Material: PVC
Color: Red
Sterilization: Eto
OD: 49.00 mm



BC-006

Bag Port, No Membrane, with Break-Off Cap and Female Luer

Material: PC
Color: Clear
Sterilization: Eto, E-Beam, Gamma, Steam
ID: 3.0 mm



BC-007

Bag Port, No Membrane, with Break-Off Cap and Female Luer

Material: PC
Color: Clear
Sterilization: Eto, E-Beam, Gamma, Steam
ID: 3.0 mm



BC-008

Bag Port, No Membrane, with Break-Off Cap and Female Luer

Material: PC
Color: White
Sterilization: Eto, Gamma, Steam
ID: 6.0 mm



BC-011

Twist-Off Spike Port with Membrane

Material: PVC
Color: Blue
Sterilization: Eto
ID: 8.0 mm



BC-012A

Bag Port with Membrane and Tear-Off Cap

Material: PVC DEHP Free
Color: Clear
Sterilization: Eto, Gamma
ID: 5.4 mm | OD: 7.1 mm



BC-016

Piercing Connector with Membrane

Material: PVC
Color: Clear
Sterilization: Eto, Gamma
ID: 5.0 mm



BC-029

Bag Port with Membrane, Tear-Off Cap and Air Filter

Material: Acrylic, PP, PVC-DEHT
Color: Clear, White
Sterilization: Eto, Gamma
ID: 5.0 mm | OD: 6.8 mm
Pore Size: 3 µm
Pressure: 300 psi



BC-054DHPF

Bag Port with Membrane and Tear-Off Cap

Material: PVC DEHP Free
Color: Clear
Sterilization: Eto, Gamma
ID: 5.5 mm | OD: 6.8 mm



BC-055

Bag Port with Membrane, Tear-Off Cap and Air Filter, Photoprotective

Material: PVC
Color: Yellow
Sterilization: Eto
OD: 6.0 mm



BC-056DEHT

Bag Port with Membrane and Tear-Off Cap

Material: PVC-DEHT
Color: Clear
Sterilization: Eto
OD: 6.8 mm



BC-057

Bag Port with Flange

Material: PVC
Color: Clear
Sterilization: Eto



BC-059

Bag Port with Flange

Material: PVC
Color: Clear
Sterilization: Eto



BC-063

Bag Port with Flange

Material: PVC
Color: Black
Sterilization: N/A
ID: 0.281"



BC-064

Bag Port with Flange

Material: PUR, PVC
Color: Clear
Sterilization: Eto
ID: 9.4 mm



BC-067

Bag Port with Membrane and Tear-Off Cap

Material: EVA
Color: Natural
Sterilization: Eto, Gamma



BC-068

Bag Port with Membrane and Tear-Off Cap

Material: TPE
Color: Natural
Sterilization: Eto, E-Beam, Gamma, Steam



BC-068E

Bag Port with Membrane and Tear-Off Cap

Material: EVA
Color: Natural
Sterilization: Eto, Gamma



BC-078

NEW

Bag Port, No Membrane - Bag Valve

Material: PVC
Color: Blue
Sterilization: Eto



BC-091

Bag Port with Male Luer

Material: MABS, PVC
DEHP Free
Color: Clear
Sterilization: N/A



BC-111

Bag Port with Membrane and Tear-Off Cap

Material: PVC-Free, SEBS
Color: Natural
Sterilization: Eto, Gamma, Steam
ID: 5.2 mm | OD: 7.1 mm



BC-111TPE

Bag Port with Membrane and Tear-Off Cap

Material: TPE
Color: Natural
Sterilization: Eto, Gamma, Steam
ID: 5.2 mm | OD: 7.1 mm



BC-112

Bag Port with Membrane and Tear-Off Cap

Material: EVA
Color: Clear
Sterilization: Eto, Gamma
ID: 5.5 mm



BC-113

Bag Port with Membrane and Tear-Off Cap

Material: PVC DEHP Free
Color: Clear
Sterilization: Eto, Steam
ID: 6.8 mm



BC-119

Bag Port with Membrane and Tear-Off Cap

Material: PVC
Color: Clear
Sterilization: Eto



CA-001

Male to Female Luer Cap,
Non-Vented

Material: PE
Color: White
Sterilization: Eto



CA-001A

Male to Female Luer Cap,
Non-Vented

Material: HDPE
Color: White
Sterilization: Eto



CA-003

Female Luer Cap,
Non-Vented

Material: HDPE
Color: Natural
Sterilization: Eto, E-Beam,
Gamma



CA-005

Female Luer Cap,
No Stem

Material: HDPE
Color: Light Blue
Sterilization: Eto, Gamma



CA-008

Spike Cap, Vented

Material: LDPE
Color: White
Sterilization: Eto, Gamma



CA-012

Spike Cap

Material: LDPE
Color: Natural
Sterilization: Eto, Gamma



CA-014

Female Luer Cap,
Non-Vented

Material: HDPE
Color: Clear
Sterilization: Eto, Gamma



CA-017

Male to Female Luer Cap,
Non-Vented

Material: PC
Color: Clear
Sterilization: Eto, Gamma



CA-022

Female Luer Cap, Vented

Material: ABS
Color: Green
Sterilization: Eto, Gamma



CA-023

Male Luer Cap,
Non-Vented

Material: LDPE
Color: Natural
Sterilization: Eto, Gamma



CA-024

Female Luer Cap,
Non-Vented

Material: ABS
Color: White
Sterilization: Eto, Gamma



CA-024PP

Female Luer Cap,
Non-Vented

Material: PP
Color: White
Sterilization: Eto, Steam



CA-025

Female Luer Cap, Vented

Material: PP
Color: Natural
Sterilization: Eto, Gamma,
Steam



CA-025BLUE

Female Luer Cap, Vented

Material: PP
Color: Blue
Sterilization: Eto, Steam



CA-025BLUEA

Female Luer Cap, Vented

Material: PP
Color: Blue
Sterilization: Eto, Steam



CA-025G

Female Luer Cap, Vented

Material: HDPE
Color: Clear
Sterilization: Eto



CA-025NG

Female Luer Cap, Vented

Material: PP
Color: Clear
Sterilization: Eto, Steam



CA-025RED

Female Luer Cap, Vented

Material: PP
Color: Red
Sterilization: Eto, Steam



CA-027

Female Luer Cap, Vented

Material: PP

Color: Blue

Sterilization: Eto, Steam



CA-030

Male to Female Luer Cap,
Non-Vented

Material: HDPE

Color: Light Blue

Sterilization: Eto



CA-030RED

Male to Female Luer Cap,
Non-Vented

Material: HDPE

Color: Red

Sterilization: Eto



CA-031

Male to Female Luer Cap,
Vented

Material: HDPE

Color: White

Sterilization: Eto



CA-032

Male Luer Cap, Vented

Material: PVC DEHP Free

Color: Clear

Sterilization: Eto



CA-032B

Male Luer Cap, Vented

Material: PVC DEHP Free

Color: Blue

Sterilization: Eto



CA-032G

Male Luer Cap, Vented

Material: PVC DEHP Free

Color: Clear

Sterilization: Eto, E-Beam,
Gamma



CA-032PP

Male Luer Cap, Vented

Material: PP

Color: Clear

Sterilization: Eto, Steam



CA-032R

Male Luer Cap, Vented

Material: PVC

Color: Red

Sterilization: Eto



CA-032RG

Male Luer Cap, Vented

Material: PVC DEHP Free

Color: Red

Sterilization: Eto, E-Beam,
Gamma



CA-033

Male Luer Cap,
Non-Vented

Material: ABS

Color: White

Sterilization: Eto, Gamma



CA-039

Male Luer Cap, Vented

Material: LDPE

Color: Clear

Sterilization: Eto, Gamma



CA-041

Male Luer Slip Cap

Material: LDPE

Color: Natural

Sterilization: Eto, Gamma



CA-042

Male Luer Slip Cap

Material: PE

Color: White

Sterilization: Eto



CA-043

Male Luer Cap, Vented

Material: LDPE

Color: Clear

Sterilization: Eto, Gamma



CA-048

Female Luer Cap, Vented

Material: HDPE

Color: White

Sterilization: Eto



CA-049W

Spike Cap, Vented

Material: LDPE

Color: White

Sterilization: Eto



CA-050

Spike Cap, Vented

Material: LDPE

Color: Clear

Sterilization: Eto, Gamma



CA-052

Female Luer Cap,
Non-Vented

Material: ABS
Color: Blue
Sterilization: Eto, Gamma



CA-052A

Female Luer Cap,
Non-Vented

Material: ABS
Color: Blue
Sterilization: Eto, Gamma



CA-053

Spike Cap, Vented

Material: LDPE
Color: White
Sterilization: Eto, Gamma



CA-054

Spike Cap, Vented

Material: LDPE
Color: Blue
Sterilization: Eto, Gamma



CA-055

Female Luer Cap,
Non-Vented

Material: PP
Color: Red
Sterilization: Eto, Steam



CA-055ABS

Female Luer Cap,
Non-Vented

Material: ABS
Color: Red
Sterilization: Eto, Gamma



CA-055B

Female Luer Cap,
Non-Vented

Material: PP
Color: Blue
Sterilization: Eto, Steam



CA-060

Female Luer Cap,
Non-Vented

Material: HDPE
Color: Yellow
Sterilization: Eto



CA-061

Female Luer Cap, Vented

Material: HDPE
Color: Yellow
Sterilization: Eto



CA-062

Female Luer Cap,
Non-Vented

Material: PP
Color: White
Sterilization: Eto, Steam



CA-066

Cardio Connector Cap,
Vented

Material: LDPE
Color: Blue
Sterilization: Eto, Gamma



CA-068

Female Luer Cap, Vented,
No Stem

Material: LDPE
Color: Natural
Sterilization: Eto, Gamma



CA-068NV

Female Luer Cap, Non-
Vented, No Stem

Material: LDPE
Color: Natural
Sterilization: Eto, Gamma



CA-070

Male Luer Cap with
Air Vent Filter

Material: Acrylic, HDPE, PP
Color: Clear, White
Sterilization: Eto
Pore Size: 3 µm



CA-080

Male / Female Luer Cap,
Non-Vented

Material: PP
Color: Blue
Sterilization: Eto, Gamma



CA-080W

Male / Female Luer Cap,
Non-Vented

Material: PP
Color: White
Sterilization: Eto, Gamma



CA-086A

Male Luer Cap, Vented
with Air Filter

Material: Acrylic, PP
Color: Blue, White
Sterilization: Eto, Gamma
Microns: 1.2 µm



CA-087

Flexible Cap, Vented

Material: TPE
Color: Red
Sterilization: Eto
ID: 1/4"



CA-087W

Flexible Cap, Vented

Material: TPE
Color: White
Sterilization: Eto
ID: 1/4"



CA-089

Flexible Cap

Material: TPE
Color: Yellow
Sterilization: Eto
ID: 1/4"



CA-091

Female Luer Cap,
Non-Vented

Material: PC
Color: Clear
Sterilization: Eto, Steam



CA-093

Female Luer Cap,
Non-Vented

Material: HDPE
Color: Blue
Sterilization: Eto, Gamma



CA-094

Male Luer Cap,
Non-Vented

Material: HDPE
Color: White
Sterilization: Eto



CA-094GS

Male Luer Cap,
Non-Vented

Material: HDPE
Color: White
Sterilization: Eto, Gamma



CA-094V

Male Luer Cap, Vented

Material: HDPE
Color: White
Sterilization: Eto, Gamma



CA-095

Male Luer Cap,
Non-Vented

Material: PP
Color: Natural
Sterilization: Eto



CA-099

Male Luer Cap, Vented

Material: LDPE
Color: Natural
Sterilization: Eto, Gamma



CA-100

Spike Cap, Vented

Material: LDPE
Color: Clear
Sterilization: Eto



CA-103

Female Luer Cap,
No Stem

Material: PP
Color: Natural
Sterilization: Eto



CA-104

Female Luer Cap,
No Stem

Material: PP
Color: Natural
Sterilization: Eto



CA-106

Female Luer Cap, Vented

Material: PP
Color: Clear
Sterilization: Eto



CA-107

Female Luer Cap, Non-
Vented, Long Stem

Material: ABS
Color: Blue
Sterilization: Eto



CA-108

Female Luer Cap,
No Stem

Material: LDPE
Color: Red
Sterilization: Eto



CA-109

Spike Cap

Material: LDPE
Color: White
Sterilization: Eto



CA-109B

Male Luer Cap, Vented

Material: HDPE
Color: Blue
Sterilization: Eto



CA-109N

Male Luer Cap, Vented

Material: HDPE
Color: Natural
Sterilization: Eto



CA-1090

Male Luer Cap, Vented

Material: HDPE
Color: Orange
Sterilization: Eto



CA-112

Female Luer Cap,
Non-Vented

Material: ABS
Color: White
Sterilization: Eto



CA-113

Spike Cap, Vented

Material: LDPE
Color: Clear
Sterilization: Eto



CA-114

Male Luer Cap,
Non-Vented

Material: PP
Color: Clear
Sterilization: Eto, E-Beam,
Gamma, Steam



CA-115

Male Luer Cap, Vented

Material: PP
Color: Clear
Sterilization: Eto, Gamma,
Steam



CA-116

Female Luer Cap,
Non-Vented

Material: HDPE
Color: Blue
Sterilization: Eto, Gamma



CA-117

Female Luer Cap,
Non-Vented

Material: HDPE
Color: Natural
Sterilization: Eto, Gamma



CA-118

Female Luer Cap,
Non-Vented

Material: HDPE
Color: White
Sterilization: Eto, Gamma



CA-119

Female Luer Cap,
Non-Vented

Material: PC
Color: Clear
Sterilization: Eto, Gamma



CA-120

Female Luer Cap,
Non-Vented

Material: HDPE
Color: Green
Sterilization: Eto, Gamma



CA-120V

Female Luer Cap, Vented

Material: HDPE
Color: Green
Sterilization: Eto, Gamma



CA-121

Female Luer Cap,
Non-Vented

Material: HDPE
Color: Red
Sterilization: Eto, Gamma



CA-122

Male Luer Cap, Vented

Material: LDPE
Color: White
Sterilization: Eto, Gamma
ID: 4.5 mm



CA-123

Male Luer Cap, Vented

Material: HDPE
Color: Blue
Sterilization: Eto



CA-124

Male Luer Cap, Vented

Material: PP
Color: Natural
Sterilization: Eto
ID: 4.5 mm



CA-125

Male Luer Cap, Vented

Material: LDPE
Color: Red
Sterilization: Eto, Gamma



CA-129

Male Luer Cap,
Non-Vented

Material: LDPE
Color: Orange
Sterilization: Eto, Gamma



CA-131

Female Luer Cap,
Non-Vented

Material: ABS
Color: White
Sterilization: Eto



CA-131V

Female Luer Cap, Vented

Material: ABS
Color: White
Sterilization: Eto



CA-132

Female Luer Cap

Material: ABS
Color: White
Sterilization: Eto, E-Beam, Gamma



CA-133

Female Luer Cap with Air Filter

Material: Acrylic, PP
Color: Clear
Sterilization: Eto
Pore Size: 3 µm



CA-133B

Female Luer Cap with Air Filter

Material: Acrylic, PP
Color: Blue
Sterilization: Eto
Pore Size: 3 µm



CA-134

Male Luer Cap with Air Filter

Material: Acrylic, PP
Color: Natural
Sterilization: Eto
Pore Size: 3 µm



CA-135

Male Luer Cap with Air Filter

Material: Acrylic, PE, PP
Color: Red
Sterilization: Eto
Pore Size: 1.2 µm



CA-136

Female Luer Cap, Vented

Material: HDPE
Color: Natural
Sterilization: Eto, Gamma



CA-148

Spike Cap, Vented

Material: PP
Color: Clear
Sterilization: Eto, Steam
ID: 7.4 mm



CA-151

Female Luer Cap, Vented, with Air Filter

Material: Acrylic, PE, PP
Color: Clear, Red
Sterilization: Eto
Pore Size: 1.2 µm



CA-151B

Female Luer Cap, Vented, with Air Filter

Material: Acrylic, PE, PP
Color: Blue, Red
Sterilization: Eto
Pore Size: 1.2 µm



CA-152

Male Luer Cap, Vented, with Air Filter

Material: Acrylic, LDPE, PP
Color: Natural, Red
Sterilization: Eto
Pore Size: 1.2 µm



CA-153

Male Luer Cap, Vented, with Air Filter

Material: Acrylic, PP
Color: Natural, Red
Sterilization: Eto, Gamma
Pore Size: 1.2 µm



CA-158

Male Luer Cap, Vented

Material: HDPE
Color: White
Sterilization: Eto



CA-160

Female Luer Cap, Non-Vented, Long Stem

Material: PP
Color: White
Sterilization: Eto



CA-161

NEW

Male to Female Luer Cap, Vented

Material: HDPE
Color: Natural
Sterilization: Eto



CA-162

NEW

Syringe Cap, Threaded

Material: HDPE
Color: Red
Sterilization: Eto



CA-164

Male Luer Cap, Vented

Material: HDPE
Color: Natural
Sterilization: Eto



CA-165

Female Luer Cap, Dead-End

Material: HDPE
Color: Natural
Sterilization: Eto, Steam



CA-167

Female Luer Cap, Vented,
No Stem

Material: PP
Color: Clear
Sterilization: Eto, Gamma



CA-168

Female Luer Cap, Tear-Off

Material: LDPE
Color: Natural
Sterilization: Eto, Gamma



CA-171

Male Luer Cap, Vented,
with Air Filter

Material: Acrylic, LDPE, PP
Color: Natural, White
Sterilization: Eto, Steam
Pore Size: 3 µm



CA-177

Spike Cap for SP-177,
Vented

Material: LDPE
Color: Natural
Sterilization: Eto



CA-200

Male Luer Cap, Vented

Material: PP
Color: Natural
Sterilization: Eto



CA-324PP

Female Luer Cap,
Non-Vented

Material: PP
Color: White
Sterilization: Eto, Gamma



CA-009

Tethered Cap, Non-Vented, for Female Luers

Material: PP
Color: Clear
Sterilization: Eto, E-Beam, Gamma



CA-016DHPF

Tethered Cap, Non-Vented, for Female Luers

Material: HDPE, Latex-Free, PVC DEHP Free
Color: White
Sterilization: Eto



CA-016ST

Cap Tether Only for CA-016DHPF

Material: PVC DEHP Free
Color: White
Sterilization: Eto



CA-021

Tethered Cap, Vented, for Female Luers, Non-Threaded

Material: LDPE
Color: Natural
Sterilization: Eto, Gamma



CA-021MP

Tethered Cap, Non-Vented, for Male Luers, Non-Threaded

Material: LDPE
Color: Natural
Sterilization: Eto, Gamma



CA-045

Tethered Cap, Vented, for Injection Sites

Material: LDPE
Color: Clear
Sterilization: Eto, Gamma



CA-046

Tethered Cap, Vented, with Air Filter for Injection Sites

Material: LDPE
Color: Clear
Sterilization: Eto, Gamma



CA-063

Tethered Cap, Vented, with Filter for Female Luers

Material: Acrylic, PP
Color: Clear, White
Sterilization: Eto
Microns: 3 µm



CA-065S

Tethered Cap, Non-Vented, for Female Luers

Material: HDPE
Color: Blue
Sterilization: Eto, Gamma



CA-065SW

Tethered Cap, Non-Vented, for Female Luers

Material: HDPE
Color: White
Sterilization: Eto, Steam



CA-127

Tethered Cap for Female Luers

Material: HDPE, PVC
Color: White
Sterilization: Eto



CA-128

Tethered Cap for Female Luers

Material: HDPE, PVC
Color: Blue
Sterilization: Eto



CA-169

Cap Tether Only for CA-133

Material: PVC-DEHT
Color: Clear
Sterilization: Eto, Gamma
ID: 6.3 mm



CA-170

Cap Tether Only for CA-133

Material: PVC-DEHT
Color: Clear
Sterilization: Eto, Gamma
ID: 6.4 mm



BB-113

NEW

Barbed T-Connector,
3/32" OD Barbs

Material: ABS
Color: White
Sterilization: Eto, Gamma



BC-040DHPF

Barbed Connector

Material: PVC DEHP Free
Color: Clear
Sterilization: Eto, E-Beam,
Gamma
ID: 6.4 mm | OD: 8.9 mm



BC-066

Stepped Connector

Material: PVC
Color: Natural
Sterilization: Eto
ID: 6.6 mm



BC-071

Stepped Connector

Material: ABS
Color: White
Sterilization: Eto, Gamma
ID: 4.1 mm



CB-003

Barbed Connector

Material: PP
Color: Natural
Sterilization: Eto
OD: 7.0 - 10.00 mm



CB-012

Barbed T-Connector

Material: PP
Color: Clear
Sterilization: Eto, Gamma
OD: 7.0 - 9.0 mm



BC-033

Enteral Feeding Connector

Material: ABS
Color: White
Sterilization: Eto, Gamma
ID: 4.1 mm | OD: 5.2 mm



BC-043DEHT

Stepped Connector

Material: PVC-DEHT
Color: Clear
Sterilization: Eto
ID: 4.1 mm | OD: 4.8 mm



BC-069

Stepped Connector

Material: ABS
Color: White
Sterilization: Eto, Gamma
OD: 6.0 mm



CB-001A

Barbed Connector

Material: PP
Color: Natural
Sterilization: Eto, Steam
OD: 7.5 - 11.00 mm



CB-006

Barbed Connector

Material: PP
Color: Natural
Sterilization: Eto
OD: 10.6 - 14.00 mm



CB-015

Barbed Y-Connector

Material: PP
Color: Natural
Sterilization: Eto
OD: 6.7 - 7.9 mm



BC-034

Stepped Connector

Material: ABS
Color: White
Sterilization: Eto, Gamma
ID: 5.5 mm | OD: 9.0 mm



BC-065

Stepped Connector

Material: PVC
Color: Natural
Sterilization: Eto
ID: 5.5 mm



BC-070

Vented Cap for BC-069

Material: LDPE
Color: Natural
Sterilization: Eto, Steam
ID: 10.00 mm



CB-002

Barbed Connector

Material: PP
Color: Natural
Sterilization: Eto, Steam
OD: 4.1 - 6.1 x
7.5 - 11.6 mm



CB-008

Barbed Connector

Material: PP
Color: Clear
Sterilization: Eto, Steam
OD: 4.0 - 6.0 mm



CB-016

Barbed Y-Connector

Material: PP
Color: Natural
Sterilization: Eto, Steam
OD: 13.00 mm



CB-022

Stepped Connector

Material: ABS
Color: White
Sterilization: Eto, Gamma
OD: 7.0 - 10.00 mm



CB-023

Barbed Connector

Material: PP
Color: Natural
Sterilization: Eto
OD: 7.5 - 10.00 mm



CB-024

Barbed Y-Connector,
3/16" x 1/4" x 5/16" OD
Barbs

Material: PP
Color: Natural
Sterilization: Eto, Gamma



CB-025

Barbed Connector, 10-32
Thread to 3/16" Barb

Material: Nylon
Color: White
Sterilization: Eto, Gamma



CB-027

Barbed T-Connector,
3/32" x 3/32" x 3/32" OD
Barbs

Material: Nylon
Color: White
Sterilization: Eto, Gamma



CB-030

Barbed Connector

Material: PS
Color: Natural
Sterilization: Eto
OD: 7.2 mm



CB-031

Barbed Connector

Material: PS
Color: Natural
Sterilization: Eto
OD: 4.0 mm



CB-032

Barbed Connector

Material: PS
Color: Natural
Sterilization: Eto
OD: 10.00 mm



CB-033

Barbed Connector

Material: PS
Color: Natural
Sterilization: Eto
OD: 4.0 - 7.2 mm



CB-034

Barbed Connector

Material: PS
Color: Natural
Sterilization: Eto
OD: 4.0 - 10.00 mm



CB-035

Barbed Connector

Material: PS
Color: Natural
Sterilization: Eto
OD: 7.0 - 10.00 mm



CB-036

Barbed Connector

Material: PP
Color: Natural
Sterilization: Eto
ID: 4.4 - 4.9 mm
OD: 7.3 mm



CB-037

Barbed Connector,
1/4" x 1/4" OD Barbs

Material: PP
Color: Clear
Sterilization: Eto, E-Beam,
Gamma



CB-038

Barbed Connector,
3/16" x 1/8" OD Barbs

Material: PP
Color: Clear
Sterilization: Eto, E-Beam,
Gamma



CB-039

Barbed Connector,
3/32" x 3/32" OD Barbs

Material: PP
Color: Natural
Sterilization: Eto, E-Beam,
Gamma, Steam



CB-040

Barbed Connector,
3/32" to 1/16" OD Barbs

Material: PP
Color: Natural
Sterilization: Eto, E-Beam,
Gamma, Steam



CB-041N

NEW

Barbed T-Connector,
1/8" OD Barbs

Material: Nylon
Color: White
Sterilization: Eto



CB-041PP

NEW

Barbed T-Connector,
1/8" OD Barbs

Material: PP
Color: Clear
Sterilization: Eto, E-Beam,
Gamma



CB-042N

NEW

Barbed Connector,
1/4" to 1/8" OD Barbs

Material: Nylon
Color: Natural
Sterilization: Eto, Gamma,
Steam



CB-042PP

NEW

Barbed Connector,
1/4" to 1/8" OD Barbs

Material: PP
Color: Clear
Sterilization: Eto, E-Beam,
Gamma



CB-043

NEW

Barbed Y-Connector,
1/8" OD Barbs

Material: Nylon
Color: White
Sterilization: Eto



CB-044N

NEW

Barbed Elbow Connector,
3/16" OD Barbs

Material: Nylon
Color: White
Sterilization: Eto



CB-044PVDF

NEW

Barbed Elbow Connector,
3/16" OD Barbs

Material: PVDF
Color: Natural
Sterilization: Eto, Gamma,
Steam



CB-045N

NEW

Barbed T-Connector,
1/8" OD Barbs

Material: Nylon
Color: White
Sterilization: Eto, Gamma



CB-046

NEW

Barbed Y-Connector,
1.6 mm (1/16") OD Barbs

Material: Nylon
Color: White
Sterilization: Eto



CB-047

NEW

Barbed Y-Connector,
2.4 mm (3/32") OD Barbs

Material: Nylon
Color: White
Sterilization: Eto



CB-048

NEW

Barbed Y-Connector,
3.2 mm (1/8") OD Barbs

Material: Nylon
Color: White
Sterilization: Eto



CB-049

NEW

Barbed Elbow Connector,
3.2 mm (1/8") OD Barbs

Material: Nylon
Color: White
Sterilization: Eto



CB-050

NEW

Barbed Elbow Connector,
4.0 mm (5/32") OD Barbs

Material: Nylon
Color: White
Sterilization: Eto



CB-051

NEW

Barbed Elbow Connector,
4.8 mm (3/16") OD Barbs

Material: Nylon
Color: White
Sterilization: Eto



CB-052

NEW

Barbed T-Connector,
1.6 mm (1/16") OD Barbs

Material: Nylon
Color: White
Sterilization: Eto



CB-053B

NEW

Barbed T-Connector,
2.4 mm (3/32") OD Barbs

Material: Nylon
Color: Blue
Sterilization: Eto



CB-054

NEW

Barbed T-Connector,
3.2 mm (3/32") OD Barbs

Material: Nylon
Color: White
Sterilization: Eto



CB-055

NEW

Barbed T-Connector,
4.0 mm (5/32") OD Barbs

Material: Nylon
Color: White
Sterilization: Eto



CB-056

NEW

Barbed T-Connector,
4.8 mm (3/16") OD Barbs

Material: Nylon
Color: White
Sterilization: Eto



CB-057

NEW

Barbed Connector,
4.0 mm (5/32") OD Barbs

Material: Nylon
Color: White
Sterilization: Eto



CB-058

NEW

Barbed Connector,
4.8 mm (3/16") OD Barbs

Material: Nylon
Color: White
Sterilization: Eto



CB-061

NEW

Barbed Connector,
2.4 mm (3/32") OD Barbs

Material: Nylon
Color: White
Sterilization: Eto



LB-022

Barbed Connector,
3/32" x 1/16" OD Barbs

Material: Nylon
Color: White
Sterilization: Eto, Gamma



TT-064

Barbed Connector,
5/32" to 3/32" OD Barbs

Material: Nylon
Color: White
Sterilization: Eto



CB-059

NEW

Barbed Connector,
6.4 mm (1/4") OD Barbs

Material: Nylon
Color: White
Sterilization: Eto



CB-062

NEW

Barbed Connector,
3.2 mm (1/8") OD Barbs

Material: Nylon
Color: White
Sterilization: Eto



LB-030

Elbow Connector,
1/8" x 1/8" OD Barbs

Material: Nylon
Color: White
Sterilization: Eto



CB-060

NEW

Barbed Connector,
1.6 mm (1/16") OD Barbs

Material: Nylon
Color: White
Sterilization: Eto



LB-008

Barbed Connector,
1/8" x 1/16" OD Barbs

Material: PP
Color: Natural
Sterilization: Eto, E-Beam,
Gamma



TC-043

Stepped Connector

Material: PVC
Color: Green
Sterilization: Eto
OD: 6.8 - 11.3 mm



CC-010

Cardio Connector, Straight

Material: PC
Color: Clear
Sterilization: Eto
OD: 1/4" x 1/4"



CC-011

Cardio Connector, Straight

Material: PC
Color: Clear
Sterilization: Eto
OD: 3/16" x 3/16"



CC-012

Cardio Connector, Straight

Material: PC
Color: Clear
Sterilization: Eto
OD: 1/4" x 3/16"



CC-013

Cardio Y-Connector

Material: PP
Color: Clear
Sterilization: Eto, Gamma,
Steam
OD: 1/8" x 1/8" x 1/8"



CC-014

Cardio Connector, Straight

Material: PC
Color: Clear
Sterilization: Eto
OD: 3/8" x 3/8"



CC-016

Cardio Connector, Straight

Material: PC
Color: Clear
Sterilization: Eto
OD: 3/8" x 1/2"



CC-017

Cardio Connector, Straight

Material: PC
Color: Clear
Sterilization: Eto
OD: 1/4" x 3/8"



CC-018

Cardio Connector, Straight

Material: PC
Color: Clear
Sterilization: Eto
OD: 1/2" x 1/2"



CC-019

Cardio Connector, Straight
with Sampling Port

Material: PC
Color: Clear
Sterilization: Eto
OD: 1/4" x 1/4"



CC-020

Cardio Connector, Straight
with Luer Port

Material: PC
Color: Clear
Sterilization: Eto
OD: 3/16" x 3/16"



CC-021

Cardio Connector, Straight
with Luer Port

Material: PC
Color: Clear
Sterilization: Eto
OD: 1/4" x 3/16"



CC-022

Cardio Connector, Straight
with Luer Port

Material: PC
Color: Clear
Sterilization: Eto
OD: 1/4" x 3/8"



CC-023

Cardio Connector, Straight
with Luer Port

Material: PC
Color: Clear
Sterilization: Eto
OD: 3/8" x 3/8"



CC-025

Cardio Connector, Straight
with Luer Port

Material: PC
Color: Clear
Sterilization: Eto
OD: 1/2" x 3/8"



CC-027

Cardio Connector, Straight
with Luer Port

Material: PC
Color: Clear
Sterilization: Eto
OD: 1/2" x 1/2"



CC-028

Cardio Y-Connector

Material: PC
Color: Clear
Sterilization: Eto
OD: 1/4" x 3/8" x 3/8"



CC-029

Cardio Y-Connector

Material: PC
Color: Clear
Sterilization: Eto
OD: 3/16" x 3/16" x 3/16"



CC-030

Cardio Y-Connector

Material: PC
Color: Clear
Sterilization: Eto
OD: 1/4" x 1/4" x 1/4"



CC-031

Cardio Y-Connector

Material: PC
Color: Clear
Sterilization: Eto
OD: 3/8" x 1/4" x 3/8"



CC-032

Cardio Y-Connector

Material: PC
Color: Clear
Sterilization: Eto
OD: 3/8" x 3/8" x 3/8"



CC-033

Cardio Y-Connector

Material: PC
Color: Clear
Sterilization: Eto
OD: 3/8" x 1/4" x 1/4"



CC-034

Cardio Y-Connector

Material: PC
Color: Clear
Sterilization: Eto
OD: 1/2" x 3/8" x 3/8"



CC-035

Cardio Y-Connector

Material: PC
Color: Clear
Sterilization: Eto
OD: 1/2" x 1/2" x 1/2"



CC-036

Cardio Y-Connector

Material: PC
Color: Clear
Sterilization: Eto
OD: 1/2" x 3/8" x 1/2"



CC-038

Cardio Y-Connector with
Luer Port

Material: PC
Color: Clear
Sterilization: Eto
OD: 3/16" x 3/16" x 3/16"



CC-039

Cardio Y-Connector with
Luer Port

Material: PC
Color: Clear
Sterilization: Eto
OD: 1/4" x 1/4" x 1/4"



CC-040

Cardio Y-Connector with
Luer Port

Material: PC
Color: Clear
Sterilization: Eto
OD: 3/8" x 1/4" x 3/8"



CC-041

Cardio Y-Connector with
Luer Port

Material: PC
Color: Clear
Sterilization: Eto
OD: 3/8" x 3/8" x 3/8"



CC-042

Cardio Y-Connector with
Luer Port

Material: PC
Color: Clear
Sterilization: Eto
OD: 3/8" x 1/4" x 1/4"



CC-043

Cardio Y-Connector with
Luer Port

Material: PC
Color: Clear
Sterilization: Eto
OD: 1/2" x 3/8" x 1/2"



CC-044

Cardio Y-Connector with
Luer Port


Material: PC
Color: Clear
Sterilization: Eto
OD: 1/2" x 1/2" x 1/2"



CC-045

Cardio Y-Connector with
Luer Port

Material: PC
Color: Clear
Sterilization: Eto
OD: 1/2" x 3/8" x 1/2"



CC-047

Cardio Connector,
Straight, 1/8" x 1/8" Barbs


Material: PP
Color: Natural
Sterilization: Eto, Gamma
OD: 1/8" x 1/8"



CC-048

Cardio Y-Connector with
Luer Port


Material: PC
Color: Clear
Sterilization: Eto
OD: 3/16" x 3/16" x 1/4"



CC-049

Cardio Y-Connector

Material: PP
Color: Clear
Sterilization: Eto, Gamma,
Steam
OD: 1/8" x 1/8" x 1/8"



BC-048DEHT

Pump Connector

Material: PVC-DEHT
Color: Clear
Sterilization: Eto
ID: 5.5 to 6.8 mm



BC-052

Barbed Connector

Material: ABS
Color: Clear
Sterilization: Eto, Gamma
ID: 4.8 mm | OD: 6.8 mm



BC-053

Barbed Connector

Material: PVC
Color: Clear
Sterilization: Eto
ID: 4.8 mm | OD: 5.5 mm



BC-131

NEW

Hansen Connector
with Cap

Material: PVC-DEHT
Color: Clear
Sterilization: Eto
ID: 6.6 mm



CH-001DHPF

Pump Connector

Material: PVC DEHP Free
Color: Clear
Sterilization: Eto
ID: 4.1 mm



CH-003

Pump Connector

Material: ABS
Color: Clear
Sterilization: Eto, Gamma
ID: 4.1 mm



CH-005

Pump Connector

Material: PVC
Color: Clear
Sterilization: Eto, E-Beam,
Gamma
ID: 6.0 mm
OD: 4.1, 8.0 mm



DY-001

Dialysis Connector, Male

Material: PVC
Color: Clear
Sterilization: E-Beam,
Gamma
ID: 6.8 mm



DY-002

Dialysis Connector, Female

Material: PVC
Color: Clear
Sterilization: Eto
ID: 6.8 mm



DY-002R

Dialysis Connector, Female

Material: PVC
Color: Red
Sterilization: Eto
ID: 6.8 mm



DY-004ADEHT

Dialysis T-Connector

Material: PVC-DEHT
Color: Clear
Sterilization: Eto
ID: 5.5 mm | OD: 6.8 mm



DY-006

Reducing Adapter

Material: PVC
Color: Clear
Sterilization: Eto
ID: 6.8 to 5.5 mm



DY-007BDEHT

Dialysis Connector, Female

Material: PVC-DEHT
Color: Blue
Sterilization: Eto
ID: 6.8 mm



DY-007BG

Dialysis Connector, Female

Material: PVC-DEHT
Color: Blue
Sterilization: Eto, E-Beam,
Gamma
ID: 6.8 mm



DY-007DEHT

Dialysis Connector, Female

Material: PVC-DEHT
Color: Clear
Sterilization: Eto
ID: 6.8 mm



DY-007G

Dialysis Connector, Female

Material: PVC DEHP Free
Color: Clear
Sterilization: Eto, E-Beam,
Gamma

ID: 6.8 mm



DY-007RDEHT

Dialysis Connector, Female

Material: PVC-DEHT
Color: Red
Sterilization: Eto
ID: 6.8 mm



DY-007RG

Dialysis Connector, Female

Material: PVC-DEHT
Color: Red
Sterilization: Eto, E-Beam,
Gamma

ID: 6.8 mm



DY-008DEHT

Dialysis T-Connector

Material: PVC-DEHT
Color: Clear
Sterilization: Eto, E-Beam,
Gamma

ID: 2.5, 6.8 mm



DY-018

Dialysis Connector, Female

Material: PVC DEHP Free
Color: Red
Sterilization: Eto, E-Beam,
Gamma

ID: 6.8 mm



DY-019

Dialysis Connector, Female

Material: PVC DEHP Free
Color: Blue
Sterilization: Eto, E-Beam,
Gamma

ID: 6.8 mm



DY-021

Dialysis Cap, Vented

Material: HDPE
Color: Natural
Sterilization: Eto, E-Beam,
Gamma



DY-022

Dialysis Cap, Non-Vented

Material: PP
Color: Natural
Sterilization: Eto, Gamma,
Steam



DY-025A

Dialysis Cap, Vented,
Threaded, for Male
Connector

Material: HDPE
Color: White
Sterilization: Eto



DY-027

Dialysis Cap, Vented

Material: LDPE
Color: Natural
Sterilization: Eto, Gamma



DY-029

Recirculating Connector,
Non-Vented, with
Closed Cap

Material: HDPE
Color: White
Sterilization: Eto, E-Beam,
Gamma



DY-030

Pump Connector

Material: PVC DEHP Free
Color: Natural
Sterilization: Eto, E-Beam,
Gamma

OD: 6.8 mm



DY-031

Arterial Pillow Chamber
with Valve

Material: PVC DEHP Free
Color: Clear
Sterilization: Eto, E-Beam
ID: 5.5 mm
Volume: 9 mL



DY-033DEHT

Arterial Pillow Chamber

Material: PVC
Color: Clear
Sterilization: Eto, E-Beam,
Gamma

ID: 6.8 mm

Volume: 4 mL



DY-035DHPF

Arterial Pillow Chamber

Material: PVC DEHP Free
Color: Clear
Sterilization: Eto, E-Beam
OD: 6.8 mm
Volume: 9 mL



DY-036DHPF

Arterial Pillow Chamber with Valve

Material: PVC DEHP Free
Color: Clear
Sterilization: Eto
ID: 5.5 mm
Volume: 9 mL



DY-039DEHT

Arterial Pillow Chamber with Valve, Service Line Connector

Material: PVC-DEHT
Color: Clear
Sterilization: Eto
ID: 6.8 mm
Volume: 4 mL



DY-041

Recirculating Connector with Vented Cap

Material: PP
Color: Natural
Sterilization: Eto, E-Beam, Gamma, Steam



DY-053

Dialysis T-Connector

Material: PVC
Color: Natural
Sterilization: Eto
ID: 5.5, 6.8 mm



DY-054

Dialysis Cap, Vented

Material: PP
Color: Natural
Sterilization: Eto, E-Beam, Gamma



DY-055

Reducing Adapter

Material: PVC
Color: Clear
Sterilization: Eto
ID: 5.5 to 6.8 mm



DY-056

Pump Connector

Material: PVC DEHP Free
Color: Natural
Sterilization: Eto, E-Beam, Gamma
ID: 4.1, 6.8, 9.54 mm



DY-057

NEW

Dialysis Cap, Vented

Material: HDPE
Color: Natural
Sterilization: Eto, E-Beam, Gamma



DY-058

Dialysis Cap, Vented

Material: LDPE
Color: White
Sterilization: Eto, E-Beam, Gamma



DY-059

Male DIN Adapter

Material: ABS
Color: Clear
Sterilization: Eto, Gamma
ID: 6.8 mm



DY-060DEHT

NEW

Dialysis Cap, Vented

Material: PVC
Color: White
Sterilization: Eto
ID: 6.8 mm



DY-061

Dialysis Cap, Non-Vented

Material: PVC
Color: Clear
Sterilization: Eto
ID: 5.9 mm



DY-143

Discharge Hook, Vented, with Cap

Material: PP
Color: White
Sterilization: Eto, Gamma, Steam



FC-006LGR

Funnel Connector

Material: PVC
Color: Light Green
Sterilization: Eto, Gamma
**ID (French Precision):
6 Fr**



FC-008LBL

Funnel Connector

Material: PVC
Color: Blue
Sterilization: Eto, Gamma
**ID (French Precision):
8 Fr**



FC-010BLK

Funnel Connector

Material: PVC
Color: Black
Sterilization: Eto, Gamma
**ID (French Precision):
10 Fr**



FC-012WHT

Funnel Connector

Material: PVC
Color: White
Sterilization: Eto, E-Beam,
Gamma
**ID (French Precision):
12 Fr**



FC-014GRE

Funnel Connector

Material: PVC
Color: Green
Sterilization: Eto
**ID (French Precision):
14 Fr**



FC-016ORA

Funnel Connector

Material: PVC
Color: Orange
Sterilization: Eto, Gamma
**ID (French Precision):
16 Fr**



FC-018RED

Funnel Connector

Material: PVC
Color: Red
Sterilization: Eto, Gamma
**ID (French Precision):
18 Fr**



FC-020YEL

Funnel Connector

Material: PVC
Color: Yellow
Sterilization: Eto, Gamma
**ID (French Precision):
20 Fr**



FC-022PIN

Funnel Connector

Material: PVC
Color: Pink
Sterilization: Eto, Gamma
**ID (French Precision):
22 Fr**



BC-021A

Large Bore Connector,
Male

Material: ABS
Color: White
Sterilization: Eto, Gamma
ID: 4.8 mm



BC-022

Large Bore Connector,
Male

Material: ABS
Color: White
Sterilization: Eto
ID: 5.5 mm | OD: 6.0 mm



BC-023

Cap for Large Bore
Female Connector, Closed

Material: ABS
Color: White
Sterilization: Eto, Gamma



BC-024

Large Bore Connector,
Female

Material: ABS
Color: White
Sterilization: Eto, Gamma
ID: 3.9 mm
OD: 7.1 mm



BC-024A

Large Bore Connector,
Female

Material: ABS
Color: White
Sterilization: Eto, Gamma
OD: 4.8 mm



BC-025

Large Bore Connector,
Female

Material: ABS
Color: White
Sterilization: Eto, Gamma
ID: 4.5 mm | OD: 7.1 mm



BC-026

Male Luer Cap,
Large Bore

Material: ABS
Color: White
Sterilization: Eto, Gamma
ID: 7.1 mm



BC-083

Large Bore Connector,
Female

Material: ABS
Color: White
Sterilization: Eto, Gamma
ID: 4.4 mm | OD: 6.0 mm



BC-084

Large Bore Connector,
Male

Material: ABS
Color: White
Sterilization: Eto
ID: 6.8 mm | OD: 8.0 mm



BC-085

Large Bore Connector

Material: ABS
Color: White
Sterilization: Gamma
OD: 8.0 mm



TC-005

Suction Connector

Material: PVC
Color: Blue
Sterilization: N/A
ID (French Precision):
30 Fr



TC-008DHPF

Suction Connector

Material: PVC DEHP Free
Color: Blue
Sterilization: Eto
ID (French Precision):
28 Fr



TC-016

Suction Connector

Material: PVC-DEHT
Color: Green
Sterilization: Eto
ID (French Precision):
24 Fr



TC-030A

Suction Connector

Material: PVC DEHP Free
Color: Blue
Sterilization: Eto, Gamma
ID (French Precision):
10 Fr



TC-031A

Suction Connector

Material: PVC DEHP Free
Color: Blue
Sterilization: Eto, Gamma
ID (French Precision):
12 Fr



TC-032A

Suction Connector

Material: PVC DEHP Free
Color: Blue
Sterilization: Eto, Gamma
ID (French Precision):
14 Fr



TC-034B

Suction Connector

Material: PVC DEHP Free
Color: Blue
Sterilization: Eto, Gamma
ID (French Precision):
18 Fr



TC-035A

Suction Connector

Material: PVC DEHP Free
Color: Blue
Sterilization: Eto, Gamma
ID (French Precision):
20 Fr



TC-037A

Suction Connector

Material: PVC DEHP Free
Color: Blue
Sterilization: Eto, Gamma
ID (French Precision):
24 Fr



TC-037G

Suction Connector

Material: PVC DEHP Free
Color: Green
Sterilization: Eto, Gamma
ID (French Precision):
24 Fr



TC-039A

Suction Connector

Material: PVC DEHP Free
Color: Blue
Sterilization: Eto, Gamma
ID (French Precision):
28 Fr



TC-040A

Suction Connector

Material: PVC DEHP Free
Color: Blue
Sterilization: Eto, Gamma
ID (French Precision):
30 Fr



TC-050

Suction Connector

Material: PVC DEHP Free
Color: Green
Sterilization: Eto
ID: 4.1 mm



TC-053DHPF

Suction Connector

Material: PVC DEHP Free
Color: Natural
Sterilization: Eto, E-Beam,
Gamma
ID: 4.1 mm



TC-057

Suction Connector

Material: PVC DEHP Free
Color: Clear
Sterilization: Eto
ID (French Precision):
28 Fr



TC-057R

Suction Connector

Material: PVC DEHP Free
Color: Red
Sterilization: Eto
ID (French Precision):
28 Fr



BB-016

Elbow Connector

Material: PUR, PVC-DEHT
Color: Natural
Sterilization: N/A
OD: 8.0 mm



CH-010

Locking Ring

Material: PP
Color: Clear
Sterilization: Eto, Steam
ID: 7.9 mm



CM-008

Tubing Manifold, 4 to 1

Material: PVC DEHP Free
Color: Clear
Sterilization: Eto
ID: 6.8 mm | OD: 8.4 mm



CM-009

Tubing Manifold, 5 to 1

Material: PVC DEHP Free
Color: Clear
Sterilization: Eto
ID: 6.8 mm | OD: 8.4 mm



CM-010

Tubing Manifold, 8 to 1

Material: PVC DEHP Free
Color: Clear
Sterilization: Eto
ID: 6.8 mm | OD: 8.4 mm



PO-004A

Tubing Port with Cap

Material: PVC DEHP Free
Color: Orange
Sterilization: Eto, Gamma
ID (French Precision): 4 Fr



PO-008A

Tubing Port with Cap

Material: PVC DEHP Free
Color: Orange
Sterilization: Eto, Gamma
ID (French Precision): 8 Fr



TC-009

Tubing Coupler

Material: PVC DEHP Free
Color: Clear
Sterilization: Eto
ID: 5.5 mm



TC-044

Tubing Connector

Material: PVC
Color: Natural
Sterilization: Eto
ID: 6.0 mm



TC-045

Elbow Connector with Slip and Port

Material: PP
Color: Natural
Sterilization: Eto
ID: 8.4 to 6.24 mm



TA-008

T-Connector

Material: ABS
Color: Clear
Sterilization: Eto
ID: 4.1 mm



TA-009

T-Connector

Material: PVC
Color: Clear
Sterilization: Eto
ID: 2.4 mm



TA-011

T-Connector

Material: PVC
Color: Blue, Clear
Sterilization: Eto, Gamma
ID: 2.1 mm



TA-011ABS

T-Connector

Material: ABS
Color: Clear
Sterilization: Eto
ID: 2.1 mm



TA-012

Rotating Hub for
T-Connectors


Material: ABS
Color: Clear
Sterilization: Eto



TA-013

T-Connector with
Luer Slip

Material: ABS
Color: Clear
Sterilization: Eto
ID: 3.5 mm



TA-014

T-Connector with Rotating
Male Luer Lock

Material: ABS
Color: Clear
Sterilization: Eto
ID: 4.1 mm



TA-014PVC

T-Connector with Rotating
Male Luer Lock

Material: ABS, PVC
Color: Clear
Sterilization: Eto
ID: 4.1 mm



TA-015

T-Connector with Rotating
Male Luer Lock

Material: ABS
Color: Clear
Sterilization: Eto
ID: 2.4 mm



TA-016

T-Connector

Material: PVC DEHP Free
Color: Clear
Sterilization: Eto, E-Beam,
Gamma
ID: 5.5 mm



TA-018

T-Connector with Rotating
Male Luer Lock

Material: ABS, PVC
Color: Clear
Sterilization: Eto
ID: 2.4 mm



TA-019

T-Connector with Rotating
Male Luer Lock

Material: ABS
Color: Clear
Sterilization: Eto
ID: 3.5 mm



TA-020

T-Connector

Material: ABS
Color: Clear
Sterilization: Eto, Gamma
ID: 4.1 mm



TA-021

T-Connector with Rotating
Male Luer Lock

Material: ABS, PVC
Color: Blue, Natural
Sterilization: Eto, Gamma
ID: 2.1 mm



TA-024

T-Connector with Rotating
Male Luer Lock and 2
Female Luers

Material: COPE, PVC
Color: Clear
Sterilization: Eto



TC-015

T-Connector with Female
Luer and Rotating Male
Luer Lock

Material: ABS, PP
Color: Clear
Sterilization: Eto, Gamma
ID: 5.5 mm



TT-010

T-Connector with 2 Female
Luers and Male Luer Slip

Material: PC
Color: Clear
Sterilization: Eto, Steam



TT-038

T-Connector with 2 Female
Luers and Rotating Male
Luer Lock

Material: PC
Color: Clear
Sterilization: Eto, Steam
ID: 2.2 mm



CY-001DEHT

Y-Connector

Material: PVC-DEHT
Color: Clear
Sterilization: Eto
ID: 2.5 mm



CY-002DEHT

Y-Connector

Material: PVC-DEHT
Color: Clear
Sterilization: Eto, E-Beam, Gamma
ID: 4.1 mm



CY-002DHPFG

Y-Connector

Material: PVC DEHP Free
Color: Clear
Sterilization: Eto, E-Beam, Gamma
ID: 4.1 mm



CY-002NGDEHT

Y-Connector

Material: PVC-DEHT
Color: Clear
Sterilization: Eto
ID: 4.1 mm



CY-003DEHT

Y-Connector, Manifold

Material: PVC-DEHT
Color: Clear
Sterilization: Eto, E-Beam, Gamma
ID: 5.5 mm



CY-003NGDEHT

Y-Connector, Manifold

Material: PVC-DEHT
Color: Clear
Sterilization: Eto
ID: 5.5 mm



CY-004DEHT

Y-Connector

Material: PVC-DEHT
Color: Clear
Sterilization: Eto, E-Beam, Gamma
ID: 6.8 mm



CY-005A

Y-Connector with Female Luer

Material: PVC
Color: Light Blue
Sterilization: Eto, Gamma
ID: 4.1 mm



CY-005DHPFG

Y-Connector with Female Luer

Material: PVC DEHP Free
Color: Light Blue
Sterilization: Eto, E-Beam, Gamma
ID: 4.1 mm



CY-006A

Y-Connector with Female Luer

Material: PVC DEHP Free
Color: Clear
Sterilization: Eto
ID: 4.1 mm



CY-006C

Y-Connector

Material: PVC-DEHT
Color: Clear
Sterilization: Eto
ID: 4.1 mm



CY-006DEHT

Y-Connector

Material: PVC-DEHT
Color: Clear
Sterilization: Eto, E-Beam, Gamma
ID: 4.1 mm



CY-007DEHT

Y-Connector

Material: PVC-DEHT
Color: Clear
Sterilization: Eto
ID: 3.0 mm



CY-008DEHT

Triple Connector

Material: PVC-DEHT
Color: Clear
Sterilization: Eto, E-Beam, Gamma
ID: 4.1 mm



CY-009DEHT

Triple Connector

Material: PVC-DEHT
Color: Clear
Sterilization: Eto, E-Beam, Gamma
ID: 4.1 mm



CY-009NGDEHT

Triple Connector

Material: PVC-DEHT
Color: Clear
Sterilization: Eto
ID: 4.1 mm



CY-010DEHT

Triple Connector

Material: PVC-DEHT
Color: Clear
Sterilization: Eto
ID: 5.5 mm



CY-010DHPF

Triple Connector

Material: PVC DEHP Free
Color: Clear
Sterilization: Eto, Gamma
ID: 5.5 mm



CY-011DHPF

Triple Connector

Material: PVC DEHP Free
Color: Clear
Sterilization: Eto, Gamma
ID: 6.8 mm



CY-012

Triple Connector

Material: PVC
Color: Light Blue
Sterilization: Eto, Gamma
ID: 6.8 mm



CY-016

Y-Connector with Male and Female Luer Locks

Material: ABS, PP
Color: Clear
Sterilization: Eto, Gamma



CY-017

Y-Connector with Female Luers to Male Luer Slip

Material: ABS
Color: Clear
Sterilization: Eto, Gamma



CY-018

Triple Connector with Rotating Male Luer Lock and 3 Female Luer Locks

Material: ABS, PP
Color: Clear
Sterilization: Eto, Gamma



CY-019

Triple Connector with 3 Female Luers to Male Luer Slip

Material: ABS
Color: Clear
Sterilization: Eto, Gamma



CY-019W

Triple Connector with 3 Female Luers to Male Luer Slip

Material: ABS
Color: White
Sterilization: Eto, Gamma



CY-020DEHT

Y-Connector

Material: PVC-DEHT
Color: Clear
Sterilization: Eto
ID: 4.1 mm



CY-020DHPF

Y-Connector

Material: PVC DEHP Free
Color: Clear
Sterilization: Eto, E-Beam, Gamma
ID: 4.1 mm



CY-020DHPFG

Y-Connector

Material: PVC DEHP Free
Color: Clear
Sterilization: Eto, E-Beam, Gamma
ID: 4.1 mm



CY-020GSDEHT

Y-Connector

Material: PVC-DEHT
Color: Clear
Sterilization: Eto, E-Beam, Gamma
ID: 4.1 mm



CY-021DEHT

Y-Connector

Material: PVC-DEHT
Color: Clear
Sterilization: Eto
ID: 5.5 mm



CY-022DHPF

Y-Connector

Material: PVC DEHP Free
Color: Clear
Sterilization: Eto, E-Beam,
Gamma

ID: 6.8 mm



CY-022DHPFGS

Y-Connector

Material: PVC DEHP Free
Color: Light Blue
Sterilization: Eto, E-Beam,
Gamma

ID: 6.8 mm



CY-023DEHPFREE

Y-Connector

Material: PVC DEHP Free
Color: Clear
Sterilization: Eto, E-Beam,
Gamma

ID: 6.8, 9.5 mm



CY-026

Y-Connector

Material: PVC
Color: Blue
Sterilization: Eto, Gamma
ID: 3.35 mm



CY-027DHPF

Y-Connector

Material: PVC DEHP Free
Color: Clear
Sterilization: Eto, E-Beam,
Gamma

ID: 4.1 mm



CY-028

Y-Connector with
Luer Slips

Material: MABS
Color: Clear
Sterilization: Eto
OD: 4.21 mm



CY-030

Y-Connector

Material: PVC
Color: Clear
Sterilization: Eto
ID: 4.1 mm



CY-031

Y-Connector

Material: PVC
Color: Clear
Sterilization: Eto
ID: 2.1 mm



CY-031ABS

Y-Connector

Material: ABS
Color: Clear
Sterilization: Eto, Gamma
ID: 2.1 mm



CY-033

Triple Connector

Material: PVC
Color: Clear
Sterilization: Eto
ID: 2.1 mm



CY-033ABS

Triple Connector

Material: ABS
Color: Clear
Sterilization: Eto
ID: 0.078"



CY-035

Y-Connector with Rotating
Male Luer Lock

Material: ABS
Color: Clear
Sterilization: Eto, Gamma
ID: 2.1 mm



CY-036

Y-Connector with Rotating
Male Luer Lock

Material: ABS
Color: Clear
Sterilization: Eto, Gamma
ID: 2.1 mm



CY-038

Y-Connector with
Male Slip

Material: ABS
Color: Clear
Sterilization: Eto
ID: 2.2 - 2.4 mm



CY-039

Y-Connector

Material: PVC DEHP Free
Color: Clear
Sterilization: Eto
ID: 6.7 mm



CY-040

Y-Connector with
Female Luer

Material: ABS
Color: Clear
Sterilization: Eto, Gamma
ID: 4.1 mm



CY-041

Y-Connector with
Female Luer

Material: ABS
Color: Clear
Sterilization: Eto, Gamma
ID: 3.0 mm



CY-041-3

Y-Connector with
Female Luer

Material: PVC
Color: Clear
Sterilization: Eto
ID: 3.0 mm



CY-041PVC

Y-Connector with
Female Luer

Material: PVC DEHP Free
Color: Clear
Sterilization: Eto
ID: 4.1 mm



CY-042

Y-Connector with
Female Luers

Material: PVC
Color: Clear
Sterilization: Eto
ID: 4.1 mm



CY-043

Y-Connector with
Female Luers

Material: PVC
Color: Clear
Sterilization: Eto
ID: 3.5 mm



CY-044

Y-Connector with
Female Luers

Material: PVC
Color: White
Sterilization: Eto
ID: 2.67 mm



CY-045

Y-Connector

Material: PVC
Color: Light Blue
Sterilization: Eto, Gamma
ID: 3.35 mm



CY-046

Y-Connector with Rotating
Male Luer Lock

Material: ABS, PVC
Color: Clear
Sterilization: Eto
ID: 2.0 mm



CY-047

Y-Connector

Material: ABS
Color: Clear
Sterilization: Eto
ID: 3.5 mm



CY-057

Y-Connector with Nut

Material: ABS
Color: Clear
Sterilization: Eto
ID: 2.4, 5.1 mm



CY-058

Quadruple Connector
with Nut

Material: ABS
Color: Clear
Sterilization: Eto
ID: 2.4 mm



CY-059

Triple Connector with Nut

Material: ABS
Color: Clear
Sterilization: Eto
ID: 2.4 mm



CY-060

Y-Connector, Parallel

Material: PVC
Color: Clear
Sterilization: Eto
ID: 4.1 mm



CY-061

Y-Connector, Adjustable

Material: Plastisol
Color: Clear
Sterilization: Eto
ID: 0.205"



CY-064

Y-Connector with Female Luer Lock

Material: ABS
Color: Clear
Sterilization: Eto, E-Beam, Gamma

ID: 4.1 mm



CY-065

NEW

Y-Connector

Material: PVC
Color: Clear
Sterilization: Eto
ID: 3.5 mm



CY-077

Triple Connector

Material: PVC
Color: Clear
Sterilization: Eto



CY-079

Y-Connector to Male Luer Slip

Material: PVC
Color: Clear
Sterilization: Eto
ID: 2.4 mm



CY-092

Y-Connector to Male Luer Slip

Material: ABS
Color: Clear
Sterilization: Eto, Gamma
ID: 2.3 mm



CY-092M

Y-Connector to Male Luer Slip

Material: ABS
Color: Clear
Sterilization: Eto
ID: 2.4 mm



CY-092PVC

Y-Connector to Male Luer Slip

Material: PVC DEHP Free
Color: Light Blue
Sterilization: Eto, Gamma
ID: 2.4 mm



CY-093PVC

Y-Connector

Material: PVC
Color: Light Blue
Sterilization: Eto, Gamma
ID: 3.5 mm



CY-101

Quadruple Connector

Material: ABS
Color: Clear
Sterilization: Eto, Gamma
ID: 4.1 mm



CY-161

Y-Connector, Adjustable

Material: Platisol
Color: Clear
Sterilization: Eto
ID: 0.205"



DC-006

Transfusion Chamber with Filter, Non-Vented, **20 Drop**

Material: ABS, PA, PP, PVC DEHP Free
Color: Clear, White
Sterilization: Eto
ID: 4.1 mm
Microns: 200 µm



DC-007DHPF

Drip Chamber with Filter, Non-Vented, **20 Drop**

Material: ABS, PES, PVC DEHP Free
Color: Clear, White
Sterilization: Eto
ID: 4.1 mm
Microns: 150 µm



DC-012

Drip Chamber with Cap and Solution Filter, Non-Vented

Material: ABS, LDPE, PA, PVC DEHP Free
Color: Clear, Natural, White
Sterilization: Eto
ID: 4.1 mm | OD: 4.8 mm
Microns: 15 µm



DC-013F

Drip Chamber Ball

Material: PP
Color: Blue
Sterilization: Eto
Dimensions: 3/8"



DC-018

Infusion Chamber, Non-Vented

Material: PVC
Color: Clear
Sterilization: Eto
ID: 4.1 mm
Volume: 9 mL
Dimensions: 16 x 60 mm



DC-024

Drip Chamber with Filter, Non-Vented

Material: PA, PVC
Color: Clear
Sterilization: Eto
ID: 4.1 mm
Microns: 15 µm



DC-025

Drip Chamber, Non-Vented

Material: PVC
Color: Clear
Sterilization: Eto
ID: 4.1 mm



DC-026

Drip Chamber, Non-Vented

Material: PVC DEHP Free
Color: Clear
Sterilization: Eto, E-Beam, Gamma
Volume: 20 mL



DC-031A

Drip Chamber, Non-Vented

Material: PVC
Color: Clear
Sterilization: Eto
ID: 4.1 mm
Volume: 13 mL



DC-031NG

Drip Chamber, Non-Vented

Material: PVC
Color: Natural
Sterilization: Eto
ID: 4.1 mm
Volume: 13 mL



DC-034

Expansion Chamber, Non-Vented

Material: PVC
Color: Clear
Sterilization: Gamma



DC-039NDGS

Drip Chamber with Cap, Non-Vented, **20 Drop**

Material: LDPE, PVC DEHP Free
Color: Clear, White
Sterilization: Eto, Gamma
OD: 3.7 mm



DC-040DEHT

Drip Chamber,
Non-Vented

Material: PVC-DEHT
Color: Clear
Sterilization: Eto
ID: 6.8 mm
Volume: 57 mL



DC-046DFG

Drip Chamber with Cap,
Non-Vented, **20 Drop**

Material: LDPE, PVC
DEHP Free
Color: Clear, White
Sterilization: Eto, Gamma



DC-051DEHT

Drip Chamber, Non-
Vented, **20 Drop**

Material: PVC-DEHT
Color: Clear
Sterilization: Eto
ID: 6.8 mm
Volume: 16 mL



DC-064

Hand Pump Barrel

Material: PVC DEHP Free
Color: Light Blue
Sterilization: Eto, Gamma



DC-079DEHPFREE

Drip Chamber with Cap,
Non-Vented, **10 Drop**

Material: ABS, LDPE, PVC
DEHP Free
Color: Blue, Clear, White
Sterilization: Eto, Gamma
ID: 4.1 mm



DC-087DEHT

Drip Chamber,
Non-Vented

Material: PVC-DEHT
Color: Clear
Sterilization: Eto
ID: 6.8 mm
Volume: 44 mL



DC-115

Drip Chamber,
Non-Vented

Material: PVC
Color: Natural
Sterilization: Eto
ID: 6.8 mm
Volume: 25 mL



DC-120

Drip Chamber with Cap,
Non-Vented, **10 Drop**

Material: ABS, LDPE, PVC
Color: Clear, White
Sterilization: Eto
ID: 4.1 mm



DC-122

Drip Chamber,
Non-Vented

Material: PVC DEHP Free
Color: Clear
Sterilization: Eto
ID: 6.8 mm
Volume: 22 mL



DC-131

Drip Chamber with Cap,
Non-Vented, **20 Drop**

Material: ABS, LDPE, PVC
DEHP Free
Color: Natural
Sterilization: Eto
ID: 4.1 mm



DC-134

Transfusion Chamber with
Cap and Blood Filter,
Non-Vented, **20 Drop**

Material: ABS, LDPE, PA,
PP, PVC
DEHP Free
Color: Clear, Natural, White
Sterilization: Eto, Gamma
ID: 3.1 mm
Microns: 200 µm



DC-140

Drip Chamber, Wall Drip,
Non-Vented

Material: PVC-DEHT
Color: Clear
Sterilization: Eto
ID: 6.8 mm



DC-142

Drip Chamber,
Non-Vented

Material: PVC DEHP Free
Color: Clear
Sterilization: Eto, Gamma
OD: 5.5 mm



DC-143

Drip Chamber,
Non-Vented

Material: PVC DEHP Free
Color: Clear
Sterilization: Eto, Gamma
ID: 6.8 mm



DC-144

Drip Chamber,
Non-Vented

Material: PVC DEHP Free
Color: Clear
Sterilization: Eto, Gamma
ID: 4.8 mm | OD: 4.1 mm



DC-145

Drip Chamber,
Non-Vented

Material: PVC DEHP Free
Color: Clear
Sterilization: Eto, Gamma
ID: 3.0 mm



DC-147

Arterial Bubble Catcher
with 2 Ports, Non-Vented

Material: PVC
Color: Clear
Sterilization: Eto
OD: 6.8 - 5.5 mm



DC-152

Drip Chamber with 2-Way
Cover, Non-Vented

Material: PVC DEHP Free
Color: Light Blue, Natural
Sterilization: Eto, E-Beam,
Gamma
OD: 2 x 5.5, 1 x 6.8 mm



DC-154

Drip Chamber with 3-Way
Cover and Conical Filter,
Non-Vented

Material: HDPE, PVC
DEHP Free
Color: Clear, Light Blue
Sterilization: Eto, E-Beam,
Gamma
OD: 2 x 5.5, 1 x 6.8 mm
Microns: 270 µm



DC-156

Drip Chamber,
Non-Vented

Material: PVC DEHP Free
Color: Clear
Sterilization: Eto, E-Beam,
Gamma
ID: 4.1 mm | OD: 4.8 mm



DC-158

Drip Chamber with Cap
and Solution Filter, Non-
Vented, **20 Drop**

Material: LDPE, PA, PVC
DEHP Free
Color: Clear, White
Sterilization: Eto, Gamma



DC-160

Hand Pump Chamber with
2 Check Balls. Non-Vented

Material: PP, PVC
Color: Blue, Natural, White
Sterilization: Eto, Gamma
ID: 4.1 mm | OD: 8.2 mm



DC-161

Drip Chamber with 3 Ports
and Conical Filter,
Non-Vented

Material: PP, PVC
DEHP Free
Color: Natural
Sterilization: Eto
ID: 2 x 5.5, 2 x 6.8 mm
Microns: 270 µm



DC-161B

NEW

Drip Chamber with 3 Ports
and Conical Filter,
Non-Vented

Material: HDPE, PVC-
DEHT
Color: Clear, Natural
Sterilization: Eto
ID: 2 x 5.5, 1 x 6.8 mm



DC-162

Drip Chamber,
Non-Vented

Material: PVC DEHP Free
Color: Natural
Sterilization: Eto, E-Beam,
Gamma

ID: 6.8 mm
Volume: 25 mL



DC-167

Drip Chamber, In-Line,
PVC-Free, Non-Vented,
20 Drop

Material: ABS, PP,
PVC-Free
Color: Clear, White
Sterilization: Eto
ID: 3.0 mm



DC-169

Blood Collection Chamber,
Non-Vented

Material: PP
Color: Clear
Sterilization: Eto
ID: 6.8 mm



DC-171

Drip Chamber,
Non-Vented

Material: PVC DEHP Free
Color: Clear
Sterilization: Eto, E-Beam,
Gamma



DC-172

Microdrip Chamber with
Cap and Solution Filter,
Non-Vented, **60 Drop**

Material: ABS, LDPE, PA,
PE, PVC
DEHP Free
Color: Clear, White
Sterilization: Eto, Gamma



DC-174

Drip Chamber with 2 Inlets
and Conical Filter,
Non-Vented

Material: PP, PVC, PVC-
DEHT
Color: Clear
Sterilization: Eto
OD: 2 x 5.5, 2 x 6.8 mm
Microns: 270 µm
Volume: 30 mL



DC-175

Drip Chamber,
Non-Vented

Material: PP
Color: Clear
Sterilization: Eto
ID: 3.0 mm



DC-177DEHT

Drip Chamber,
Non-Vented

Material: PVC-DEHT
Color: Clear
Sterilization: Eto
ID: 6.8 mm
Volume: 62 mL



DC-179

Drip Chamber with Cap
and Solution Filter, Non-
Vented, **20 Drop**

Material: LDPE, PA, PE,
PVC DEHP Free
Color: Clear, White
Sterilization: Eto, Gamma



DC-180

Drip Chamber, Non-
Vented, **20 Drop**

Material: ABS, PVC
DEHP Free
Color: Clear, White
Sterilization: Eto, Gamma
ID: 4.1 mm



DC-186

Drip Chamber with
Solution Filter,
Non-Vented

Material: ABS, PA, PVC
DEHP Free
Color: Clear, White
Sterilization: Eto
ID: 4.1 mm | OD: 4.8 mm
Microns: 15 µm



DC-187

Drip Chamber Assembly
with Cap, Non-Vented,
20 Drop

Material: ABS, LDPE, PVC
DEHP Free
Color: Clear, White
Sterilization: Eto



DC-188

Drip Chamber,
Non-Vented

Material: PVC-DEHT
Color: Clear
Sterilization: Eto
ID: 4.1 mm



DC-192

Expansion Chamber,
Non-Vented

Material: PVC-DEHT
Color: Clear
Sterilization: Eto
ID: 1 x 5.5, 2 x 6.8 mm
Volume: 25 mL



DC-193

NEW

Drip Chamber with Cap,
Non-Vented, **20 Drop**

Material: ABS, PVC
DEHP Free
Color: Clear, White
Sterilization: Eto, E-Beam,
Gamma
ID: 4.1 mm



DC-196

NEW

Drip Chamber, Non-
Vented, with 4 Ports and
Conical Filter

Material: HDPE, PVC
DEHT
Color: Clear
Sterilization: Eto
ID: 3 x 5.5, 1 x 6.8 mm



DC-228G

Drip Chamber,
Non-Vented

Material: PVC DEHP Free
Color: Clear
Sterilization: Eto, Gamma
ID: 4.1 mm



DC-004

Drip Chamber with Air and Solution Filters, Vented, **20 Drop**

Material: LDPE, Nylon, PA, PP, PVC DEHP Free
Color: Clear, Natural, Red, White
Sterilization: Eto, Gamma
ID: 2.7 mm
Microns: 3 µm (Air Filter)
15 µm (Sol. Filter)



DC-015

Drip Chamber with Cap, Air and Solution Filters, Vented, **20 Drop**

Material: Acrylic, LDPE, PE
Color: Clear, Natural, White
Sterilization: Eto, Gamma
Microns: 1.2 µm (Air Filter)
15 µm (Sol. Filter)



DC-020DEHT

Drip Chamber with Cap and Air Filter, Vented, **20 Drop**

Material: ABS, LDPE, PP, PVC-DEHT, PVDF
Color: Clear, White
Sterilization: Eto
OD: 4.1 mm
Microns: 3 µm



DC-021DEHT

Drip Chamber with Cap, Air and Solution Filters, Vented, **20 Drop**

Material: ABS, Acrylic, LDPE, PA, PP, PVDF, PVC-DEHT
Color: Clear, White
Sterilization: Eto
ID: 4.1 mm
Microns: 3 µm (Air Filter)
15 µm (Sol. Filter)



DC-032

Burette, **20 Drop**

Material: ABS, COPE, PE, PIP, PP, PVDF
Color: Blue, Clear, White
Sterilization: Eto
Volume: 150 mL



DC-033

Burette

Material: ABS, Nylon, PE, PES, PETG, PP, PVC, PVC DEHP Free, Rubber, Silicone
Color: Blue, Clear, Natural, White
Sterilization: Eto
Volume: 150 mL



DC-054

Drip Chamber with Cap, Air and Solution Filters, Vented, **20 Drop**

Material: Acrylic, LDPE, PA, PP, PVC DEHP Free
Color: Blue, Clear, Natural, White
Sterilization: Eto, Gamma
ID: 3.0 mm
Microns: 3 µm (Air Filter)
15 µm (Sol. Filter)



DC-054G

Drip Chamber with Cap, Air and Solution Filters, Vented, **20 Drop**

Material: ABS, LDPE, PA, PP, PVC DEHP Free, PVDF
Color: Blue, Clear, Natural, White
Sterilization: Eto, Gamma
Microns: 3 µm (Air Filter)
15 µm (Sol. Filter)



DC-069

Drip Chamber with Cap and Air Filter, Vented, **20 Drop**

Material: ABS, LDPE, PP, PVC, PVDF
Color: Blue, Clear, White
Sterilization: Eto
OD: 4.1 mm
Microns: 3 µm



DC-095

Drip Chamber with Cap and Air Filter, Vented, **15 Drop**

Material: ABS, Acrylic, HDPE, LDPE, PVC DEHP Free
Color: Clear, White
Sterilization: Eto, E-Beam, Gamma
ID: 4.1 mm
Microns: 1.2 µm



DC-097GDF20

Drip Chamber with Cap and Air Filter, Vented, **20 Drop**

Material: Acrylic, LDPE, PP, PVC DEHP Free
Color: Clear, Red, White
Sterilization: Eto, Gamma
OD: 2.7 mm
Microns: 3 µm



DC-126

Drip Chamber with Cap and Air Filter, Vented, **20 Drop**

Material: Acrylic, LDPE, PP, PVC DEHP Free
Color: Clear, Natural, Red, White
Sterilization: Eto, Gamma
ID: 3.0 mm
Microns: 3 µm



DC-130

Drip Chamber with Cap and Air Filter, Vented, **10 Drop**

Material: ABS, Acrylic, LDPE, Nylon, PVC

Color: Clear, Natural, White

Sterilization: Eto

ID: **4.1 mm**

Microns: 0.8 µm



DC-133

Transfusion Chamber, High Flow, with Cap, Air and Blood Filters, Vented, **20 Drop**

Material: ABS, LDPE, PA, PP, PVC

Color: Clear, Green, White

Sterilization: Eto

ID: **3.0 mm**

Microns: 0.3 µm (Air Filter)
200 µm (Blood Filter)



DC-137

Drip Chamber with Cap and Air Filter, Vented, **20 Drop**

Material: ABS, LDPE, Nylon, PVC
DEHP Free

Color: Blue, Clear, White

Sterilization: Eto

ID: **4.1 mm**

Microns: 0.5 µm



DC-138

Drip Chamber with Cap and Air Filter, Vented, **60 Drop**

Material: ABS, PE, PVC
DEHP Free

Color: Blue, Clear, White

Sterilization: Eto

ID: **4.1 mm**

Microns: 0.5 µm



DC-159

Drip Chamber with Cap, Vented, **60 Drop**

Material: ABS, LDPE, PVC DEHP Free

Color: Clear, White

Sterilization: Eto, Gamma



DC-163

Drip Chamber with Cap, Air and Solution Filters, Vented, **20 Drop**

Material: Acrylic, LDPE, Nylon, PA, PP, PVC DEHP Free

Color: Clear, Red, White

Sterilization: Eto, Gamma

Microns: 3 µm (Air Filter)
15 µm (Sol. Filter)



DC-165

Drip Chamber with Cap, Air and Solution Filters, Vented, **20 Drop**

Material: ABS, PA, PE, PP, PVC-Free, PVDF

Color: Clear, Green, Natural, White

Sterilization: Eto

ID: **3.0 mm**

Microns: 0.8 µm (Air Filter)
15 µm (Sol. Filter)



DC-166

Drip Chamber with Cap, Air and Solution Filters, Vented, **60 Drop**

Material: ABS, PE, PP, PVC-Free

Color: Green, Natural, White

Sterilization: Eto

ID: **3.0 mm**

Microns: 0.8 µm (Air Filter)
15 µm (Sol. Filter)



DC-181

Microdrip Chamber with Cap, Air and Solution Filters, Vented, **60 Drop**

Material: ABS, Acrylic, LDPE, Nylon, PA, PP, PVC
DEHP Free

Color: Clear, Red, White

Sterilization: Eto, Gamma

Microns: 3 µm (Air Filter)
15 µm (Sol. Filter)



DC-195

NEW

Drip Chamber with Cap and Air Filter, Vented, **10 Drop**

Material: ABS, Glass Fiber, LDPE, PVC
DEHP Free

Color: Blue, Clear, White

Sterilization: Eto

Microns: 0.5 µm



DC-197-60

Microdrip Chamber with Cap and Air Filter, Vented, **60 Drop**

Material: ABS, Acrylic, LDPE, Nylon, PE, PP, PVC
DEHP Free

Color: Clear, Red, White

Sterilization: Eto, Gamma

Microns: 3 µm



SP-042

Drip Chamber with Cap and Air Filter, Vented, **10 Drop**

Material: ABS, Acrylic, HDPE, LDPE, PVC DEHP Free

Color: Clear, Natural, White

Sterilization: Eto, Gamma

Microns: 1.2 µm



CA-155

Burette Cap with Air Vent,
Needle-Free

Material: ABS, PVC
Color: Blue, White
Sterilization: Eto



DC-008ADF

Drip Chamber Cover

Material: PVC DEHP Free
Color: Natural
Sterilization: Eto
ID: 4.1 mm



DC-008DF

Drip Chamber Cover,
20 Drop

Material: PVC DEHP Free
Color: Natural
Sterilization: Eto, E-Beam,
Gamma
ID: 4.1 mm



DC-008NG

Drip Chamber Cover,
20 Drop

Material: PVC DEHP Free
Color: Natural
Sterilization: Eto
ID: 4.1 mm



DC-062

Hand Pump Cover,
Outflow, White Check Ball

Material: PP, PVC
Color: Natural, White
Sterilization: Eto, Gamma
ID: 4.1 mm



DC-063

Hand Pump Cover,
Inflow, Blue Check Ball

Material: PP, PVC
Color: Blue, Clear
Sterilization: Eto, Gamma
ID: 4.1 mm



DC-068DEHT

Drip Chamber Cover
with 4 Ports

Material: PVC-DEHT
Color: Clear
Sterilization: Eto
ID: 3 x 5.5, 1 x 6.8 mm



DC-076

Drip Chamber Cover
with 2 Ports

Material: PVC
Color: Natural
Sterilization: Eto
ID: 4.1 mm



DC-100

Drip Chamber Cover

Material: PVC DEHP Free
Color: White
Sterilization: Eto, Gamma
ID: 2.0 mm | OD: 3.0 mm



DC-123

Drip Chamber Cover
with 2 Ports

Material: PVC DEHP Free
Color: Purple
Sterilization: Eto, Steam
ID: 4.1 mm



DC-127

Drip Chamber Cover
with 3 Ports

Material: PVC
Color: Purple
Sterilization: Eto
ID: 4.1, 5.5, 6.8 mm



DC-128DFG

Drip Chamber Cover

Material: PVC DEHP Free
Color: Light Blue
Sterilization: Eto, Gamma
ID: 4.1 mm | OD: 5.5 mm



DC-128DFNG

Drip Chamber Cover

Material: PVC DEHP Free
Color: Clear
Sterilization: Eto
ID: 4.1 mm | OD: 5.5 mm



DC-141

Drip Chamber Cover
with 3 Ports

Material: PVC
Color: Purple
Sterilization: Eto
ID: 2 x 5.5, 1 x 6.8 mm



DC-146

Drip Chamber Cover
with 3 Ports

Material: PVC DEHP Free
Color: Clear
Sterilization: Eto, Gamma
ID: 2 x 5.5, 1 x 6.8 mm



DC-168

Drip Chamber Cover,
One-Way

Material: PVC
Color: Natural
Sterilization: Eto
OD: 4.1 mm



DC-176

Drip Chamber Cover

Material: ABS
Color: White
Sterilization: Eto
ID: 3.0 mm



DC-178

Drip Chamber Cover

Material: ABS
Color: White
Sterilization: Eto, Gamma
ID: 4.1 mm



DC-189

Drip Chamber Cover
with 3 Ports

Material: PVC DEHP Free
Color: Clear
Sterilization: Eto
ID: 2 x 5.5, 1 x 6.8 mm



DC-190

Drip Chamber Cover
with 3 Ports

Material: PVC DEHP Free
Color: Clear
Sterilization: Eto
ID: 2 x 5.5, 1 x 6.8 mm



DC-198

NEW

Drip Chamber Cover

Material: PVC DEHP Free
Color: Light Blue
Sterilization: Eto, E-Beam,
Gamma
ID: 6.8 mm



BC-028

ENFit System Bottle
Adapter, Vented

Material: HDPE, LDPE, PS,
Silicone
Color: Purple, White
Sterilization: Eto, E-Beam,
Gamma

ID: 3.7 mm

Pressure: ≤ 8 mbar



BC-036

ENFit System Cap, Vented

Material: LDPE
Color: Natural
Sterilization: Eto, Gamma
ID: 10.00 mm



EN-007

ENFit System Cap Tether

Material: PE
Color: Clear
Sterilization: Eto, Gamma



EN-0080

ENFit System Female Luer

Material: PVC
Color: Orange
Sterilization: Eto
ID: 3.5 mm



EN-008P

ENFit System Female Luer

Material: PVC
Color: Purple
Sterilization: Eto
ID: 3.5 mm



EN-0090

ENFit System Male Luer

Material: PVC
Color: Orange
Sterilization: Eto
ID: 3.5 mm



EN-0100R

ENFit System Female Luer

Material: PVC
Color: Orange
Sterilization: Eto
ID: 4.1 mm



EN-010P

ENFit System Female Luer

Material: PVC
Color: Purple
Sterilization: Eto
ID: 4.1 mm



EN-010W

ENFit System Female Luer

Material: PVC
Color: White
Sterilization: Eto
ID: 4.1 mm



EN-0110R

ENFit System Male Luer

Material: PVC
Color: Orange
Sterilization: Eto
ID: 4.1 mm



EN-011P

ENFit System Male Luer

Material: PVC
Color: Purple
Sterilization: Eto
ID: 4.1 mm



EN-0150R

ENFit System Female Luer

Material: PVC
Color: Orange
Sterilization: Eto
ID: 2.1 mm



EN-015P

ENFit System Female Luer

Material: PVC
Color: Purple
Sterilization: Eto
ID: 2.0 mm



EN-0160R

ENFit System Male Luer

Material: PVC
Color: Orange
Sterilization: Eto
ID: 2.0 mm



EN-016P

ENFit System Male Luer

Material: PVC
Color: Purple
Sterilization: Eto
ID: 2.0 mm



EN-018

ENFit System Tethered
Cap for Male Luers

Material: PVC
Color: Clear, Purple
Sterilization: Eto, Gamma
ID: 5.5 mm



EN-018B

ENFit System Tethered
Cap for Male Luers

Material: PVC DEHP Free
Color: Blue, Clear
Sterilization: Eto, Gamma
ID: 5.5 mm



EN-019P

ENFit System Cap for
Female Luers, Vented

Material: MABS
Color: Purple
Sterilization: Eto, Gamma
OD: 5.4 mm



EN-019PA

ENFit System Cap for
Female Luers, Vented

Material: ABS
Color: Purple
Sterilization: Eto, Gamma
OD: 5.4 mm



EN-020P

ENFit System Male
Luer Cap

Material: PVC DEHP Free
Color: Purple
Sterilization: Eto



EN-020PA

ENFit System Cap for
Male Luers, Vented

Material: ABS
Color: Purple
Sterilization: Eto, Gamma



EN-020PT

ENFit System Cap on
Strap for Male Luers

Material: ABS, PE
Color: Purple
Sterilization: Eto



EN-0210

ENFit System Female Luer

Material: PVC
Color: Orange
Sterilization: Eto
ID: 2.4 mm



EN-021P

ENFit System Female Luer

Material: PVC
Color: Purple
Sterilization: Eto
ID: 2.4 mm



EN-023

ENFit System Cap for
Female Luers, Vented

Material: HDPE
Color: Natural
Sterilization: Eto



EN-024

ENFit System Cap for
Female Luers, Non-Vented

Material: HDPE
Color: Natural
Sterilization: Eto



EN-026

ENFit System Transition /
Syringe Adapter

Material: ABS
Color: Purple
Sterilization: Eto



EN-029

ENFit System Stepped
Adapter

Material: ABS
Color: Purple
Sterilization: Eto, Gamma



EN-033

ENFit System Connector,
Male Luer Lock to Port

Material: ABS
Color: White
Sterilization: Eto, Gamma
ID: 5.5 mm



EN-035

ENFit System Female
Connector with Cap

Material: HDPE, PVC
DEHP Free
Color: Purple
Sterilization: Eto
ID: 4.1 mm



EN-043

ENFit System IV Spike
with Cap, Non-Vented

Material: ABS, LDPE
Color: Clear, Purple
Sterilization: Eto
ID: 4.1 mm



DY-011DEHT

Transducer Protector with Male Luer Lock and Luer Slip

Material: PTFE, PVC-DEHT
Color: Blue, Clear
Sterilization: Eto, E-Beam, Gamma

Microns: 0.2 µm



DY-012

Transducer Protector with Male and Female Luer Locks

Material: PVC
Color: Clear
Sterilization: Eto
Microns: 0.2 µm



DY-013ADHPF NEW

Transducer Protector, Female Luer Lock to Luer Slip

Material: PTFE, PVC
DEHP Free
Color: Clear, Natural, White
Sterilization: Eto, E-Beam, Gamma
Microns: 0.2 µm



DY-014DEHT

Transducer Protector with Male Luer Slips

Material: PTFE, PVC-DEHT
Color: Clear, White
Sterilization: Eto
Microns: 0.2 µm



DY-015

Transducer Protector with Male Luer Slip and Female Luer Lock

Material: PTFE, PVC
Color: Clear, White
Sterilization: Eto
Microns: 0.2 µm



DY-045DEHT

Transducer Protector with Drip Chamber and Male Luer Slip

Material: PTFE, PVC-DEHT
Color: Clear
Sterilization: Eto
ID: 3.5, 6.0 mm
Microns: 0.2 µm



DY-045DHPF

Transducer Protector with Drip Chamber and Male Luer Slip

Material: PTFE, PVC
DEHP Free
Color: Clear
Sterilization: Eto
ID: 3.5 mm
Microns: 0.2 µm



DY-047

Transducer Protector with Male Luer Slips

Material: PVC
Color: Clear
Sterilization: Eto
Microns: 0.2 µm



DY-047A

Transducer Protector with Luer Slips, Positive to Negative

Material: PVC, PVC
DEHP Free
Color: Clear
Sterilization: Eto
Microns: 0.2 µm
Pressure: 87 psi



DY-047DHPFG

Transducer Protector with Luer Slips, Positive to Negative

Material: PVC DEHP Free
Color: Light Blue, Natural
Sterilization: Eto, Gamma
Pressure: 87 psi



DY-052

Transducer Protector with Female Luer Lock and Tubing Port

Material: PTFE, PVC
DEHP Free
Color: Natural, White
Sterilization: Eto
ID: 5.5 mm
Microns: 0.2 µm



DY-062DEHT

Transducer Protector with Male Luer Slip and Female Luer Lock

Material: PTFE, PVC-DEHT
Color: Clear, White
Sterilization: Eto
Microns: 0.2 µm



DY-110

Transducer Protector with Male and Female Luer Locks

Material: PTFE, PVC
Color: Clear, Natural
Sterilization: Eto, Gamma
Microns: 0.2 μ m



DY-110NG

Transducer Protector with Male and Female Luer Locks

Material: PTFE, PVC
Color: Natural
Sterilization: Eto
Microns: 0.2 μ m



FF-003

Inline IV Filter with Female and Rotating Male Luer Locks

Material: Acrylic, PES, PTFE
Color: Clear
Sterilization: Eto, Gamma
Microns: 0.2 μ m
Pressure: 45 psi



FF-003A

Inline IV Filter, Vented, with Female and Rotating Male Luer Locks

Material: Acrylic, PES, PTFE
Color: Clear, White
Sterilization: Eto, Gamma
Microns: 0.2 μ m
Pressure: 46.4 psi



FF-004

Inline IV Filter, Adult

Material: Acrylic, PES, PTFE
Color: Blue, Clear
Sterilization: Eto, Gamma
ID: 4.1 mm
Microns: 0.03 μ m (Vent) 1.2 μ m (Filter)
Pressure: 45 psi



FF-005

Inline IV Filter, Adult

Material: Acrylic, PES, PTFE
Color: Clear, Natural
Sterilization: Eto, Gamma
ID: 3.6 mm
Microns: 5 μ m
Pressure: 46.4 psi



FF-007

Inline IV Filter, Adult

Material: Acrylic, PES, PTFE
Color: Clear
Sterilization: Eto, Gamma
ID: 4.1 mm
Microns: 0.03 μ m (Vent) 0.2 μ m (Filter)
Pressure: 45 psi



FF-008

Inline IV Filter, Pediatric

Material: Acrylic, ePTFE, PES
Color: Clear, Natural
Sterilization: N/A
ID: 2.2 mm
Microns: 0.2 μ m



FF-010

Small Volume Filter with Male Luer Lock, Hydrophobic

Material: ePTFE, PVC
DEHP Free
Color: Clear, Yellow
Sterilization: Eto, Gamma
Microns: 0.2 μ m



FF-030

Inline IV Filter, Adult

Material: Acrylic, PES, PTFE
Color: Clear
Sterilization: Eto, Gamma
ID: 4.1 mm
Microns: 0.03 μ m (Vent) 1.2 μ m (Filter)
Pressure: 3.2 bar



FF-034

Syringe Filter with Male Luer Lock

Material: PES
Color: Blue, Natural
Sterilization: N/A
Microns: 0.2 μ m



FF-034C

Syringe Filter with Male Luer Lock

Material: Acrylic, PES
Color: Clear, White
Sterilization: Eto, Gamma
Microns: 0.22 μ m
Pressure: 80 psi



FF-037

Inline IV Filter, Neonatal,
Positive Charge

Material: Acrylic, PES,
PTFE

Color: Clear, Natural

Sterilization: Eto, Gamma

ID: 2.0 mm

**Microns: 0.03 µm (Vent)
0.2 µm (Filter)**

Pressure: 75.4 psi



FF-038

Inline IV Filter, Neonatal,
Positive Charge

Material: Acrylic, PES,
PTFE

Color: Clear, White

Sterilization: Eto, Gamma

ID: 2.0 mm

Microns: 5 µm

Pressure: 75.4 psi



FF-051


Tubular Filter

Material: ABS, PA

Color: Natural

Sterilization: N/A

Microns: 300 µm



FF-052

Tubular Filter

Material: Nylon, PP

Color: Natural

Sterilization: Eto, E-Beam,
Gamma

Microns: 200 µm



FF-053

Disk Filter

Material: LDPE, PA

Color: White

Sterilization: Eto, Gamma

Microns: 15 µm



FF-054A

Tubular Filter

Material: ABS, PE

Color: Natural

Sterilization: Eto

Microns: 150 µm



FF-055

Tubular Filter

Material: ABS, PA

Color: White

Sterilization: Eto, Gamma

Microns: 200 µm



FF-057

Inline IV Filter

Material: Acrylic, PES,
PTFE

Color: Red, White

Sterilization: Eto, Gamma

ID: 2.0 mm

Microns: 0.2 µm

Pressure: 3.1 bar - 45 psi
Max



FF-059

Air Vent Filter

Material: Acrylic, PP

Color: Red

Sterilization: Eto, Gamma,
Steam

Microns: 3 µm



FF-060

Syringe Filter with Male
Luer Slip

Material: Acrylic, Nylon

Color: Blue, Clear

Sterilization: Eto, Gamma

Microns: 0.2 µm



FF-061

Tubular Filter

Material: PA, PP

Color: Natural

Sterilization: Eto, Gamma

Microns: 200 µm



FF-064

Conical Filter

Material: PP

Color: Natural

Sterilization: Eto, Gamma

Microns: 270 µm



FF-066

Tubular Filter

Material: Nylon, PA, PP
Color: White
Sterilization: Eto
Microns: 263 µm



FF-077

Air Vent Filter

Material: Acrylic, PP
Color: White
Sterilization: Eto, Gamma, Steam
Microns: 3 µm



FF-077B

NEW

Air Vent Filter

Material: Acrylic, PP
Color: Blue, Natural
Sterilization: Eto, Gamma, Steam
Microns: 3 µm



FF-083

Drip Chamber Filter /
Nutrition Filter, Vented,
1.2 µm

Material: COPE, PES,
PVDF
Color: Clear, Yellow
Sterilization: Eto, Gamma
**Microns: 0.1 µm (Vent)
1.2 µm (Filter)**
Dimensions: 3.3 mm



FF-084

Inline IV Filter, Adult,
Female Luer Lock to
Luer Slip

Material: Acrylic, PES,
PTFE
Color: Clear, White
Sterilization: Eto, Gamma
Microns: 0.2 µm
Pressure: 46.4 psi /3.2 bar
Dimensions: 10 cm²
Filtration Area



FF-085

Tubular Filter

Material: ABS, Nylon
Color: Clear
Sterilization: Eto
Microns: 263 µm



FF-086

Inline IV Filter, Epidural,
with Female and Rotating
Male Luer Locks

Material: Acrylic, PES, PP
Color: Clear, White
Sterilization: Eto, Gamma
Microns: 0.2 µm



FF-087

Tubular Filter

Material: ABS, Polyester
Color: White
Sterilization: Eto, Gamma
Microns: 150 µm



FF-088

Disk Filter

Material: LDPE, PA
Color: Natural, White
Sterilization: Eto, Gamma
Microns: 15 µm
Dimensions: 16.2 mm



FF-089

Tubular Filter

Material: Nylon, PA
Color: Natural
Sterilization: Eto, Gamma
Microns: 40 µm



FF-093

Disk Filter

Material: ABS, Nylon
Color: White
Sterilization: Eto, Gamma
Microns: 15 µm



FF-094

Disk Filter

Material: Nylon, PP
Color: Natural, White
Sterilization: Eto, Gamma
Microns: 15 µm
Dimensions: 13 mm



FF-095

Syringe Filter with Male Luer Lock

Material: PP, PTFE
Color: White
Sterilization: N/A
Microns: 0.45 µm



FF-096

Tubular Filter

Material: ABS, Nylon
Color: White
Sterilization: Eto, E-Beam, Gamma
Microns: 170 µm



FF-097

Conical Filter

Material: HDPE
Color: Natural
Sterilization: Eto, E-Beam, Gamma
Microns: 270 µm



FF-098

Inline IV Filter with Male and Female Luer Locks

Material: Acrylic, PVC
DEHP Free
Color: Light Blue, Clear
Sterilization: Eto, Gamma
ID: 4.1 mm
Microns: 0.2 µm



FF-099

Inline IV Filter with Male and Female Luer Locks

Material: Acrylic, PE, PTFE
Color: Blue, Clear
Sterilization: Eto
Microns: 1.2 µm



FF-101

Inline IV Filter with Male and Female Luers

Material: Acrylic, PE, PTFE
Color: Clear
Sterilization: Eto, Gamma
Microns: 1.2 µm



FF-102

Inline IV Filter, Adult

Material: Acrylic, PES, PTFE
Color: Clear
Sterilization: Eto, Gamma
Microns: 0.03 µm (Vent) 1.2 µm (Filter)
Pressure: 3.2 bar



FF-103

Inline IV Filter, Adult, Positive Charge

Material: Acrylic, PES, PTFE
Color: Clear
Sterilization: Eto, Gamma
Microns: 0.03 µm (Vent) 0.2 µm (Filter)
Pressure: 3.2 bar



FF-104

Inline IV Filter, Infant

Material: Acrylic, PES, PTFE
Color: Clear
Sterilization: Eto, Gamma
ID: 3.0 mm
Microns: 0.2 µm



FF-105

Inline IV Filter, Infant

Material: Acrylic, PES, PTFE
Color: Clear
Sterilization: Eto, Gamma
ID: 3.0 mm
Microns: 1.2 µm



FF-106

Inline IV Filter, Infant

Material: Acrylic, PES, PTFE
Color: Clear
Sterilization: Eto, Gamma
ID: 3.0 mm
Microns: 0.2 µm



FF-107

Disk Filter

Material: ABS, Nylon
Color: White
Sterilization: Eto
Microns: 15 µm



FF-108

NEW

Inline IV Filter, Female
Luer to Rotating Male
Luer Lock

Material: Acrylic, PES, PP
Color: Clear, Natural
Sterilization: Eto, Gamma
Microns: 0.2 µm
Pressure: 115 psi



FF-110

Inline IV Filter

Material: COPE, PES,
PVDF
Color: Clear, Natural
Sterilization: Eto, Gamma
ID: 3.3 mm
Microns: 0.1 µm (Vent)
0.2 µm (Filter)
Pressure: 45 psi
Priming Volume: 3 mL



FF-111

Inline IV Filter,
Non-Vented

Material: MBS, PTFE
Color: Clear, White
Sterilization: Eto, Gamma
ID: 2.5 mm
Microns: 0.2 µm



FF-112

Inline IV Filter

Material: Acrylic, PTFE
Color: Clear
Sterilization: Eto
Microns: 0.2 µm



FF-116

Inline IV Filter

Material: Acrylic, PES,
PTFE
Color: Clear, White
Sterilization: Eto, Gamma
ID: 2.7 mm
Microns: 0.2 µm



FF-120

Air Vent Filter, Barbed

Material: Glass, PP
Color: Clear, Green
Sterilization: N/A
OD: 5.9 - 8.0 mm



FF-122

Inline IV Filter, Pediatric

Material: Acrylic, ePTFE,
PES
Color: Clear, Natural
Sterilization: Eto
ID: 1.97 mm
Microns: 0.2 µm



FF-123

Inline IV Filter, Pediatric

Material: Acrylic, ePTFE,
PES
Color: Clear, Natural
Sterilization: Eto
ID: 1.97 mm
Microns: 1.2 µm



FF-124

Inline IV Filter

Material: Acrylic, PES,
Polyester, PTFE
Color: Clear, Natural
Sterilization: Eto, Gamma
ID: 1.45mm | OD: 3.42mm
Microns: 0.2 µm
Pressure: 3.17 bar / 45 psi
(Max)
Priming Volume: 2.9 mL
Dimensions: 36 x 72 mm



FF-127

NEW

Tubular Filter

Material: Nylon, PS
Color: Natural
Sterilization: Eto
Microns: 200 µm



FF-125

Inline IV Filter

Material: Acrylic, PES,
Polyester, PTFE
Color: Blue, Clear, Natural
Sterilization: Eto, Gamma
ID: 1.45mm | OD: 3.42mm
Microns: 1.2 µm
Pressure: 0.3 bar / 5 psi
(Max)
Priming Volume: 2.9 mL
Dimensions: 36 x 72 mm



FF-128

NEW

Inline IV Filter, Adult,
Positive Charge

Material: MABS, PES,
PTFE
Color: Clear, White
Sterilization: Eto
Microns: 0.03 µm (Vent)
0.2 µm (Filter)
Pressure: 3.6 bar (52 psi)



FF-129

NEW

Inline IV Filter, Pediatric,
Positive Charge

Material: MABS, PES,
PTFE

Color: Clear, White

Sterilization: Eto

**Microns: 0.03 µm (Vent)
0.2 µm (Filter)**

Pressure: 3.6 bar (52 psi)



FF-130

Air Vent Filter, Bacterial

Material: PE, PET, PVDF

Color: White

Sterilization: Eto, Gamma

Microns: 3 µm



FF-131

NEW

Inline IV Filter, Infant,
Positive Charge

Material: MABS, PES,
PTFE

Color: Clear, White

Sterilization: Eto

**Microns: 0.03 µm (Vent)
0.2 µm (Filter)**

Pressure: 3.6 bar (52 psi)



FF-141

NEW

Air Vent Filter

Material: Acrylic, LDPE, PP

Color: Natural, Red

Sterilization: Eto

Microns: 1.2 µm



FF-178

NEW

Air Vent Filter

Material: Acrylic, HDPE

Color: Green, Natural

Sterilization: Eto, E-Beam,
Gamma

Microns: 3 µm



FF-224

Air Vent Filter,
Hydrophobic

Material: PP, PVDF

Color: Blue

Sterilization: Eto, Gamma

Microns: 3 µm



FF-226

Air Vent Filter,
Hydrophobic

Material: ABS, Acrylic

Color: White

Sterilization: Eto, Gamma

Microns: 3 µm



FF-301

Air Vent Filter

Material: Acrylic, LDPE,
Nylon

Color: Blue, Natural

Sterilization: Eto, Gamma

Microns: 3 µm



FF-303

Air Vent Filter

Material: Acrylic, PP

Color: Natural

Sterilization: Eto, Gamma

Microns: 3 µm



FF-304

Air Vent Filter with Cap

Material: Acrylic, Nylon,
PP

Color: White

Sterilization: Eto

Microns: 0.3 µm



RF-003

Respiratory Filter,
Bacterial / Viral

Material: K-Resin

Color: Clear

Sterilization: N/A

ID: 15.00 mm

OD: 22.00 mm



RW-009


Syringe Filter with
Luer Lock

Material: PC, PET

Color: Purple

Sterilization: Eto

Microns: 5 µm



RW-010

Syringe Filter with Luer Lock

Material: MABS, PET
Color: Clear
Sterilization: Eto
Microns: 5 µm
Dimensions: 25 mm



RW-017

Inline IV Filter with Luer Slip, Positive Charge

Material: Nylon, PC, PTFE
Color: Clear
Sterilization: Eto
Microns: 0.2 µm



RW-018

Inline IV Filter with Caps, Male and Female Luer Locks

Material: MABS, PET, PTFE
Color: Clear, White
Sterilization: N/A
Microns: 1.2 µm



RW-019

Syringe Filter

Material: MABS, PET
Color: Clear
Sterilization: Eto
Microns: 5 µm
Pressure: 5 bar
Dimensions: 18 mm



RW-020

Inline IV Filter, Pediatric, with Male and Female Luer Locks

Material: MABS, PET, PTFE
Color: Clear
Sterilization: Eto
Microns: 0.2 µm



RW-021

Inline IV Filter, Epidural, with Male and Female Luers

Material: PC, PET
Color: Clear
Sterilization: Eto, Gamma
Microns: 0.2 µm



RW-025

Inline IV Filter

Material: PC, PET, PTFE
Color: Clear, Natural
Sterilization: Eto
ID: 2.0 mm
Microns: 1.2 µm



RW-122

Syringe Filter with Male and Female Luer Locks

Material: MABS, PET
Color: Clear
Sterilization: Eto
Microns: 0.2 µm
Pressure: 4.5 bar



SP-140B

Air Vent Filter with Closed Cap

Material: PP, PVDF
Color: Blue
Sterilization: Eto, Gamma
Microns: 3 µm



FC-010

Flow Indicator

Material: Acrylic, PP
Color: Blue, Clear
Sterilization: Eto
ID: 1/4" - 5/16"



FC-025

NEW

Flow Indicator

Material: Acrylic, PP
Color: Blue, Clear
Sterilization: Gamma
ID: 3/8"



FR-001

Flow Regulator, On / Off

Material: ABS, PC, POM,
Silicone
Color: Blue, Clear
Sterilization: Eto



FR-002

Flow Regulator, Dual Scale

Material: ABS, SEBS,
Silicone
Color: Blue, Clear, White
Sterilization: Eto, Gamma
ID: 3.0 mm
Pressure: 7.25 psi
Flow Rate: 20-250 mL/hr



FR-009

Flow Regulator

Material: ABS, SEBS
Color: Blue, Clear, White
Sterilization: Eto, Gamma
OD: 3.0 mm
Flow Rate: 20-250 mL/hr



FR-014

Flow Regulator

Material: PC, TPE
Color: Blue, White
Sterilization: N/A
Flow Rate: 5-350 mL/hr



FR-016

Flow Regulator

Material: ABS, TPE
Color: Clear, White
Sterilization: Eto
OD: 3.0 mm
Flow Rate: 0-250 mL/hr



FR-017

Flow Regulator, Horizontal

Material: PC, POM, PP
Color: Black, White
Sterilization: Eto, Gamma
ID: 3.0 mm | OD: 4.1 mm
Flow Rate: 10-360 mL/hr



FR-018

Flow Regulator, Dual Scale

Material: ABS, Silicone,
TPE
Color: Green, White
Sterilization: Eto
ID: 3.0 mm | OD: 4.1 mm
Flow Rate: 10-250 mL/hr



FR-019

Flow Regulator

Material: ABS, SEBS,
Silicone
Color: Blue, White
Sterilization: Eto, Gamma
Pressure: 0.5 bar
Flow Rate: 20-250 mL/hr



FR-021

Flow Regulator

Material: ABS, Silicone
Color: Blue, White
Sterilization: Eto, Gamma
ID: 3.0 mm | OD: 4.1 mm
Flow Rate: 5-250 mL/hr



FR-025

Flow Regulator

Material: ABS, Silicone
Color: Blue, White
Sterilization: Eto, Gamma
ID: 3.0 mm | OD: 4.1 mm
Pressure: 50 kPa / 7.252 psi
Flow Rate: 5-300 mL/hr



FR-030

Precision Flow Restrictor

Material: Nylon, PC,
Silicone
Color: Tinted, White
Sterilization: N/A
Volume: 100 mL
Pressure: 55 kPa / 7.98 psi
Flow Rate: 0.44-0.61 mL/hr



IS-001

Y-Injection Site, Dual Access

Material: MABS, PIP
Color: Clear, Natural
Sterilization: Eto
ID: 4.1 mm



IS-004

Y-Injection Site

Material: PVC
Color: Clear, Natural
Sterilization: Eto, Gamma
ID: 3.4 mm



IS-006

T-Injection Site

Material: PVC
Color: Red
Sterilization: Eto
ID: 6.5 mm



IS-010

T-Injection Site with Rotating Male Luer Lock

Material: IR, MABS, PVC
Color: Clear, Yellow
Sterilization: Eto, Gamma
ID: 2.4 mm



IS-010C

T-Injection Site with Cap and Rotating Male Luer Lock

Material: ABS, HDPE, IR, PVC
Color: Clear, Natural, White
Sterilization: Eto
ID: 2.4 mm



IS-011

Injection Stopper

Material: ABS, PIP
Color: Yellow
Sterilization: Eto
ID: 4.4 mm



IS-011CL

Injection Stopper

Material: ABS, IR
Color: Clear, Natural
Sterilization: Eto



IS-014

T-Injection Site with Rotating Male Luer Lock

Material: ABS, IR, PVC
Color: Clear, Yellow
Sterilization: Eto
ID: 4.1 mm



IS-014S

T-Injection Site with Vented Cap and Male Luer Slip

Material: ABS, HDPE, IR, PVC
Color: Clear, White, Yellow
Sterilization: Eto
ID: 4.1 mm



IS-018

Injection Site, Microbore

Material: ABS, IR
Color: Clear, Yellow
Sterilization: Eto, Gamma
ID: 2.1 mm



IS-019

Injection Site

Material: IR, PVC
Color: Clear, Yellow
Sterilization: Eto
ID: 2.4 mm



IS-022

Injection Site

Material: PC, PIP
Color: Clear, Natural
Sterilization: Eto, Gamma, Steam
ID: 3.6 mm | OD: 7.1 mm



IS-025

Y-Injection Site

Material: ABS, PIP
Color: Clear, Natural
Sterilization: Eto, Gamma
ID: 4.1 mm



IS-025A

Y-Injection Site

Material: ABS, PIP
Color: Clear, Natural
Sterilization: Eto, Gamma
ID: 4.1 mm



IS-030

Y-Injection Site

Material: IR, PVC
Color: Clear, Yellow
Sterilization: Eto, Gamma
ID: 2.1 mm



IS-032

Y-Injection Site

Material: ABS, PIP
Color: Natural
Sterilization: Eto, Gamma
ID: 4.1 mm



IS-033

T-Injection Site with Rotating Male Luer Lock and 2 Ports


Material: ABS
Color: Clear, Yellow
Sterilization: Eto, Gamma
ID: 2.4 mm



IS-036DEHT

T-Injection Site

Material: PVC-DEHT, TPE
Color: Blue, Clear, Natural
Sterilization: Eto
ID: 6.8 mm



IS-039

T-Injection Site

Material: PP, PVC, TPE
Color: Clear, Natural, White
Sterilization: Eto
ID: 4.1 mm



IS-039G

T-Injection Site

Material: PP, PVC, TPE
Color: Clear, Natural, White
Sterilization: Eto, Gamma
ID: 4.1 mm



IS-041DEHT

T-Injection Site


Material: PP, PVC-DEHT, TPE
Color: Red, Clear, Natural
Sterilization: Eto, Gamma
ID: 6.8 mm



IS-045

T-Injection Site

Material: PIP, PP, PVC
Color: Clear, Natural, Red
Sterilization: Eto
ID: 5.5, 6.8 mm



IS-046

T-Injection Site with Female Port

Material: ABS
Color: Clear, Natural, White
Sterilization: Eto
ID: 5.5 mm



IS-047

Injection Site with Tear-Off Cap and Male Luer Slip

Material: PC, PIP, PP
Color: Clear, Natural
Sterilization: Eto, E-Beam, Gamma
ID: 3.0 mm | OD: 4.8 mm



IS-048

Injection Site with Tear-Off Cap

Material: PC, PIP, PP
Color: Clear, Natural
Sterilization: Eto, E-Beam, Gamma, Steam
ID: 5.0 mm | OD: 6.0 mm



IS-049

Injection Site with Tear-Off Cap and Male Luer Slip

Material: PIP, PP, PVC
Color: Clear, Natural
Sterilization: Eto
OD: 4.8 mm



IS-052

Injection Site

Material: PC, PP, TPE
Color: Clear, Natural
Sterilization: Eto, Gamma, Steam
OD: 5.0 mm



IS-052A

Injection Site

Material: PC, PP
Color: Clear, Natural
Sterilization: Eto
OD: 5.0 mm



IS-054

Injection Site

Material: PC, PP, TPE
Color: Clear, Natural
Sterilization: Eto, Steam
ID: 4.3 mm | OD: 6.8 mm



IS-054B

Injection Site

Material: PC, PP, TPE
Color: Blue, Clear, Natural
Sterilization: Eto, Gamma
OD: 6.0 mm



IS-054R

Injection Site

Material: PC, PP, TPE
Color: Clear, Natural, Red
Sterilization: Eto, Gamma, Steam
ID: 4.3 mm | OD: 6.0 mm



IS-055

Injection Site, Barbed

Material: PC, PP, TPE
Color: Blue, Clear, Natural
Sterilization: Eto, Gamma, Steam
ID: 3.0 mm | OD: 5.0 mm



IS-056

Injection Site

Material: BIIR, PC, PP
Color: Black, Clear, Natural
Sterilization: Eto, Steam
OD: 6.0 - 6.8 mm



IS-057

Injection Site

Material: PIP, PP, PVC
Color: Clear, Natural
Sterilization: Eto
OD: 5.0 mm



IS-058

Injection Site with Male Luer Lock

Material: ABS, TPE
Color: Clear, Natural
Sterilization: Eto, Gamma
ID: 2.6 mm



IS-060

Injection Site with Male Luer Slip

Material: ABS, TPE
Color: Clear, Natural
Sterilization: Eto, Gamma
ID: 2.5 mm | OD: 5.0 mm



IS-062

Injection Site


Material: ABS, TPE
Color: Clear, Yellow
Sterilization: Eto, Gamma
ID: 6.0 mm



IS-063

T-Injection Site with Cap

Material: IR, PP, PVC
DEHP Free
Color: Blue, Natural
Sterilization: Eto
OD: 6.8 mm



IS-065

Y-Injection Site, Pre-Slit

Material: ABS, PIP
Color: Clear, Natural
Sterilization: Eto, Gamma
ID: 4.1 mm



IS-069

T-Injection Site with Cap

Material: IR, PP, PVC
DEHP Free
Color: Natural, Red
Sterilization: Eto
OD: 6.8 mm



IS-092

Injection Spike

Material: Acrylic
Color: Clear
Sterilization: Eto, E-Beam, Gamma



IS-100

Injection Site

Material: PC, PIP, PVC
Color: Clear, Natural
Sterilization: Eto, Gamma, Steam
OD: 6.8 mm



IS-101

Injection Site

Material: PC, PIP
Color: Clear, Natural
Sterilization: Eto, Steam
OD: 6.0 mm



IS-108

Y-Injection Site, Pre-Slit

Material: IR, PC
Color: Clear, Natural
Sterilization: N/A



IS-109

Y-Injection Site with Male Luer Slip

Material: ABS, IR, PVC
Color: Clear, Yellow
Sterilization: Eto
ID: 2.4 mm | OD: 4.1 mm



IS-110

Y-Injection Site with Male Luer Slip

Material: ABS, IR, PVC
Color: Clear, Yellow
Sterilization: Eto
ID: 2.0 mm



IS-112

T-Injection Site with Vented Cap and Male Luer Lock

Material: ABS, HDPE, IR, PVC
Color: Clear, Natural, White
Sterilization: Eto
ID: 2.4 mm



IS-112S

T-Injection Site with Vented Cap and Male Luer Slip

Material: ABS, HDPE, IR, PVC
Color: Clear, Natural, White
Sterilization: Eto
ID: 2.4 mm



IS-113

T-Injection Site with 2 Ports and Male Luer Slip

Material: ABS, HDPE
Color: Clear, Yellow
Sterilization: Eto
ID: 2.2 mm



IS-128

Injection Site with Male Luer Lock

Material: IR, PC
Color: Clear, Natural
Sterilization: Eto



IS-130

Y-Injection Site

Material: ABS
Color: Clear, Natural
Sterilization: Eto



IS-131

Y-Injection Site

Material: IR, PVC
Color: Clear, Natural
Sterilization: Eto
ID: 2.4 mm



IS-149

NEW

Standard Injection Site with Male Luer Slip

Material: ABS, TPE
Color: Clear, Natural
Sterilization: Eto, Gamma



IS-151

Injection Site with Male Luer Slip

Material: PC, PIP, PP
Color: Clear, Natural
Sterilization: Eto, Steam
ID: 2.7 mm | OD: 4.1 mm



IS-154C

Injection Site

Material: IR, PC, PP
Color: Clear, Natural
Sterilization: Eto, Gamma, Steam
ID: 6.0 mm



IS-211

Injection Stopper / Heparin Cap

Material: ABS, PIP
Color: Yellow
Sterilization: Gamma



IS-005

Negative Needleless
Injection Site with
Male Luer Lock

Material: COPE, PC,
Silicone
Color: Blue, Clear
Sterilization: Eto
ID: 2.5 mm



IS-007A

Negative Needleless
Y Injection Site

Material: ABS, PC, Silicone
Color: Blue, Clear
Sterilization: Eto, Gamma
ID: 4.1 mm



IS-008A

Negative Needleless
Injection Site

Material: ABS, PC, Silicone
Color: Blue, Clear
Sterilization: Eto, Gamma
ID: 4.1 mm



IS-009A

Negative Needleless
Injection Site with
Male Luer Lock

Material: PC, Silicone
Color: Blue, Clear
Sterilization: Eto, Gamma
ID: 4.1 mm



IS-015

Negative Needleless
Injection Site with
Luer Slip

Material: PC, Silicone
Color: Blue, Clear
Sterilization: Eto, Gamma
ID: 2.1 mm
Pressure: 30 psi



IS-016

Negative Needleless
Injection Site

Material: PC, Silicone
Color: Blue, Clear
Sterilization: Eto, Gamma
ID: 4.1 mm
Pressure: 30 psi



IS-067

Negative Needleless
Injection Site

Material: PC, Silicone
Color: Natural
Sterilization: Eto, Gamma
OD: 5.5 mm
Pressure: 30 psi



IS-068

Negative Needleless
T Injection Site

Material: PC, Silicone
Color: Clear, Natural, Yellow
Sterilization: Eto, E-Beam,
Gamma
ID: 6.8 mm
Flow Rate: 250 mL/min



IS-071

Negative Needleless
Injection Site with
Male Luer Slip

Material: COPE, PC,
Silicone
Color: Clear
Sterilization: Eto
ID: 2.5 mm



IS-079

Negative Needleless
Injection Site with Male
Luer Lock, Photoprotective

Material: COPE, PC,
Silicone
Color: Brown
Sterilization: Eto
ID: 2.5 mm



IS-080

Negative Needleless
Injection Site,
Photoprotective

Material: COPE, Silicone
Color: Brown
Sterilization: Eto
ID: 2.3 mm
Priming Volume: 0.11 mL



IS-081

Negative Needleless
Y Injection Site,
Photoprotective

Material: COPE, PC,
Silicone
Color: Green
Sterilization: Eto
ID: 4.1 mm



IS-082

Negative Needleless
Injection Site with
Male Luer Lock

Material: COPE, PC,
Silicone
Color: Clear
Sterilization: Eto, Gamma
ID: 2.5 mm



IS-083

Negative Needleless
Injection Site

Material: COPE, Silicone
Color: Clear
Sterilization: Eto
ID: 4.1 mm



IS-094

Negative Needleless
Injection Site

Material: COPE, Silicone
Color: Clear
Sterilization: Eto
ID: 2.5 mm



IS-095

Negative Needleless
Y Injection Site

Material: COPE, PC,
Silicone
Color: Clear
Sterilization: Eto
ID: 4.1 mm



IS-096

Negative Needleless
Y Injection Site

Material: PC, Silicone
Color: Clear, Tinted
Sterilization: Eto, Gamma
ID: 4.1 mm



IS-097

Negative Needleless
Injection Site with Rotating
Male Luer Lock

Material: PC, Silicone
Color: Clear, Purple
Sterilization: Eto, Gamma
Pressure: 325 psi



IS-133

Negative Needleless
T Injection Site

Material: PC, PVC
DEHP Free
Color: Blue, Clear
Sterilization: Eto
ID: 6.8 mm



IS-138

Negative Needleless
T Injection Site

Material: PC, Silicone
Color: Clear, Natural, Red
Sterilization: Eto, E-Beam,
Gamma
ID: 6.8 mm
Flow Rate: 250 mL/min



IS-139

Negative Needleless
T Injection Site

Material: PC, Silicone
Color: Blue, Clear, Natural
Sterilization: Eto, E-Beam,
Gamma
ID: 6.8 mm
Flow Rate: 250 mL/min



IS-141

NEW

Negative Needleless
Injection Site, Swabbable
Valve

Material: PC, Silicone
Color: Blue, Tinted
Sterilization: Eto, Gamma
Pressure: 30 psi (Backflow)
Priming Volume: <0.2 mL



IS-143

Negative Needleless
T Injection Site

Material: PC, Silicone
Color: Blue, Clear
Sterilization: Eto, Gamma
ID: 6.8 mm
Pressure: 30 psi



IS-146

NEW

Negative Needleless
Injection Site with Male
Luer for CT Syringes

Material: PC, Silicone
Color: Natural, Tinted
Sterilization: N/A



IS-147

NEW

Negative Needleless
Injection Site for
CT Syringes

Material: PC, Silicone
Color: Natural, Tinted
Sterilization: N/A



IS-148

NEW

Negative Needleless
Injection Site, One-Way

Material: MABS, PC,
Silicone
Color: Blue, Natural,
Tinted, White
Sterilization: N/A
ID: 2.3 mm



IS-245

Negative Needleless
Y Injection Site

Material: COPE, Silicone
Color: Blue, Clear
Sterilization: Eto, Gamma
ID: 4.1 mm | OD: 6.8 mm
Pressure: 30 psi



CA-157

Cap, Olav, for CV-157

Material: TPE
Color: White
Sterilization: N/A



CV-157

Neutral Needleless
Injection Site with
Female Luer Lock, Olav

Material: PC, PSU, Silicone,
TPE
Color: Clear, Orange
Sterilization: Eto, Gamma



IS-026

Neutral Needleless
Injection Site, Closed
Connector, Female to Male

Material: ABS, PC, Silicone
Color: Blue, Clear
Sterilization: Eto, Gamma
ID: 4.1 mm



IS-028

Neutral Needleless
Injection Site

Material: ABS, PP, Silicone
Color: Blue, Clear
Sterilization: Eto
ID: 2.3 mm



IS-066

Neutral Needleless
Injection Site, Closed,
with Male Luer Lock

Material: HDPE, MABS,
PC, PP, Silicone
Color: Clear, Natural,
Yellow
Sterilization: Eto



IS-104

Neutral Needleless
Injection Site with
Male Luer Lock

Material: PC, Silicone
Color: Clear, Purple
Sterilization: Eto



IS-105

Neutral Needleless
Injection Site with
Male Luer Lock, Pediatric

Material: PC, Silicone
Color: Clear, Purple
Sterilization: Eto



IS-106

Neutral Needleless
Injection Site

Material: PC, Silicone
Color: Clear, Purple
Sterilization: Eto
ID: 4.1 mm



IS-107

Neutral Needleless
Y-Injection Site

Material: PC, Silicone
Color: Clear, Purple
Sterilization: Eto
ID: 4.1 mm



IS-118

Neutral Needleless
T-Injection Site with
Male Luer Slip

Material: ABS, PC, Silicone
Color: Clear, Natural
Sterilization: Eto, Gamma,
Steam


ID: 2.3 mm



IS-119

Neutral Needleless
Injection Site with
Male Luer Lock

Material: ABS, PC, Silicone
Color: Clear, Natural
Sterilization: Eto, Gamma,
Steam



IS-120

Neutral Needleless
Y-Injection Site

Material: ABS, PC, Silicone
Color: Clear, Natural
Sterilization: Eto, Gamma,
Steam

ID: 2.3 mm



IS-121

Neutral Needleless
Injection Site

Material: ABS, PC, Silicone
Color: Clear, Natural
Sterilization: Eto, Gamma,
Steam

ID: 4.1 mm



IS-122

Neutral Needleless
Y-Injection Site

Material: ABS, PC, Silicone
Color: Clear, Natural
Sterilization: Eto, Gamma,
Steam

ID: 3.8 mm



IS-124

Neutral Needleless
Y-Injection Site

Material: ABS, PC, Silicone
Color: Natural, Tinted
Sterilization: Eto, Gamma,
Steam

ID: 4.1 mm



IS-125

Neutral Needleless
Injection Site with
Male Luer Lock

Material: ABS, PC, Silicone
Color: Clear
Sterilization: Eto
Pressure: 350 psi



IS-126

Neutral Needleless
Y-Injection Site

Material: ABS, PP, Silicone
Color: Clear
Sterilization: Eto
ID: 4.1 mm
Pressure: 350 psi



IS-127

Neutral Needleless
Injection Site

Material: ABS, PP, Silicone
Color: Clear
Sterilization: Eto
ID: 2.9 mm
Pressure: 350 psi



IS-129

Neutral Needleless
Injection Site with Male
Luer Lock and Bag Port

Material: HDPE, MABS, PC,
PVC DEHP Free,
Silicone
Color: Clear, White, Yellow
Sterilization: Eto



IS-132

Neutral Needleless
Injection Site

Material: MABS, PC,
Silicone
Color: Clear, Natural
Sterilization: Eto
ID: 2.3 mm



IS-134

Neutral Needleless
Y Injection Site

Material: ABS, PP, Silicone
Color: Blue, Clear
Sterilization: Eto



IS-136

Neutral Needleless Injection
Site with Male Luer Lock

Material: COPE, MABS,
Silicone
Color: Clear, Light Blue,
Natural
Sterilization: Eto, E-Beam,
Gamma, Steam
Volume: 0.07 mL
Pressure: 29 psi



IS-137

Neutral Needleless
Y-Injection Site

Material: COPE, MABS,
Silicone
Color: Clear, Light Blue,
Natural
Sterilization: Eto, E-Beam,
Gamma, Steam
ID: 4.0 mm
Volume: 0.11 mL
Pressure: 29 psi



IS-150B

NEW

Neutral Needleless
Injection Site

Material: PC, PP, Silicone
Color: Blue, Natural,
Tinted
Sterilization: Eto



Neutral Needleless Injection Sites

67

IS-150G

NEW

Neutral Needleless
Injection Site

Material: PC, PP, Silicone
Color: Green, Natural,
Tinted
Sterilization: Eto



IS-150OR

NEW

Neutral Needleless
Injection Site

Material: PC, PP, Silicone
Color: Natural, Orange,
Tinted
Sterilization: Eto



IS-150R

NEW

Neutral Needleless
Injection Site

Material: PC, PP, Silicone
Color: Natural, Red, Tinted
Sterilization: Eto



IS-150W

NEW

Neutral Needleless
Injection Site

Material: PC, PP, Silicone
Color: Natural, Tinted
White
Sterilization: Eto



DirectMed

+1 (516) 656-3377 • sales@directmed.com
www.directmed.com

DC-013K

Non-Vented IV Spike for
Drip Chambers

Material: PVC
Color: White
Sterilization: Eto
ID: 5.5 mm
Drops: 20



RW-026

Non-Vented IV Spike with
Snap Cap and Filter

Material: ABS, PTFE
Color: Green, White
Sterilization: Eto
Microns: 0.1 µm



RW-027

Non-Vented IV Spike with
Snap Cap and Filter

Material: ABS, PET, PTFE
Color: Blue, White
Sterilization: Eto
Microns: 0.1 µm



SP-002

Non-Vented Micro IV Spike

Material: PVC
Color: White
Sterilization: Eto



SP-009

Non-Vented IV Spike,
One-Way, with Vented Cap

Material: ABS, LDPE
Color: Natural, White
Sterilization: Eto, E-Beam,
Gamma
ID: 4.1 mm



SP-011

Non-Vented IV Spike for
Drip Chambers

Material: ABS
Color: White
Sterilization: Eto, Gamma
Drops: 20



SP-012

Non-Vented IV Spike

Material: ABS
Color: White
Sterilization: Eto, Gamma
ID: 5.0 mm



SP-013

Non-Vented IV Spike,
One-Way

Material: ABS
Color: White
Sterilization: Eto, Gamma
ID: 5.5 mm



SP-014

Non-Vented IV Spike

Material: ABS
Color: White
Sterilization: Eto, Gamma
ID: 4.1 mm



SP-014CC

Non-Vented IV Spike
with Vented Cap

Material: ABS, LDPE
Color: Clear, White
Sterilization: Eto, Gamma
ID: 4.1 mm



SP-014CLEAR

Non-Vented IV Spike
with Vented Cap

Material: LDPE, PVC
DEHP Free
Color: Clear
Sterilization: Eto
ID: 4.1 mm



SP-015

Non-Vented IV Spike
with Female Luer

Material: ABS
Color: White
Sterilization: Eto, Gamma



SP-015A

Non-Vented IV Spike with Female Luer Lock

Material: ABS
Color: White
Sterilization: Eto, Gamma



SP-015WC

Non-Vented IV Spike with Female Luer Lock and Cap

Material: ABS, LDPE
Color: White
Sterilization: Eto, Gamma



SP-020

Non-Vented IV Spike

Material: ABS
Color: White
Sterilization: Eto, Gamma
ID: 5.5 mm



SP-021

Non-Vented IV Spike

Material: ABS
Color: White
Sterilization: Eto, Gamma
ID: 6.8 mm



SP-023

Non-Vented IV Spike, Compounder

Material: ABS
Color: White
Sterilization: Eto, Gamma
ID: 6.8 mm



SP-025

Non-Vented Chemo IV Spike with Snap Cap, Air and Solution Filters, and Female Luer Lock

Material: MABS, PTFE
Color: Clear, Red
Sterilization: Eto
Microns: 0.1 µm (Air Filter)
5 µm (Sol. Filter)
Pressure: 1 bar



SP-029

Non-Vented IV Spike for Drip Chambers

Material: ABS
Color: White
Sterilization: Eto
Drops: 20



SP-033

Non-Vented Mini IV Spike

Material: ABS
Color: White
Sterilization: Eto
ID: 5.0 mm | OD: 6.0 mm



SP-033C

Spike Cap for SP-033

Material: PP
Color: Natural
Sterilization: Eto
ID: 6.5 mm



SP-034

Non-Vented IV Spike with Cap

Material: ABS, LDPE
Color: White
Sterilization: Eto, Gamma
ID: 4.1 mm



SP-035

Non-Vented IV Spike with Cap and Female Luer Port

Material: ABS, LDPE
Color: Clear, White
Sterilization: Eto, Gamma



SP-040

Non-Vented Micro IV Spike with Cap and Female Luer Lock

Material: PC, PP
Color: Clear, Purple
Sterilization: Eto



SP-047

Non-Vented IV Spike,
One-Way

Material: ABS
Color: White
Sterilization: Eto, Gamma



SP-050

Non-Vented IV Spike with
Cap, Injection Site, and
Bag Port

Material: EVA, MABS, PC,
PE, PVC, Silicone
Color: Clear, Natural,
Tinted, White
Sterilization: Eto



SP-052

Non-Vented IV Spike,
One-Way, with Vented Cap

Material: ABS, LDPE
Color: White
Sterilization: Eto
ID: 4.1 mm | OD: 4.8 mm



SP-057

Non-Vented IV Spike,
One-Way with Cap

Material: ABS, PE
Color: Natural, White
Sterilization: Eto
ID: 4.1 mm | OD: 4.8 mm



SP-060

NEW

Non-Vented IV Spike, Dual
Flow, with Needleless Site

Material: MABS, PC, HDPE,
PP, Silicone
Color: Natural, Tinted,
White, Yellow
Sterilization: Eto
ID: 3.88 mm



SP-061

NEW

Non-Vented IV Spike, Dual
Flow, with Needleless Site

Material: MABS, PC, HDPE,
PP, Silicone
Color: Natural, Tinted,
White, Yellow
Sterilization: Eto
ID: 3.44 mm



SP-079

Non-Vented IV Spike for
Drip Chambers

Material: ABS
Color: White
Sterilization: Eto, Gamma
Drops: 10



SP-087

Non-Vented IV Spike with
Cap, for Drip Chambers

Material: ABS, PE
Color: Clear, Natural
Sterilization: N/A



SP-230

Non-Vented IV Spike with
Barbed Connector

Material: ABS
Color: White
Sterilization: Eto, Gamma,
Steam
OD: 1/4"



RW-011

Vented IV Spike with Snap Cap and Filter

Material: ABS, PET, PTFE
Color: Blue, White
Sterilization: Eto
Microns: 0.1 µm



RW-012

Vented IV Spike with Snap Cap, Air and Solution Filters

Material: MABS, PET, PTFE
Color: Clear, Red
Sterilization: Eto
Microns: 0.1 µm (Air Filter)
5 µm (Sol. Filter)



RW-014

Vented IV Spike with Snap Cap and Filter

Material: ABS, PTFE
Color: Green, White
Sterilization: Eto, Gamma
ID: 4.1 mm
Microns: 0.1 µm



SP-003

Vented IV Spike with Cap and Filter

Material: ABS, Acrylic, HDPE, LDPE
Color: White, Natural
Sterilization: Eto
ID: 4.1 mm | OD: 4.8 mm
Microns: 3 µm



SP-006

Vented IV Spike with Cap and Filter

Material: ABS, Acrylic, LDPE, PE
Color: Natural, Red, White
Sterilization: Eto
ID: 4.1 mm
Microns: 3 µm



SP-016

Vented IV Spike with Cap and Filter

Material: ABS, LDPE, PP, PVDF
Color: White
Sterilization: Eto, Gamma
ID: 4.1 mm | OD: 4.8 mm
Microns: 3 µm



SP-017

Vented IV Spike with Cap, Filter, and Female Luer Lock

Material: ABS, LDPE, PP, PVDF
Color: White
Sterilization: Eto, Gamma
Microns: 3 µm



SP-017A

Vented IV Spike with Female Luer Lock

Material: ABS
Color: White
Sterilization: Eto, Gamma



SP-024

Vented Mini IV Spike with Female Luer Lock

Material: ABS
Color: Natural, White
Sterilization: Eto



SP-024C

Spike Cap for SP-024

Material: HDPE
Color: Natural
Sterilization: Eto



SP-031

Vented IV Spike with Snap Cap and Filter

Material: MABS, PP, PTFE
Color: Clear, Green
Sterilization: Eto
Microns: 0.1 µm



SP-032

Vented Mini Transfer Spike with Cap and Female Luer

Material: ABS, ePTFE, HDPE, PET
Color: Natural, White
Sterilization: Eto
ID: 4.1 mm
Microns: 0.2 µm



SP-038

Vented Mini IV Spike with Cap and Filter

Material: ABS, LDPE, PP
Color: Natural, White
Sterilization: Eto
ID: 3.6 mm
Microns: 0.45 µm



SP-044

Vented IV Spike with Cap

Material: ABS, ePTFE, LDPE, PET, PVC DEHP Free
Color: Blue, Natural, White
Sterilization: Eto
ID: 4.0 mm
Microns: 1 µm



SP-045

Vented IV Spike with Cap and Filter

Material: ABS, LDPE, PP, PVDF
Color: Clear, White
Sterilization: Eto
ID: 4.1 mm
Microns: 1.2 µm



SP-048

Vented IV Spike with Cap and Filter

Material: ABS, HDPE, LDPE, PE, PP, PVDF
Color: Natural, White
Sterilization: Eto
ID: 4.1 mm | OD: 4.8 mm
Microns: 0.8 µm



SP-049

Vented IV Spike with Cap and Needleless Injection Site

Material: Acrylic, EVA, MABS, Nylon, PC, PE, PP, Silicone
Color: Clear, Natural, White
Sterilization: Eto
Microns: 3 µm



SP-056

Vented IV Spike, Hyperal Piercing

Material: ABS, LDPE, PET, PVDF, Silicone, SS
Color: Natural, White
Sterilization: Eto
OD: 4.3 mm
Microns: 3 µm



SP-058

Vented IV Spike with Cap and Air Filter

Material: ABS, Acrylic, HDPE, LDPE
Color: Natural, White
Sterilization: Eto
ID: 4.1 mm
Microns: 1.2 µm



SP-059

NEW

Vented IV Spike with Cap, Female Luer Lock, and Air Filter

Material: ABS, HDPE, LDPE, PP, PVDF
Color: Blue, Natural, White
Sterilization: Eto, E-Beam, Gamma
Microns: 3 µm



SP-039

Vial Adapter with Cap, Air Filter, and Female Luer Lock

Material: ABS, ePTFE, HDPE, PET
Color: Natural, White
Sterilization: Eto
Microns: 0.2 µm



SP-128

Vial Adapter with Cap and Dual Filters

Material: ABS, Acrylic, LDPE, Nylon
Color: Clear, White
Sterilization: Eto
Microns: 0.22 µm



VA-01313

Vial Adapter

Material: PC
Color: Tinted
Sterilization: Eto, Gamma
ID: 13 mm to 13 mm



VA-01320

Vial Adapter

Material: PC
Color: Tinted
Sterilization: Eto, Gamma
ID: 13 mm to 20 mm



VA-013F

Vial Adapter with Female Luer

Material: PC
Color: Clear, Purple
Sterilization: Eto, Gamma
ID: 13 mm



VA-013IS

Vial Adapter, Vented, to Injection Site

Material: ePTFE, MABS, PC, PP, Silicone
Color: Natural, Tinted, White
Sterilization: Eto
ID: 13 mm



VA-014

Vial Adapter, Spike with Cap to Injection Site

Material: HDPE, MABS, PC, PP, Silicone
Color: Natural, Tinted
Sterilization: Eto
Dimensions: 12 mm



VA-015

Vial Adapter, Vented Spike with Cap to Injection Site

Material: ePTFE, HDPE, MABS, PC, PP, Silicone
Color: Natural, Tinted, White
Sterilization: Eto
ID: 20 mm to 13 mm



VA-02020

Vial Adapter

Material: PC
Color: Tinted
Sterilization: Eto, Gamma
ID: 20 mm to 20 mm



VA-020F

Vial Adapter with Female Luer

Material: PC
Color: Clear, Purple
Sterilization: Eto, Gamma
ID: 20 mm



VA-020IS

Vial Adapter, Vented, to Injection Site

Material: ePTFE, MABS, PC, PP, Silicone
Color: Natural, Tinted, White
Sterilization: Eto
ID: 20 mm



BH-011

NEW

Bag Handle with Bed Hook

Material: PP
Color: Blue
Sterilization: Eto



BH-013

Bag Handle

Material: PP
Color: Blue
Sterilization: Eto, E-Beam, Gamma



BH-015

NEW

Bag Hook

Material: HDPE
Color: Blue
Sterilization: Eto, Gamma



BH-016

Bag Hanger with Pins

Material: PA
Color: Blue
Sterilization: N/A



BH-017

Hook for BH-016

Material: PA
Color: Blue
Sterilization: N/A



FH-001

Forceps, Wide Grip

Material: Nylon
Color: Blue
Sterilization: N/A



FH-002

Forceps, Wide Grip

Material: Nylon
Color: Red
Sterilization: N/A



FH-009

Forceps

Material: POM
Color: Blue
Sterilization: Eto, Steam



GS-011

Silicone Gasket

Material: Silicone
Color: Natural
Sterilization: N/A
ID: 3.0 mm



GW-010

Guide Wire Holder

Material: Brass, PC
Color: Clear
Sterilization: Eto, Steam



GW-011

Guide Wire Holder

Material: HDPE, PP
Color: Blue, White
Sterilization: Eto, Steam
ID: 0.046"



GW-012

Guide Wire Holder

Material: ABS, Brass, HDPE
Color: White
Sterilization: N/A
ID: 0.96 mm



GW-013

Torquer

Material: PC
Color: Natural
Sterilization: Eto, Steam
ID: 0.50 mm



GW-014

Guide Wire Holder

Material: ABS, Brass, HDPE
Color: Blue, White
Sterilization: N/A
ID: 0.96 mm



GW-014B

Guide Wire Holder Body

Material: HDPE
Color: White
Sterilization: N/A



GW-014C

Guide Wire Holder Cap

Material: ABS
Color: Blue
Sterilization: N/A



GW-015

Guide Wire Holder

Material: PP
Color: Natural
Sterilization: Eto, Steam
ID: 0.046"



GW-017

Guide Wire Holder Cap

Material: PC
Color: Clear
Sterilization: N/A
ID: 0.99 mm



GW-018

Guide Wire Holder Cap

Material: Brass, PC
Color: Clear
Sterilization: Eto, Steam



SC-009

Scalpel with #10 Blade

Material: PS, Stainless Steel
Color: Green, Natural
Sterilization: Gamma



SC-011

Scalpel with #11 Blade

Material: PS, Stainless Steel
Color: Blue, Natural
Sterilization: Gamma



SC-015

Scalpel with #15 Blade

Material: PS, Stainless Steel
Color: Natural, Yellow
Sterilization: Gamma



LM-102

Male Luer, Microbore

Material: PC
Color: Clear
Sterilization: Eto, Steam
ID: 0.029"



LM-052

Male Luer Lock, Microbore

Material: PC
Color: Clear
Sterilization: Eto, Steam
ID: 0.031"



LM-054

Male Luer Lock, Microbore

Material: PC
Color: Clear
Sterilization: Eto, Steam
ID: 0.034"



LM-055

Male Luer Lock, Microbore

Material: PC
Color: Clear
Sterilization: Eto, Steam
ID: 0.056"



LM-029

Male Luer

Material: ABS
Color: Clear
Sterilization: Eto, Gamma
ID: 2.0 mm



LM-029A

Male Luer

Material: ABS
Color: Clear
Sterilization: Eto, Gamma
ID: 2.0 mm



LM-141

Male Luer, Rotating, with Filtered Cap

Material: ABS, Acrylic, LDPE, PP
Color: Clear, Natural, White
Sterilization: Eto
ID: 2.0mm | Microns: 3µm



LM-141B

Male Luer, Rotating, with Filtered Cap

Material: ABS, Acrylic, LDPE, PP
Color: Blue, Clear, Natural, White
Sterilization: Eto
ID: 2.0mm | Microns: 3µm



LM-145

Male Luer, Rotating, with Vented Cap and Filter

Material: MABS, PE, PP, PVDF
Color: Clear, Natural
Sterilization: Eto
ID: 2.0 mm
Microns: 0.8 µm



LM-145B

Male Luer, Rotating, with Vented Cap and Filter

Material: MABS, PE, PP, PVDF
Color: Blue, Clear, Natural
Sterilization: Eto
ID: 2.0 mm
Microns: 0.8 µm



LM-146

Male Luer Lock, Rotating, with Vented Cap and Filter

Material: ABS, PP, PTFE
Color: Blue, Clear
Sterilization: Eto
ID: 2.0 mm



LM-148B

NEW

Male Luer Lock, Rotating

Material: MABS
Color: Blue, Clear, Natural
Sterilization: Eto, E-Beam, Gamma
ID: 2.0 mm



LM-017

Male Luer Lock

Material: PVC
Color: Light Blue
Sterilization: Eto
ID: 2.1 mm



LM-088GS

Male Luer

Material: PVC
Color: Clear
Sterilization: Eto, Gamma
ID: 2.1 mm



LM-022

Male Luer Lock

Material: ABS
Color: Clear
Sterilization: Eto, Gamma
ID: 2.2 mm | OD: 4.3 mm



LM-121

Male Luer, Rotating, with Slide-Off Hub and Cap

Material: ABS, PC
Color: Clear
Sterilization: Eto, Gamma
ID: 2.2 mm



LM-035

Male Luer, Rotating

Material: ABS, PVC
Color: Clear
Sterilization: Eto
ID: 2.4 mm



LM-089

Male Luer

Material: PVC
Color: Clear
Sterilization: Eto
ID: 2.4 mm



LM-089GS

Male Luer Lock

Material: PVC
Color: Blue
Sterilization: Eto, Gamma
ID: 2.4 mm



LM-108

Male Luer, Rotating

Material: ABS, PVC
Color: Clear
Sterilization: Eto
ID: 2.4 mm



LM-117

Male Luer

Material: PVC
Color: Light Blue
Sterilization: Eto
ID: 2.4 mm



LM-120

Male Luer, Rotating

Material: ABS, PVC
Color: Clear
Sterilization: Eto
ID: 2.4 mm



LM-151

NEW

Male Luer Lock

Material: PVC
Color: Clear
Sterilization: Eto
ID: 2.4 mm



LM-152

NEW

Male Luer, Rotating Shaft

Material: PC, Silicone
Color: Clear, Red
Sterilization: N/A
ID: 2.4 mm



LM-028

Male Luer

Material: MABS
Color: Clear
Sterilization: Eto, Gamma
ID: 2.5 mm



LM-028A

Male Luer

Material: ABS
Color: Clear
Sterilization: Eto, Gamma
ID: 2.5 mm



LM-031ABS

Male Luer Lock, Rotating

Material: ABS, PP
Color: Clear
Sterilization: Eto, Gamma
ID: 2.5 mm



LM-031CAPDEHT

Male Luer, Rotating, with Cover

Material: LDPE, PP, PVC-DEHT
Color: Clear
Sterilization: Eto
ID: 2.5 mm



LM-031DEHT

Male Luer, Rotating

Material: PP, PVC-DEHT
Color: Clear
Sterilization: Eto
ID: 2.5 mm



LM-104

Male Luer, Photoprotective

Material: PC
Color: Brown
Sterilization: Eto, Gamma
ID: 2.5 mm



LM-124

Male Luer, Rotating, with Slide-Off Hub and Cap

Material: ABS, PP
Color: Clear
Sterilization: Eto
ID: 2.5 mm



LM-064

Male Luer, Rotating

Material: ABS
Color: Clear
Sterilization: Eto
ID: 2.79 mm



LM-114

Male Luer

Material: PVC
Color: Clear
Sterilization: Eto
ID: 2.8 mm



LM-114GR

Male Luer

Material: PVC
Color: Green
Sterilization: Eto
ID: 2.8 mm



LM-114GY

Male Luer

Material: PVC
Color: Grey
Sterilization: Eto
ID: 2.8 mm



LM-027

Male Luer

Material: ABS
Color: Clear
Sterilization: Eto, Gamma
ID: 3.0 mm



LM-027A

Male Luer

Material: ABS
Color: Clear
Sterilization: Eto, Gamma
ID: 3.0 mm



LM-033

Male Luer, Rotating

Material: ABS, PP
Color: Clear
Sterilization: Eto, Gamma
ID: 3.0 mm



LM-090

Male Luer

Material: ABS
Color: Clear
Sterilization: Eto, Gamma
ID: 3.0 mm



LM-122

Male Luer Lock

Material: ABS
Color: Clear
Sterilization: Eto
ID: 3.0 mm



LM-131

Male Luer Lock, Rotating, with Cover

Material: ABS, LDPE, PVC
DEHP Free
Color: Clear, Natural
Sterilization: Eto, Gamma
ID: 3.0 mm



LM-147

Male Luer Lock, Rotating, with Vented Cap and Filter

Material: ABS, HDPE, PTFE
Color: Clear, White
Sterilization: Eto
OD: 3.0 mm



LM-119

Male Luer Lock

Material: PVC
Color: Clear
Sterilization: Eto
ID: 3.2 mm



LM-015

Male Luer Lock

Material: ABS
Color: Clear
Sterilization: Eto, Gamma
ID: 3.3 mm



LM-106

Male Luer, Photoprotective

Material: PC
Color: Brown
Sterilization: Eto, Gamma
ID: 3.3 mm



LM-137

Male Luer Lock

Material: PC
Color: Clear
Sterilization: Eto, Steam
ID: 3.3 mm



LM-010

Male Luer, Rotating

Material: ABS, PVC
Color: Clear
Sterilization: N/A
ID: 3.5 mm



LM-093

Male Luer

Material: ABS
Color: Clear
Sterilization: Eto, Gamma
ID: 3.5 mm



LM-115

Male Luer

Material: PVC
Color: Clear
Sterilization: Eto
OD: 3.5 mm



LM-118

Male Luer

Material: PVC
Color: Light Blue
Sterilization: Eto
ID: 3.5 mm



LM-111

Male Luer, Rotating, with Cap

Material: ABS, Acrylic, PP
Color: Blue, Clear
Sterilization: Eto
ID: 3.6 mm
Microns: 2 µm



LM-134

Male Luer, Rotating, with Vented Cap

Material: ABS, PP
Color: Clear
Sterilization: Eto
OD: 4.0 mm



LM-014

Male Luer Lock

Material: PC
Color: Clear
Sterilization: Eto, Gamma, Steam
ID: 4.1 mm



LM-014-CAP-A

Cap for LM-014

Material: PP
Color: Natural
Sterilization: Eto



LM-014-CAP-B

Cap for LM-014

Material: PP
Color: White
Sterilization: Eto



LM-016

Male Luer Slip with Rotating Hub

Material: ABS
Color: Clear
Sterilization: N/A
ID: 4.1 mm



LM-026

Male Luer

Material: ABS
Color: Clear
Sterilization: Eto, Gamma
ID: 4.1 mm



LM-032

Male Luer Lock, Rotating

Material: ABS, PP
Color: Clear
Sterilization: Eto, Gamma
ID: 4.1 mm



LM-032A

Male Luer Lock, Rotating

Material: ABS, PP
Color: Blue, Clear
Sterilization: Eto, Gamma
ID: 4.1 mm



LM-034

Male Luer, Rotating

Material: ABS, PP
Color: Clear
Sterilization: Eto, Gamma
ID: 4.1 mm



LM-034-11

Male Luer, Rotating,
Self-Ejecting

Material: PC, PVC
DEHP Free

Color: Clear
Sterilization: Eto, E-Beam,
Gamma

ID: 4.1 mm



LM-034B

Male Luer, Rotating

Material: ABS, PP
Color: Clear
Sterilization: Eto
ID: 4.1 mm



LM-034BODY

Male Luer Stem, Rotating

Material: ABS
Color: Clear
Sterilization: Eto, Gamma
ID: 4.1 mm | OD: 6.8 mm



LM-041

Male Luer Lock, Rotating,
Self-Ejecting

Material: PC, PVC
Color: Clear
Sterilization: Eto

ID: 4.1 mm



LM-059

Male Luer Lock, Rotating,
with Cover

Material: ABS, PE, PP
Color: Clear
Sterilization: Eto
ID: 4.1 mm



LM-068

Male Luer Lock, Rotating

Material: ABS, PP
Color: Natural
Sterilization: Eto, Gamma
ID: 4.1 mm



LM-091

Male Luer

Material: PC
Color: Clear
Sterilization: Eto
ID: 4.1 mm



LM-098

Male Luer Lock, Rotating,
with Cover

Material: PC, PVC
Color: Clear
Sterilization: Eto
ID: 4.1 mm



LM-099

Male Luer, Rotating

Material: PP, PVC
Color: Clear
Sterilization: Eto
ID: 4.1 mm



LM-099ABS

Male Luer, Rotating

Material: ABS
Color: Clear
Sterilization: Eto
ID: 4.1 mm



LM-105

Male Luer, Photoprotective

Material: PC
Color: Brown
Sterilization: Eto, Gamma
ID: 4.1 mm



LM-110

Male Luer, Rotating,
with Cap

Material: ABS, PP
Color: Clear
Sterilization: Eto
ID: 4.1 mm



LM-125

Male Luer, Rotating, with
Slide-Off Hub, Cap, and
Filter

Material: ABS, PP, PTFE
Color: Blue, Clear
Sterilization: Eto
ID: 4.1 mm
Microns: 1.2 µm



LM-129

Male Luer Lock, Rotating,
with Cover

Material: MABS, PE, PP
Color: Clear, Natural
Sterilization: Eto
ID: 4.1 mm



LM-129NC

Male Luer Lock, Rotating

Material: MABS, PP
Color: Clear
Sterilization: Eto
ID: 4.1 mm



Male Luers - 4.1 to 5.5 mm ID

81

LM-129PVC

Male Luer Lock, Rotating

Material: PP, PVC
DEHP Free
Color: Clear
Sterilization: Eto
ID: 4.1 mm



LM-138

Male Luer Lock

Material: ABS
Color: Clear
Sterilization: Eto
OD: 4.1 mm



LM-139

NEW

Male Luer Lock with Cap

Material: MABS, PP
Color: Blue, Clear
Sterilization: Eto
ID: 4.1 mm



LM-140

Male Luer, Rotating, with
Filtered Cap

Material: ABS, Acrylic,
LDPE, PP
Color: Clear, Natural, White
Sterilization: Eto
ID: 4.1 mm
Microns: 3 µm



LM-140B

Male Luer, Rotating, with
Filtered Cap

Material: ABS, Acrylic,
LDPE, PP
Color: Blue, Clear, Natural,
White
Sterilization: Eto
ID: 4.1mm | Microns: 3µm



LM-142

Male Luer Lock, Rotating

Material: MABS, PVC
Color: Clear
Sterilization: Eto
ID: 4.1 mm



LM-144

Male Luer, Rotating, with
Vented Cap

Material: MABS, PE, PP,
PVDF
Color: Clear, Natural
Sterilization: Eto
ID: 4.1 mm
Microns: 0.8 µm



LM-144B

Male Luer, Rotating, with
Vented Cap

Material: MABS, PE, PP,
PVDF
Color: Blue, Clear, Natural
Sterilization: Eto
ID: 4.1 mm
Microns: 0.8 µm



LM-150

NEW

Male Luer, Rotating Shaft

Material: PC, Silicone
Color: Clear, Red
Sterilization: N/A
ID: 4.1 mm



LM-154

Male Luer

Material: PVC DEHP Free
Color: Natural
Sterilization: Eto
ID: 4.1 mm
Pressure: 4 bar



LM-123

Male Luer Lock

Material: PVC
Color: Clear
Sterilization: Eto
ID: 4.4 mm



LM-039

Male Luer Lock

Material: ABS
Color: Clear
Sterilization: Eto, Gamma
ID: 4.5 mm



LM-050

Male Luer Lock

Material: Acrylic
Color: Clear
Sterilization: Eto
ID: 0.188"



LM-063

Male Luer

Material: MABS
Color: Clear
Sterilization: Eto, Gamma
ID: 5.0 mm



LM-025

Male Luer

Material: ABS
Color: Clear
Sterilization: Eto, Gamma
ID: 5.5 mm | OD: 6.8 mm



DirectMed

+1 (516) 656-3377 • sales@directmed.com
www.directmed.com

LM-037

Male Luer, Rotating

Material: ABS, PP
Color: Clear, White
Sterilization: Eto, Gamma
ID: 5.5 mm



LM-037B

Male Luer, Rotating

Material: ABS, PP
Color: Blue, Clear
Sterilization: Eto, Gamma
ID: 5.5 mm



LM-038

Male Luer, Rotating

Material: ABS, PP
Color: Clear, Red
Sterilization: Eto, Gamma
ID: 5.5 mm



LM-067

Male Luer Lock, Rotating,
Self-Ejecting

Material: MABS, PVC
DEHP Free
Color: Blue, Natural
Sterilization: Eto
ID: 5.5 mm



LM-024

Male Luer Lock

Material: ABS
Color: Clear
Sterilization: Eto, Gamma
ID: 6.8 mm



LM-030

Male Luer, Rotating,
Self-Ejecting

Material: PC, PVC
DEHP Free
Color: Natural
Sterilization: Eto, E-Beam,
Gamma
ID: 6.8 mm



LM-030B

Male Luer, Rotating,
Self-Ejecting

Material: PC, PVC
DEHP Free
Color: Blue, Natural
Sterilization: Eto, Gamma
ID: 6.8 mm



LM-030R

Male Luer, Rotating,
Self-Ejecting

Material: PC, PVC
DEHP Free
Color: Natural, Red
Sterilization: Eto, E-Beam,
Gamma
ID: 6.8 mm



LM-036

Male Luer, Rotating

Material: PC, PVC
DEHP Free
Color: Clear
Sterilization: Eto
ID: 6.8 mm



LM-066

Male Luer Lock, Rotating,
Self-Ejecting

Material: MABS, PVC
DEHP Free
Color: Natural, Red
Sterilization: Eto
OD: 6.8 mm



LM-097

Male Luer Lock, Rotating

Material: MABS, PVC
Color: Blue, Natural
Sterilization: Eto
ID: 6.8 mm



LM-126

Male Luer Lock, Rotating

Material: ABS, PVC
DEHP Free
Color: Blue, Clear
Sterilization: Eto, Gamma
ID: 6.8 mm



LM-128

Male Luer Lock, Rotating

Material: ABS
Color: Blue, Clear
Sterilization: Eto, Gamma
ID: 6.8 mm



LB-014

Barbed Male Luer,
1/16" ID Barb

Material: PC
Color: Tinted
Sterilization: Eto, Gamma



LB-131A

Barbed Male Luer Slip,
1/16" OD Barb

Material: PC
Color: Light Blue
Sterilization: Eto, Gamma



LB-036

NEW

Barbed Male Luer,
1.6 mm (1/16") OD Barb

Material: Nylon
Color: White
Sterilization: Eto



LB-148

Barbed Male Luer,
3/32" OD Barb

Material: PC
Color: Clear
Sterilization: Eto, Steam



LB-037

NEW

Barbed Male Luer,
2.4 mm (3/32") OD Barb

Material: Nylon
Color: White
Sterilization: Eto



LB-015

Barbed Male Luer,
1/8" OD Barb

Material: PC
Color: Tinted
Sterilization: Eto, Gamma



LB-117A

Barbed Male Luer,
1/8" OD Barb

Material: Nylon
Color: White
Sterilization: Eto



LB-118A

Barbed Male Luer,
1/8" OD Barb

Material: PP
Color: Clear
Sterilization: Eto, E-Beam,
Gamma, Steam



LB-132A

Barbed Male Luer Slip,
1/8" OD Barb

Material: PC
Color: Purple
Sterilization: N/A



LB-147

Barbed Male Luer,
Rotating, 1/8" OD Barb

Material: PC
Color: Tinted
Sterilization: Eto, Gamma,
Steam



LB-150

Barbed Male Luer,
1/8" OD Barb

Material: PC
Color: Clear
Sterilization: Eto, Steam



LB-227

Barbed Male Luer,
Rotating, 1/8" OD Barb

Material: Nylon
Color: White
Sterilization: Eto



LB-038

NEW

Barbed Male Luer,
3.2 mm (1/8") OD Barb

Material: Nylon
Color: White
Sterilization: Eto



LB-029

Barbed Male Luer,
5/32" OD Barb

Material: Nylon
Color: Black
Sterilization: Eto, Gamma,
Steam



LB-145

Barbed Male Luer,
Rotating, 5/32" OD Barb

Material: PC
Color: Purple
Sterilization: Eto, Gamma,
Steam



LB-146

Barbed Male Luer,
Rotating, 5/32" OD Barb

Material: ABS
Color: Clear
Sterilization: Eto



LM-084

Barbed Male Luer Lock,
5/32" OD Barb

Material: Nylon
Color: White
Sterilization: Eto, Gamma,
Steam



LB-039

NEW

Barbed Male Luer,
4.0 mm (5/32") OD Barb

Material: Nylon
Color: White
Sterilization: Eto



LB-108

Barbed Male Luer,
Rotating, 3/16" OD Barb

Material: ABS
Color: Clear
Sterilization: Eto



LB-120A

Barbed Male Luer,
3/16" OD Barb

Material: Nylon
Color: White
Sterilization: Eto, Gamma



LB-121A

Barbed Male Luer,
3/16" OD Barb

Material: PP
Color: Clear
Sterilization: Eto



LB-122A

Barbed Male Luer,
3/16" OD Barb

Material: PC
Color: Clear
Sterilization: Eto, Steam



LB-130A

Barbed Male Luer,
Rotating, 3/16" OD Barb

Material: PC
Color: Tinted
Sterilization: Eto, Gamma,
Steam



LB-006

Barbed Male Luer Lock,
1/4" OD Barb

Material: PP
Color: Natural
Sterilization: Eto, Gamma



LB-006PC

Barbed Male Luer Lock,
1/4" OD Barb

Material: PC
Color: Clear
Sterilization: Eto, Gamma,
Steam



LB-133A

Barbed Male Luer Slip,
1/4" OD Barb

Material: PC
Color: Tinted
Sterilization: Eto, Gamma



LF-007

Female Luer, Microbore

Material: PC
Color: Clear
Sterilization: Eto, Steam
ID: 0.74 mm



LF-079

Female Luer, Microbore

Material: PC
Color: Blue
Sterilization: Eto, Steam
ID: 0.78 mm



LF-079CL

Female Luer, Microbore

Material: PC
Color: Natural
Sterilization: Eto, Steam
ID: 0.78 mm



LF-012

Female Luer, Microbore

Material: PC
Color: Clear
Sterilization: Eto, Steam
ID: 0.9 mm



LF-013

Female Luer Lock, Microbore

Material: PC
Color: Clear
Sterilization: Eto, Steam
ID: 0.056"



LF-091GS

Female Luer Lock, Microbore

Material: PC
Color: Clear
Sterilization: Eto, Gamma
ID: 0.057"



LF-059

Female Luer, Microbore

Material: PC
Color: Clear
Sterilization: Eto, Steam
ID: 1.45 mm



LF-178

Female Luer

Material: COPE
Color: Clear
Sterilization: N/A
ID: 0.078"



LF-018

Female Luer

Material: PVC
Color: Clear
Sterilization: Eto
ID: 2.0 mm



LF-021DEHT

Female Luer

Material: PVC-DEHT
Color: Clear
Sterilization: Eto, E-Beam, Gamma
ID: 2.0 mm



LF-028

Female Luer

Material: PVC
Color: Clear
Sterilization: Eto
ID: 2.0 mm



LF-135

Female Luer Lock

Material: PVC DEHP Free
Color: Light Blue
Sterilization: Eto, E-Beam, Gamma
ID: 2.0 mm



LF-088

Female Luer

Material: PVC DEHP Free
Color: Light Blue
Sterilization: Eto, Gamma
ID: 2.1 mm



LF-094

Female Luer

Material: PVC
Color: Clear
Sterilization: Eto
ID: 2.4 mm



LF-094GS

Female Luer

Material: PVC DEHP Free
Color: Light Blue
Sterilization: Eto, Gamma
ID: 2.4 mm



LF-112

Female Luer

Material: PVC
Color: Clear
Sterilization: Eto
ID: 2.4 mm



LF-118

Female Luer

Material: PVC
Color: Brown
Sterilization: Eto
ID: 2.4 mm



LF-022DEHT

Female Luer

Material: PVC-DEHT
Color: Clear
Sterilization: Eto, E-Beam, Gamma
ID: 2.5 mm



LF-029A

Female Luer

Material: PVC
Color: Clear
Sterilization: Eto
ID: 2.5 mm



LF-029DF

Female Luer

Material: PVC DEHP Free
Color: Clear
Sterilization: Eto, E-Beam, Gamma
ID: 2.5 mm



LF-097

Female Luer

Material: PC
Color: Clear
Sterilization: Eto, E-Beam, Gamma, Steam
ID: 2.5 mm



LF-106

Female Luer

Material: PVC DEHP Free
Color: Clear
Sterilization: Eto
ID: 2.5 mm



LF-113

Female Luer

Material: PVC
Color: Clear
Sterilization: Eto
ID: 2.8 mm



LF-023DEHT

Female Luer

Material: PVC-DEHT
Color: Clear
Sterilization: Eto, E-Beam, Gamma
ID: 3.0 mm



LF-023DHPF

Female Luer

Material: PVC DEHP Free
Color: Clear
Sterilization: Eto, E-Beam, Gamma
ID: 3.0 mm



LF-048

Female Luer

Material: Acrylic
Color: Clear
Sterilization: Eto, E-Beam, Gamma
ID: 3.0 mm



LF-066

Female Luer

Material: PVC
Color: Clear
Sterilization: Eto, Gamma
OD: 3.0 mm



LF-066B

Female Luer

Material: PVC
Color: Blue
Sterilization: Eto, Gamma
OD: 3.0 mm



LF-066R

Female Luer

Material: PVC
Color: Red
Sterilization: Eto, Gamma
OD: 3.0 mm



LF-071

Female Luer Lock, High Pressure

Material: PVC DEHP Free
Color: Clear
Sterilization: Eto
ID: 3.0 mm



LF-034

Female Luer

Material: PC
Color: Clear
Sterilization: Eto
ID: 0.125"



LF-003

Female Luer

Material: PVC DEHP Free
Color: Natural
Sterilization: Eto
ID: 3.5 mm



LF-005

Female Luer

Material: PVC
Color: Clear
Sterilization: Eto
ID: 3.5 mm



LF-093

Female Luer

Material: PVC DEHP Free
Color: Light Blue
Sterilization: Eto, Gamma
ID: 3.5 mm



LF-114

Female Luer

Material: PVC
Color: Clear
Sterilization: Eto
ID: 3.5 mm



LF-039

Female Luer

Material: PC
Color: Clear
Sterilization: Eto, Gamma
ID: 4.0 mm



Female Luers - 4.0 to 5.0 mm ID

87

LF-042

Female Luer Lock

Material: ABS
Color: Clear
Sterilization: Eto, E-Beam,
Gamma

ID: 4.0 mm



LF-001

Female Luer

Material: PC
Color: Clear
Sterilization: Eto, E-Beam,
Gamma, Steam

ID: 4.1 mm



LF-002

Female Luer

Material: MABS
Color: Clear
Sterilization: Eto, E-Beam,
Gamma

ID: 4.1 mm



LF-002PVCG

Female Luer

Material: PVC DEHP Free
Color: Light Blue
Sterilization: Eto, E-Beam,
Gamma

ID: 4.1 mm



LF-002PVCNG

Female Luer

Material: PVC DEHP Free
Color: Clear
Sterilization: Eto
ID: 4.1 mm



LF-024

Female Luer

Material: PVC
Color: Clear
Sterilization: Eto
ID: 4.1 mm



LF-024DHPF

Female Luer

Material: PVC DEHP Free
Color: Clear
Sterilization: Eto
ID: 4.1 mm



LF-074

Female Luer Lock,
High Pressure

Material: PVC DEHP Free
Color: Clear
Sterilization: Eto
ID: 4.1 mm



LF-095

Female Luer

Material: PC
Color: Clear
Sterilization: Eto, Gamma
ID: 4.1 mm



LF-102

Female Luer Lock

Material: MABS
Color: Clear
Sterilization: Eto, Gamma
ID: 4.1 mm



LF-107

Female Luer

Material: PVC DEHP Free
Color: Clear
Sterilization: Eto, Gamma
ID: 4.1 mm



LF-115

Female Luer

Material: PVC
Color: Clear
Sterilization: Eto
ID: 4.1 mm



LF-122

Female Luer Lock

Material: PVC DEHP Free
Color: Clear
Sterilization: Eto
ID: 4.1 mm



LF-041

Female Luer Lock

Material: PC
Color: Clear
Sterilization: Eto, Steam
ID: 3/16"



LF-041PVC

Female Luer Lock

Material: PVC DEHP Free
Color: Light Blue
Sterilization: Eto, Gamma
ID: 3/16"



LF-116

Female Luer

Material: PVC
Color: Natural
Sterilization: Eto
ID: 4.8 mm



LF-016

Female Luer

Material: PVC DEHP Free
Color: Clear
Sterilization: Eto
ID: 5.0 mm



LF-016DHPFG

Female Luer

Material: PVC DEHP Free
Color: Clear
Sterilization: Eto, E-Beam,
Gamma
ID: 5.0 mm



DirectMed

+1 (516) 656-3377 • sales@directmed.com
www.directmed.com

LF-064

Female Luer

Material: MABS
Color: Clear
Sterilization: Eto, Gamma
ID: 5.0 mm



LF-125

Female Luer

Material: PVC
Color: Clear
Sterilization: Eto
ID: 5.0 mm



LF-117

Female Luer

Material: PVC
Color: Clear
Sterilization: Eto
ID: 5.3 mm



LF-004

Female Luer Lock

Material: PC
Color: Clear
Sterilization: Eto, Gamma
ID: 5.5 mm



LF-008

Female Luer

Material: PVC
Color: Clear
Sterilization: Eto
ID: 5.5 mm



LF-026DEHT

Female Luer Lock

Material: PVC-DEHT
Color: Clear
Sterilization: Eto, E-Beam, Gamma
ID: 5.5 mm



LF-026RNG

Female Luer Lock

Material: PVC
Color: Red
Sterilization: Eto
ID: 5.5 mm



LF-031

Female Luer

Material: PVC
Color: Clear
Sterilization: Eto
ID: 5.5 mm



LF-103

NEW

Female Luer Lock

Material: PVC DEHP Free
Color: Clear
Sterilization: Eto
ID: 5.5 mm



LF-155

Female Luer

Material: PVC DEHP Free
Color: Clear
Sterilization: Eto, E-Beam, Gamma
ID: 5.5 mm



LF-010

Female Luer

Material: PC
Color: Clear
Sterilization: Eto, E-Beam, Gamma, Steam
ID: 6.8 mm



LF-027A

Female Luer Lock

Material: PVC
Color: Clear
Sterilization: Eto
ID: 6.8 mm



LF-027BDEHT

Female Luer Lock

Material: PVC-DEHT
Color: Blue
Sterilization: Eto, E-Beam, Gamma
ID: 6.8 mm



LF-027DEHT

Female Luer Lock

Material: PVC-DEHT
Color: Clear
Sterilization: Eto, E-Beam, Gamma
ID: 6.8 mm



LF-027DF

Female Luer Lock

Material: PVC DEHP Free
Color: Clear
Sterilization: Eto, E-Beam, Gamma
ID: 6.8 mm



LF-027DHPF

Female Luer Lock

Material: PVC DEHP Free
Color: Clear
Sterilization: Eto
ID: 6.8 mm



LF-027RDEHT

Female Luer Lock

Material: PVC-DEHT
Color: Red
Sterilization: Eto, E-Beam, Gamma
ID: 6.8 mm



LF-011

Female Luer

Material: PVC DEHP Free
Color: Clear
Sterilization: Eto
ID: 7.0 mm



LB-040

NEW

Barbed Female Luer,
1.6 mm (1/16") OD Barb

Material: Nylon
Color: White
Sterilization: Eto



LB-011

Barbed Female Luer,
1/16" OD Barb

Material: PC
Color: Tinted
Sterilization: Eto, Gamma



LB-041

NEW

Barbed Female Luer,
2.4 mm (3/32") OD Barb

Material: Nylon
Color: White
Sterilization: Eto



LB-149

Barbed Female Luer,
3/32" OD Barb

Material: PP
Color: Clear
Sterilization: Eto, E-Beam,
Gamma, Steam



LB-010

Barbed Female Luer,
1/8" OD Barb | ID: 3.0 mm

Material: Nylon
Color: White
Sterilization: Eto, Gamma,
Steam



LB-135A

Barbed Female Luer,
1/8" OD Barb

Material: PC
Color: Clear
Sterilization: Eto, Gamma



LB-135APP

Barbed Female Luer,
1/8" OD Barb

Material: PP
Color: Natural
Sterilization: Eto, Gamma



LF-092

Barbed Female Luer,
1/8" OD Barb | ID: 2.25 mm

Material: Nylon
Color: Blue
Sterilization: Eto



LB-028

Barbed Female Luer,
5/32" OD Barb

Material: PC
Color: Clear
Sterilization: Eto, E-Beam,
Gamma, Steam



LB-142A

Barbed Female Luer Lock,
5/32" OD Barb

Material: PP
Color: Natural
Sterilization: Eto, E-Beam,
Gamma, Steam



LB-143A

Barbed Female Luer Lock,
5/32" OD Barb

Material: PC
Color: Clear
Sterilization: Eto, Gamma,
Steam



LB-012C

Barbed Female Luer,
3/16" OD Barb

Material: PP
Color: Natural
Sterilization: Eto, E-Beam,
Gamma



LB-013A

Barbed Female Luer Lock,
3/16" OD Barb

Material: Nylon
Color: White
Sterilization: Eto, Gamma



LB-018A

Barbed Female Luer,
3/16" OD Barb

Material: PC
Color: Clear
Sterilization: Eto, E-Beam,
Gamma, Steam



LB-042

NEW

Barbed Female Luer,
4.8 mm (3/16") OD Barb

Material: Nylon
Color: White
Sterilization: Eto



LB-034

Barbed Female Luer Lock,
High Flow, 1/4" Barb

Material: PC
Color: Clear
Sterilization: Eto, Steam



LF-014

Male to Female Luer Lock

Material: PP
Color: White
Sterilization: Eto, Gamma



LF-087

Female to Female Luer

Material: PC
Color: Tinted
Sterilization: Eto, Gamma



LF-119

Female to Female Luer

Material: PP
Color: Clear
Sterilization: Eto, Steam
ID: 4.1 mm



LF-119GS

Female to Female Luer

Material: PP
Color: Clear
Sterilization: Eto, E-Beam, Gamma, Steam
ID: 4.1 mm



LF-120B

Female to Female Luer

Material: PVC DEHP Free
Color: Clear
Sterilization: Eto



LF-123

NEW

Female to Female Luer, Double

Material: PP
Color: Red
Sterilization: Eto



LL-002

Male to Female Luer Lock

Material: ABS
Color: Clear
Sterilization: Eto, Gamma



LL-002PC

Male to Female Luer Lock

Material: PC
Color: Clear
Sterilization: Eto, E-Beam, Gamma



LL-003

Female to Female Luer

Material: ABS
Color: Natural
Sterilization: Eto
ID: 3.0 mm



LL-006

NEW

Luer, Elbow, Male Slip to Female Slip

Material: PC
Color: Clear
Sterilization: Eto, Gamma, Steam



LL-008

Male to Female Luer, Elbow, Rotating

Material: ABS, PP
Color: Clear
Sterilization: Eto, Gamma
ID: 2.6 mm



LL-009

Male Luer, Elbow, Rotating

Material: ABS, PP
Color: Clear
Sterilization: Eto, Gamma
ID: 5.5 mm



LL-014

Female Luer Lock, Elbow

Material: ABS
Color: Clear
Sterilization: Eto, Gamma
ID: 4.1 mm



LL-015

Female Luer Lock, Elbow, to Male Luer Slip

Material: ABS
Color: Clear
Sterilization: Eto, Gamma
ID: 2.6 mm



LL-018

NEW

Male to Female Luer Lock

Material: PVC
Color: Light Blue
Sterilization: Eto, Gamma



LL-019

NEW

Female Luer to
Male Luer Slip

Material: PC
Color: Clear
Sterilization: Eto, E-Beam,
Gamma, Steam

ID: 2.6 mm



LL-026

Male to Female Luer Lock,
Rotating

Material: ABS
Color: Clear
Sterilization: Eto, Gamma
ID: 2.7 mm



LM-002

Male to Female Luer,
Rotating, High-Pressure

Material: EPDM, PC
Color: Clear, Red
Sterilization: Eto
ID: 2.7 mm
Pressure: 1200 psi



LM-003

Male to Female Luer Lock

Material: PC
Color: White
Sterilization: Eto
ID: 1.3 mm



LM-004

Male Luer, Rotating,
High-Pressure

Material: EPDM, PC
Color: Clear, Red
Sterilization: Eto
ID: 4.6 mm
Pressure: 1200 psi



LM-005

Male Luer, Rotating,
Double, High-Pressure

Material: EPDM, PC
Color: Clear, Red
Sterilization: Eto



LM-061

Male Luer with
Tubing Clip

Material: PC
Color: Clear
Sterilization: Eto, Steam
ID: 2.6 mm



LM-074

Male Luer Lock, Double


Material: PC
Color: Clear
Sterilization: Eto, Gamma,
Steam



LM-074R

Male Luer Lock, Double

Material: PC
Color: Red
Sterilization: Eto, Steam
ID: 2.35 mm



LM-081

Male Luer Lock, Double

Material: PP
Color: Natural
Sterilization: Eto, Gamma



LM-082

Male Luer Lock, Rotating,
Double

Material: ABS
Color: Clear
Sterilization: Eto



LM-082A

Male Luer Lock, Rotating,
Double

Material: ABS, PVC
Color: Blue
Sterilization: N/A



LM-082B

Male Luer Lock, Rotating,
Double

Material: ABS, PC
Color: Clear, Purple
Sterilization: N/A



LM-116

Male to Female Luer Lock

Material: PC
Color: Tinted
Sterilization: Eto, Gamma
ID: 2.6, 2.8 mm



LM-174

Male Luer Lock, Double

Material: PP
Color: Natural
Sterilization: Eto, Gamma



LC-011

Luer Collar

Material: ABS
Color: Clear
Sterilization: Eto
ID: 5.0 mm



LC-012

Luer Collar

Material: PC
Color: Clear
Sterilization: Eto
ID: 5.0 mm



LC-013

Luer Collar

Material: ABS
Color: Clear
Sterilization: Eto
ID: 5.0 mm



LC-014

Luer Collar

Material: ABS
Color: White
Sterilization: Eto, Gamma,
ID: 5.0 mm



LC-015

Luer Collar

Material: PC
Color: Natural
Sterilization: Eto, Gamma,
Steam
ID: 5.0 mm



LC-016

Luer Collar

Material: PP
Color: Clear
Sterilization: Eto
ID: 5.0 mm



LC-018

Luer Collar

Material: PP
Color: Natural
Sterilization: Eto, Gamma,
Steam
ID: 6.4 mm



LM-042HUB

Luer Collar

Material: PS
Color: Natural
Sterilization: Eto, E-Beam,
Gamma
ID: 5.9 mm



LF-032

Female Luer to Luer Slip

Material: ABS
Color: White
Sterilization: Eto, Gamma
ID: 3.5 - 4.1 mm



LM-011

Male Luer Slip

Material: PVC
Color: Clear
Sterilization: Eto
ID: 1.5 mm



LS-001A

Male Luer Slip, Double

Material: ABS
Color: White
Sterilization: Eto, Gamma
ID: 2.3 mm



LS-001AC

Male Luer Slip, Double

Material: ABS
Color: Clear
Sterilization: Eto, Gamma
ID: 2.3 mm



LS-002

Male Luer Slip

Material: ABS
Color: Clear
Sterilization: Eto, Gamma
ID: 2.5 mm



LS-003

Male Luer Slip

Material: ABS
Color: Clear
Sterilization: Eto, Gamma
ID: 4.1 mm



LS-004

Male Luer Slip to
Female Luer

Material: HDPE
Color: Natural
Sterilization: Eto



LS-005DEHT

Female Luer to
Male Luer Slip

Material: PVC-DEHT
Color: Clear
Sterilization: Eto
ID: 3.5 mm



LS-007

Male Luer Slip

Material: ABS
Color: Clear
Sterilization: Eto, Gamma
ID: 6.8 mm



NH-014

NEW

Needle Hub

Material: PC
Color: Clear
Sterilization: Eto
ID: 14 Gauge



NH-015

NEW

Needle Hub

Material: PC
Color: Clear
Sterilization: Eto
ID: 15 Gauge



NH-016

NEW

Needle Hub

Material: PC
Color: White
Sterilization: Eto
ID: 16 Gauge



NH-017

NEW

Needle Hub

Material: PC
Color: Red
Sterilization: Eto
ID: 17 Gauge



NH-018

NEW

Needle Hub

Material: PC
Color: Pink
Sterilization: Eto
ID: 18 Gauge



NH-019

NEW

Needle Hub

Material: PC
Color: Beige
Sterilization: Eto
ID: 19 Gauge



NH-020

NEW

Needle Hub

Material: PC
Color: Yellow
Sterilization: Eto
ID: 20 Gauge



NH-021

NEW

Needle Hub

Material: PC
Color: Light Green
Sterilization: Eto
ID: 21 Gauge



NH-022

NEW

Needle Hub

Material: PC
Color: Black
Sterilization: Eto
ID: 22 Gauge



NA-001

NRFit Check Valve with Female Connector

Material: Acrylic, Silicone
Color: Clear, Yellow
Sterilization: Eto, Gamma
ID: 3.6 mm



NA-002

NRFit Inline IV Filter

Material: Acrylic, LDPE, MABS, PES
Color: Clear, Yellow
Sterilization: Eto, Gamma
Microns: 0.2 µm



NA-003

NRFit Rotating Hub for NA-028 and NA-029

Material: ABS
Color: Clear
Sterilization: Eto
ID: 3.7 mm



NA-004

NRFit Female Luer Lock

Material: PVC DEHP Free
Color: Clear
Sterilization: Eto, E-Beam, Gamma
ID: 3.0 mm



NA-005

NRFit Female Luer Lock

Material: PVC DEHP Free
Color: Clear
Sterilization: Eto, E-Beam, Gamma
ID: 2.5 mm



NA-010

NRFit Female Luer


Material: PC
Color: Clear
Sterilization: Eto, Gamma
ID: 2.0 mm



NA-011

NRFit Female Luer

Material: PVC
Color: Clear
Sterilization: Eto
ID: 2.4 mm



NA-020

NRFit Rotating Male Luer

Material: ABS
Color: Clear
Sterilization: Eto, Gamma
ID: 2.5 mm



NA-021

NRFit Rotating Male Luer

Material: ABS
Color: Clear
Sterilization: Eto, Gamma
ID: 3.5 mm



NA-023

NRFit Male Luer

Material: ABS
Color: Clear
Sterilization: Eto, Gamma
ID: 2.0 mm



NA-025

NRFit Male Luer

Material: ABS
Color: Clear
Sterilization: Eto, Gamma
ID: 2.5 mm



NA-026

NRFit Male Luer

Material: ABS
Color: Clear
Sterilization: Eto, Gamma
ID: 3.6 mm



NA-027

NRFit Male Luer
Material: PVC
Color: Clear
Sterilization: Eto
ID: 2.4 mm



NA-028

NRFit Male Luer Slip
Material: PVC
Color: Clear
Sterilization: Eto
ID: 2.4 mm



NA-029

NRFit Male Luer Slip
Material: PVC
Color: Clear
Sterilization: Eto
ID: 4.1 mm



NA-030

NRFit Non-Vented Cap for
Female Connectors
Material: ABS
Color: White
Sterilization: Eto, Gamma



NA-031

NRFit Non-Vented Cap for
Male Connectors
Material: PC
Color: Yellow
Sterilization: Eto, Gamma



NA-032

NRFit Vented Cap with
Filter for Male Luers
Material: Acrylic, HDPE
Color: Natural
Sterilization: N/A
Microns: 3 µm



NA-034

NRFit Non-Vented Cap for
Female Luers
Material: ABS
Color: White
Sterilization: Eto, Gamma



NA-036

NRFit Non-Vented Cap for
Female Luers
Material: ABS
Color: Yellow
Sterilization: Eto, Gamma



NA-037

NRFit Non-Vented Cap for
Female Luers
Material: HDPE
Color: Clear
Sterilization: Eto



NA-038

NRFit Vented Cap for
Male Luers
Material: LDPE
Color: Clear
Sterilization: Eto



ST-006

1-Way Stopcock with Luer Slip and Female Luer

Material: HDPE, PC
Color: Blue, Clear
Sterilization: Eto, Gamma



ST-006LL

1-Way Stopcock with Rotating Male Luer and Female Luer

Material: HDPE, PC, Silicone
Color: Blue, Clear, Natural
Sterilization: Eto, Gamma



ST-011NC

1-Way Stopcock with Luer Slip

Material: HDPE, PC
Color: Clear, Natural
Sterilization: Eto, Gamma



ST-033

4-Way Stopcock with 2 Female Luers

Material: HDPE, PC, Silicone
Color: Blue, Clear, Natural
Sterilization: Eto, Gamma
ID: 5.5 mm



ST-034

4-Way Stopcock with Rotating Male Luer and 2 Female Luers

Material: HDPE, PC, Silicone
Color: Clear, Natural, Yellow
Sterilization: Eto, Gamma



ST-035

4-Way Stopcock with Rotating Male Luer and 2 Female Luers

Material: HDPE, PC, Silicone
Color: Blue, Clear, Natural
Sterilization: Eto, Gamma



ST-036

4-Way Stopcock with Luer Slip and 2 Female Luers

Material: HDPE, PC, Silicone
Color: Blue, Clear, Natural
Sterilization: Eto, Gamma



ST-037

4-Way Stopcock with Caps and 2 Female Luers

Material: HDPE, PC, PP, Silicone
Color: Blue, Clear
Sterilization: Eto, Gamma
ID: 4.1 mm



ST-037NC

4-Way Stopcock with 2 Female Luers

Material: HDPE, PC, Silicone
Color: Blue, Clear
Sterilization: Eto, Gamma
ID: 4.1 mm



ST-039

3-Way Stopcock with Caps, 2 Female Luers and Male Luer Lock

Material: HDPE, LDPE, PC, PP, Silicone
Color: Blue, Clear, Natural, White
Sterilization: Eto, Gamma



ST-059

3-Way Stopcock with Caps, Rotating Male Luer Lock, and 2 Female Luer Locks

Material: HDPE, LDPE, PC, PP, Silicone
Color: Clear, Natural, White
Sterilization: Eto, Gamma
ID: 4.1 mm



ST-064

3-Way Stopcock with Male Luer Lock and 2 Female Luers

Material: HDPE, PC
Color: Clear, White
Sterilization: Eto, Gamma
ID: 1.7 mm



ST-066

3-Way Stopcock with Caps, Male Luer Lock and 2 Female Luers

Material: HDPE, PC
Color: Natural, Red, White
Sterilization: Eto, Gamma



ST-067

4-Way Stopcock, Photosensitive

Material: PC
Color: Clear, Natural, Yellow
Sterilization: Eto, E-Beam, Gamma



ST-070

3-Way Stopcock with 2 Female Luer Locks

Material: HDPE, PC
Color: Clear, White
Sterilization: Eto, Gamma



ST-071

3-Way Stopcock with 2 Female Luer Locks

Material: HDPE, PC, PP
Color: Black, Clear, Natural, White
Sterilization: Eto, Gamma
ID: 4.0 mm



ST-072

3-Way Stopcock, High Pressure, with Female Luer and Rotating Male Luer Lock

Material: PC, POM
Color: Blue, Clear
Sterilization: Eto
Pressure: 1200 psi



ST-077B

3-Way Stopcock with Caps, Male Luer, and 2 Female Luers

Material: HDPE, PC, PP
Color: Blue, Clear, Tinted
Sterilization: Eto
Pressure: 200 kpa
Flow Rate: 2,600 mL/10 mins



ST-077R

3-Way Stopcock with Caps, Male Luer, and 2 Female Luers

Material: HDPE, PC, PP
Color: Clear, Red, Tinted
Sterilization: Eto
Pressure: 200 kpa
Flow Rate: 2,600 mL/10 mins



ST-077Y

3-Way Stopcock with Caps, Male Luer, and 2 Female Luers

Material: ABS, PC, POM
Color: Clear, Green, Yellow
Sterilization: Eto



ST-078

1-Way Stopcock, Female Luer to Male Luer

Material: HDPE, PC, Silicone
Color: Blue, Clear
Sterilization: Eto, Gamma



ST-079

1-Way Stopcock

Material: HDPE, PC, Silicone
Color: Clear, Natural, White
Sterilization: Eto



ST-081

3-Way Stopcock, High Flow, Large Bore

Material: HDPE, PC, PP, PSU
Color: Clear, Tinted, White
Sterilization: Eto, Gamma



ST-082

NRFit 4-Way Stopcock

Material: HDPE, PC, Silicone
Color: Tinted, Yellow
Sterilization: Eto, Gamma



ST-083

NEW

3-Way Stopcock with Cap,
Male Luer Slip, and
2 Female Luers

Material: HDPE, LDPE, PC
Color: Blue, Clear, White
Sterilization: Eto, Steam
ID: 1.9 mm



ST-087

NEW

4-Way Stopcock with
Needleless Port

Material: HDPE, PC,
Silicone
Color: Blue, Tinted
Sterilization: N/A



ST-088

NEW

3-Way Stopcock with
3 Female Luers,
High Pressure

Material: PC, POM
Color: Clear, White
Sterilization: Eto
Pressure: 200 psi



ST-154

3-Way Stopcock with Cap,
2 Female Luers and
Male Luer Lock

Material: ABS, HDPE, PC
Color: Blue, Clear
Sterilization: Eto, Gamma



ST-162

1-Way Stopcock with
Female Luer and Rotating
Male Luer Lock

Material: HDPE, PC
Color: Natural, Purple,
White
Sterilization: Eto, Gamma



SR-001

3-Gang Stopcock Manifold
with Female Luers

Material: HDPE, PP
Color: Blue, Natural, Red,
White
Sterilization: Eto, Gamma



SR-004

3-Gang Stopcock Manifold
with Female Luers

Material: HDPE, PC
Color: Blue, Clear, Red,
White
Sterilization: Eto



SR-005

5-Gang Stopcock Manifold
with Female Luers

Material: HDPE, PC
Color: Blue, Clear, Green,
Red, White, Yellow
Sterilization: Eto



SR-024

3-Gang Stopcock Manifold,
Male to Female Luer

Material: ABS, COPE, PP
Color: Blue, Clear
Sterilization: Eto



SR-025

3-Gang Stopcock Manifold

Material: PC
Color: Blue, Clear, Red,
White
Sterilization: N/A



SR-026

4-Gang Stopcock Manifold
with Female Luers

Material: HDPE, PC
Color: Blue, Clear, Red,
White, Yellow
Sterilization: Eto



SR-027

5-Gang Stopcock Manifold

Material: HDPE, PC
Color: Blue, Clear, Green,
Red, White, Yellow
Sterilization: Eto



SY-026

Syringe

Material: ABS, PC, PIP
Color: Black, Clear, Red
Sterilization: Eto
Volume: 10cc



SY-027

Syringe

Material: ABS, PC, PIP
Color: Black, Clear, White
Sterilization: Eto
Volume: 20cc



SY-032

Syringe, LOR

Material: PP
Color: Black, Clear, White
Sterilization: Eto, Steam
Volume: 10cc



SY-034

Syringe

Material: ABS, COPE, Silicone
Color: Black, Blue, Clear
Sterilization: Eto, E-Beam, Gamma
Volume: 1cc



SY-035

Syringe

Material: ABS, COPE, Silicone
Color: Black, Clear, Yellow
Sterilization: Eto, E-Beam, Gamma
Volume: 1cc



SY-037

Syringe

Material: PP, Rubber
Color: Black, Clear, Natural
Sterilization: Eto
Volume: 30 mL



SY-045

Syringe

Material: ABS, PC, Silicone
Color: Black, Clear, White
Sterilization: Eto
Volume: 1 mL



SY-046

Syringe

Material: ABS, COPE, Silicone
Color: Black, Clear, White
Sterilization: Eto
Volume: 1 mL



SY-047

Syringe

Material: ABS, PC, Silicone
Color: Black, Clear, White
Sterilization: Eto
Volume: 1 mL



SY-063

Control Syringe

Material: ABS, IR, PC, Silicone
Color: Black, Clear, Red, White
Sterilization: N/A
Volume: 12 mL



SY-960

NEW

Syringe, Eccentric Luer

Material: PP, Rubber, Silicone
Color: Black, Clear, Natural
Sterilization: Eto
Volume: 60 mL



VAC160P

VACLOK Syringe

Material: ABS, PC, Silicone
Color: Black, Clear, White
Sterilization: Eto, Gamma
Volume: 60 mL



TT-002A

Tubing

Material: PVC DEHP Free
Color: Clear
Sterilization: Eto, E-Beam
ID: 3.0 mm | OD: 4.1 mm
Dimensions: 1,000 m Length

TT-006A

Tubing

Material: PVC DEHP Free
Color: Clear
Sterilization: Eto, E-Beam
ID: 16 mm | OD: 18 mm
Dimensions: 6 cm Length

TT-007A

Tubing

Material: PVC DEHP Free
Color: Clear
Sterilization: Eto, E-Beam
ID: 3.0 mm | OD: 4.1 mm
Dimensions: 10 cm (4") Length

TT-008A

Tubing

Material: PVC DEHP Free
Color: Clear
Sterilization: Eto, E-Beam
ID: 3.0 mm | OD: 4.1 mm
Dimensions: 15 cm Length

TT-041

Tubing

Material: PVC DEHP Free
Color: Clear
Sterilization: Eto, E-Beam
ID: 3.0 mm | OD: 4.1 mm
Dimensions: 18 cm (7") Length

TT-042

Tubing

Material: PVC DEHP Free
Color: Clear
Sterilization: Eto, E-Beam
ID: 3.0 mm | OD: 4.1 mm
Dimensions: 40.6 cm Length

TT-046

Tubing

Material: PVC DEHP Free
Color: Clear
Sterilization: Eto, E-Beam,
Gamma
ID: 1.0 mm | OD: 2.5 mm

TT-048

Tubing

Material: PVC DEHP Free
Color: Clear
Sterilization: Eto, E-Beam
ID: 3.0 mm | OD: 4.1 mm
Dimensions: 61 cm Length

TT-059

Tubing

Material: PVC DEHP Free
Color: Clear
Sterilization: N/A
ID: 3.0 mm | OD: 4.1 mm
Dimensions: 1,500 m Length

TT-061DEHT

Tubing

Material: PVC-DEHT
Color: Clear
Sterilization: Eto, E-Beam,
Gamma
ID: 1.0 mm | OD: 2.5 mm

TT-062DEHT

Tubing

Material: PVC-DEHT
Color: Clear
Sterilization: Eto, E-Beam,
Gamma
ID: 3.5 mm | OD: 5.5 mm
Dimensions: 700 m Length

TT-063DEHT

Tubing

Material: PVC-DEHT
Color: Clear
Sterilization: Eto, E-Beam,
Gamma
ID: 4.3 mm | OD: 6.8 mm

TT-070

Tubing

Material: PVC DEHP Free
Color: Clear
Sterilization: Eto
ID: 1.27 mm | OD: 2.4 mm

TT-072

Tubing

Material: PVC DEHP Free
Color: Clear
Sterilization: Eto, E-Beam
ID: 3.0 mm | OD: 4.1 mm
Dimensions: 2.5 cm Length

TT-073

Tubing

Material: PVC DEHP Free
Color: Clear
Sterilization: Eto, E-Beam
ID: 3.0 mm | OD: 4.1 mm

TT-074

Tubing

Material: PVC DEHP Free
Color: Clear
Sterilization: Eto, E-Beam,
Gamma
ID: 2.3 mm | OD: 4.1 mm

TT-075

Tubing

Material: PVC DEHP Free
Color: Clear
Sterilization: Eto, E-Beam,
Gamma
ID: 2.6 mm | OD: 4.1 mm

TT-076

Tubing

Material: PVC DEHP Free
Color: Clear
Sterilization: Eto, E-Beam,
Gamma
ID: 2.3 mm | OD: 4.1 mm

TT-077

Tubing

Material: PVC DEHP Free
Color: Clear
Sterilization: Eto, E-Beam
ID: 3.0 mm | OD: 4.1 mm
Dimensions: 5 cm (2") Length

TT-080

Tubing

Material: PVC DEHP Free
Color: Clear
Sterilization: Eto, E-Beam
ID: 3.0 mm | OD: 4.1 mm
Dimensions: 182.9 cm (6 ft) Length

TT-081 Tubing Material: PVC DEHP Free Color: Clear Sterilization: Eto, E-Beam ID: 3.0 mm OD: 4.1 mm Dimensions: 304.8 cm (10 ft) Length	TT-088 Tubing Material: PVC DEHP Free Color: Clear Sterilization: Eto, E-Beam, Gamma ID: 2.6 mm OD: 4.1 mm Dimensions: 2.7 cm Length	TT-089 Tubing Material: PVC DEHP Free Color: Clear Sterilization: Eto, E-Beam, Gamma ID: 2.6 mm OD: 4.1 mm Dimensions: 3.3 cm Length	TT-099 Tubing, 57SHA Material: PVC-DEHT Color: Clear Sterilization: Eto, E-Beam, Gamma ID: 6.35 mm OD: 9.52 mm Dimensions: 45.72 cm Length
TT-100 Tubing, 57SHA Material: PVC-DEHT Color: Clear Sterilization: Eto, E-Beam, Gamma ID: 6.35 mm OD: 9.52 mm Dimensions: 121.92 cm Length	TT-101 Tubing, 70SHA Material: PVC-DEHT Color: Clear Sterilization: Eto, E-Beam, Gamma ID: 6.35 mm OD: 9.52 mm Dimensions: 274.32 cm Length	TT-102 Tubing, 70SHA Material: PVC-DEHT Color: Clear Sterilization: Eto, E-Beam, Gamma ID: 6.35 mm OD: 9.52 mm Dimensions: 248.92 cm Length	TT-103 Tubing, 57SHA Material: PVC-DEHT Color: Clear Sterilization: Eto, E-Beam, Gamma ID: 3.17 mm OD: 6.35 mm Dimensions: 121.92 cm Length

DirectMed offers nearly every conceivable size and style of PVC tubing. All tubing is available in standard PVC or DEHP-Free.

Our extrusion process can be quickly and easily adjusted to create the exact tubing specification that your application calls for.

We can help you identify which PVC Tubing will work best for your application and we look forward to consulting with you regarding your project.

- PVC-DEHT or DEHP-Free
- All standard ID/OD available
- Custom colors, custom lengths
- Matched connectors



DY-048

Tubing Clip

Material: PVC
Color: Natural
Sterilization: Eto
ID: 6.0 mm



DY-049

Tubing Clip

Material: PVC
Color: Clear
Sterilization: Eto
ID: 8.0 mm



DY-050DHPF

Tubing Clip

Material: PVC DEHP Free
Color: Clear
Sterilization: Eto, E-Beam, Gamma
ID: 10.5 mm



DY-051

Tubing Clip

Material: PP
Color: Natural
Sterilization: N/A
ID: 11.6 mm



DY-077M

Tubing Clip

Material: HDPE
Color: Natural
Sterilization: N/A



TC-006

Tubing Clip, "J" with Male Luer Slip

Material: PC
Color: Clear
Sterilization: N/A



TC-117

Tubing Clip, 4-Cavity

Material: HDPE
Color: Natural
Sterilization: N/A
OD: 0.152"



TC-143

Tubing Clip, 4-Cavity

Material: HDPE
Color: White
Sterilization: N/A
OD: 0.152"



APHC0750

NEW

Hose Clamp

Material: Nylon
Color: White
Sterilization: Eto
ID: 0.73" - 0.825"



CP-002ABS

Pinch Clamp for 2.0 to 5.5 mm OD Tubing

Material: ABS
Color: White
Sterilization: Eto, Gamma



CP-004A

Pinch Clamp for 5.5 to 6.8 mm OD Tubing

Material: PP
Color: White
Sterilization: Eto, Gamma, Steam



CP-004AB

Pinch Clamp for 5.5 to 6.8 mm OD Tubing

Material: PP
Color: Blue
Sterilization: Eto, Gamma, Steam



CP-004AY

Pinch Clamp for 5.5 to 6.8 mm OD Tubing

Material: PP
Color: Yellow
Sterilization: Eto, Gamma, Steam



CP-004R

Pinch Clamp for 5.5 to 6.8 mm OD Tubing

Material: PP
Color: Red
Sterilization: Eto, Gamma, Steam



CP-006A

Pinch Clamp for 5.5 to 6.8 mm OD Tubing

Material: PP
Color: White
Sterilization: Eto, Gamma, Steam



CP-006AB

Pinch Clamp for 5.5 to 8.0 mm OD Tubing

Material: PP
Color: Blue
Sterilization: Eto, Gamma, Steam



CP-007

Pinch Clamp for 5.5 to 6.8 mm OD Tubing

Material: PP
Color: White
Sterilization: Eto, Steam



CP-008

Pinch Clamp, Non-Reopening for 5.5 to 6.8 mm OD Tubing

Material: PP
Color: White
Sterilization: Eto, Gamma, Steam



CP-008B

Pinch Clamp, Non-Reopening for 5.5 to 6.8 mm OD Tubing

Material: PP
Color: Blue
Sterilization: Eto, Gamma, Steam



CP-011A

Pinch Clamp for 0.75" OD Tubing

Material: Polyester
Color: White
Sterilization: Gamma



CP-011AB

Pinch Clamp for 0.75" OD Tubing

Material: Polyester
Color: Blue
Sterilization: Eto, Gamma, Steam



CP-014

Pinch Clamp, Non-Reopening

Material: PP
Color: Natural
Sterilization: Eto



CP-015

Umbilical Clamp, Non-Reopening

Material: POM
Color: White
Sterilization: Eto, E-Beam, Gamma, Steam



CP-017

Pinch Clamp for 3.5 to 6.8 mm OD Tubing

Material: PP
Color: White
Sterilization: Eto, E-Beam, Gamma, Steam



CP-017B

Pinch Clamp for 3.5 to 6.8 mm OD Tubing

Material: PP
Color: Blue
Sterilization: Eto, E-Beam, Gamma, Steam



CP-017R

Pinch Clamp for 3.5 to 6.8 mm OD Tubing

Material: PP
Color: Red
Sterilization: Eto, E-Beam, Gamma, Steam



CP-020

Pinch Clamp for 5.5 to 6.8 mm OD Tubing

Material: PP
Color: White
Sterilization: Eto, Gamma, Steam



CP-024

Pinch Clamp for 2.0 to 5.5 mm OD Tubing

Material: PP
Color: Blue
Sterilization: Eto



CP-025

Pinch Clamp for 2.0 to 5.5 mm OD Tubing

Material: PP
Color: Blue
Sterilization: Eto, Gamma, Steam



CP-025P

Pinch Clamp for 2.0 to 5.5 mm OD Tubing

Material: PP
Color: Purple
Sterilization: Eto, Gamma, Steam



CP-025W

Pinch Clamp for 2.0 to 5.5 mm OD Tubing

Material: PP
Color: White
Sterilization: Eto, Gamma, Steam



CP-036BL

Pinch Clamp for 5.5 to 6.8 mm OD Tubing

Material: PP
Color: Black
Sterilization: Eto, Steam



CP-036BR

Pinch Clamp for 5.5 to 6.8 mm OD Tubing

Material: PP
Color: Brown
Sterilization: Eto



CP-036GR

Pinch Clamp for 5.5 to 6.8 mm OD Tubing

Material: PP
Color: Green
Sterilization: Eto



CP-036YE

Pinch Clamp for 5.5 to 6.8 mm OD Tubing

Material: PP
Color: Yellow
Sterilization: Eto



CP-038

Pinch Clamp for 4.1 to 5.5 mm OD Tubing

Material: PP
Color: White
Sterilization: Eto, Steam



CP-039

Pinch Clamp for 5.5 to 6.8 mm OD Tubing

Material: PP
Color: White
Sterilization: Eto, Gamma, Steam



CP-040

Pinch Clamp for 2.5 to 4.1 mm OD Tubing

Material: PP
Color: White
Sterilization: Eto, Gamma, Steam



CP-042

Pinch Clamp for 2.5 to 4.1 mm OD Tubing

Material: POM
Color: Natural
Sterilization: Eto



CP-043

Pinch Clamp for 2.5 to 4.1 mm OD Tubing

Material: PP
Color: Purple
Sterilization: Eto, Steam



CP-044

Pinch Clamp for 2.5 to 4.1 mm OD Tubing

Material: PP
Color: Green
Sterilization: Eto, Steam



CP-045

Pinch Clamp for 2.5 to 4.1 mm OD Tubing

Material: PP
Color: White
Sterilization: Eto, Steam



CP-045GR

Pinch Clamp for 2.5 to 4.1 mm OD Tubing

Material: PP
Color: Green
Sterilization: Eto, Steam



CP-045YE

Pinch Clamp for 2.5 to 4.1 mm OD Tubing

Material: PP
Color: Yellow
Sterilization: Eto, Steam



CP-046

Pinch Clamp for 2.5 to 4.1 mm OD Tubing

Material: PP
Color: White
Sterilization: Eto



CP-047

Pinch Clamp for 2.5 to 4.1 mm OD Tubing

Material: PP
Color: Blue
Sterilization: Eto



CP-049

Pinch Clamp for 5.5 to 6.8 mm OD Tubing

Material: PP
Color: Red
Sterilization: Eto



CP-050

Pinch Clamp for 5.5 to 6.8 mm OD Tubing

Material: PP
Color: Blue
Sterilization: Eto



CP-051

Pinch Clamp for 5.5 to 6.8 mm OD Tubing

Material: PP
Color: White
Sterilization: Eto, Steam



CP-052

Pinch Clamp for 5.5 to 6.8 mm OD Tubing

Material: PS
Color: Blue
Sterilization: Eto, Steam



CP-053

Pinch Clamp for 5.5 to 6.8 mm OD Tubing

Material: PP
Color: Blue
Sterilization: Eto, E-Beam,
Gamma, Steam



CP-054

Pinch Clamp for 5.5 to 6.8 mm OD Tubing

Material: PP
Color: White
Sterilization: Eto, Steam



CP-055

Pinch Clamp for 5.5 to 6.8 mm OD Tubing

Material: PP
Color: White
Sterilization: Eto, Gamma,
Steam



CP-056

Open Side Pinch Clamp
for 2.0 to 4.1 mm OD
Tubing

Material: POM
Color: Red
Sterilization: Eto, Steam



CP-057

Open Side Pinch Clamp
for 2.0 to 4.1 mm OD
Tubing

Material: POM
Color: Blue
Sterilization: Eto, Steam



CP-058

Umbilical Clamp,
Non-Reopening

Material: POM
Color: Blue
Sterilization: Eto, Gamma



CP-059

Pinch Clamp for 2.5 to
4.1 mm OD Tubing

Material: POM
Color: White
Sterilization: Eto



CP-061

Pinch Clamp for 4.1 mm
OD Tubing

Material: POM
Color: Red
Sterilization: Eto



CP-063

Pinch Clamp for 4.1 mm
OD Tubing

Material: POM
Color: Blue
Sterilization: Eto



CP-064

Pinch Clamp for 4.1 mm
OD Tubing

Material: PP
Color: Blue
Sterilization: Eto



CP-067

Pinch Clamp, Ratchet, for
6.8 to 10 mm OD Tubing

Material: PP
Color: White
Sterilization: Eto, E-Beam,
Gamma, Steam



CP-068

Pinch Clamp

Material: PP
Color: White
Sterilization: Eto, Gamma,
Steam



CP-070

Pinch Clamp for 6.8 to
8.0 mm OD Tubing

Material: PP
Color: Red
Sterilization: Eto, E-Beam,
Gamma, Steam



CP-071

Pinch Clamp for 5.5 mm
OD Tubing

Material: PP
Color: White
Sterilization: Eto



CP-088B

Pinch Clamp for 2.5 to
4.1 mm OD Tubing

Material: POM
Color: Blue
Sterilization: Eto, Steam



CP-088G

Pinch Clamp for 2.5 to
4.1 mm OD Tubing

Material: POM
Color: Green
Sterilization: Eto, Steam



CP-088PPW

Pinch Clamp for 2.5 to
4.1 mm OD Tubing

Material: PP
Color: White
Sterilization: Eto, Gamma



CP-088R

Pinch Clamp for 2.5 to
4.1 mm OD Tubing

Material: POM
Color: Red
Sterilization: Eto, Steam



CP-088W

Pinch Clamp for 2.5 to 4.1 mm OD Tubing

Material: POM
Color: White
Sterilization: Eto, Steam



CP-106B

Pinch Clamp for 5.5 to 6.8 mm OD Tubing

Material: PP
Color: Blue
Sterilization: Eto, Gamma, Steam



CP-106R

Pinch Clamp for 5.5 to 6.8 mm OD Tubing

Material: PP
Color: Red
Sterilization: Eto, Gamma, Steam



CP-113

Pinch Clamp

Material: PP
Color: White
Sterilization: Eto



CP-120

Pinch Clamp for 5.5 to 6.8 mm OD Tubing

Material: PP
Color: White
Sterilization: Eto, Gamma, Steam



CP-301

Screw Tubing Clamp

Material: POM
Color: White
Sterilization: Eto



CR-002

Roller Clamp for 6.8 mm OD Tubing

Material: ABS
Color: White
Sterilization: Eto, E-Beam, Gamma



CR-002B

Roller Clamp for 5.5 to 6.8 mm OD Tubing

Material: ABS
Color: White
Sterilization: Eto



CR-002BP

Roller Clamp for 6.8 mm OD Tubing

Material: ABS
Color: Purple
Sterilization: Eto



CR-002BW

Wheel for CR-002B

Material: ABS
Color: White
Sterilization: Eto
OD: 6.8 mm



CR-002W

Wheel for CR-002

Material: ABS
Color: White
Sterilization: Eto, E-Beam, Gamma



CR-003

Roller Clamp for 4.1 mm OD Tubing

Material: ABS
Color: Orange
Sterilization: Eto



CR-003W

Wheel for CR-003

Material: ABS
Color: Natural
Sterilization: Eto



CR-004

Roller Clamp for 6.0 to 8.5 mm OD Tubing

Material: ABS, SAN
Color: White
Sterilization: Eto



CR-005

Roller Clamp for 4.1 mm OD Tubing

Material: HDPE
Color: Blue, White
Sterilization: Eto, Gamma



CR-005W

Wheel for CR-005

Material: HDPE
Color: Blue
Sterilization: Eto, Gamma
OD: 4.1 mm



CR-005WW

Wheel for CR-005

Material: HDPE
Color: White
Sterilization: Eto, Gamma
OD: 4.1 mm



CR-006

Roller Clamp for 5.5 mm OD Tubing

Material: PP
Color: White
Sterilization: Eto, Steam



CR-008

Roller Clamp for 1/4" OD Tubing

Material: ABS, PP
Color: White
Sterilization: Eto



CR-009

Roller Clamp for 4.1 mm OD Tubing

Material: ABS
Color: White
Sterilization: Eto, Gamma



CR-009ASE

Roller Clamp Assembly for 4.1 mm OD Tubing

Material: ABS
Color: White
Sterilization: Eto, Gamma



CR-009W

Wheel for CR-009

Material: ABS
Color: White
Sterilization: Eto, Gamma
OD: 4.1 mm



CR-013

Roller Clamp for 6.8 mm OD Tubing

Material: ABS
Color: White
Sterilization: Eto



CR-013W

Wheel for CR-013

Material: ABS
Color: White
Sterilization: Eto
OD: 6.8 mm



CR-014

Roller Clamp for 4.1 mm OD Tubing

Material: HDPE, PS
Color: Natural, Purple
Sterilization: Eto



CR-015

Roller Clamp for 4.1 mm OD Tubing

Material: ABS
Color: White
Sterilization: Eto



CR-017

Roller Clamp for 5.5 mm OD Tubing

Material: PP
Color: White
Sterilization: Eto, Gamma



CR-017W

Wheel for CR-017

Material: PP
Color: White
Sterilization: Eto, Gamma



CR-017ASE

Roller Clamp Assembly for 5.5 mm OD Tubing

Material: PP
Color: White
Sterilization: Eto, Gamma



CR-019BB

Roller Clamp for 4.1 mm OD Tubing

Material: ABS
Color: White
Sterilization: Eto



CR-019BBGW

Wheel for CR-019BB

Material: ABS
Color: Green
Sterilization: Eto
OD: 4.1 mm



CR-019BBWW

Wheel for CR-019BB

Material: ABS
Color: White
Sterilization: Eto
OD: 4.1 mm



CR-019BLUE

Roller Clamp for 4.1 mm OD Tubing

Material: ABS
Color: Blue
Sterilization: Eto



CR-019G

Roller Clamp for 4.1 mm OD Tubing

Material: ABS
Color: Green
Sterilization: Eto



CR-019GW

Wheel for CR-019G

Material: ABS
Color: Green
Sterilization: Eto
OD: 4.1 mm



CR-019WHEEL

Wheel for CR-019BLUE

Material: ABS
Color: White
Sterilization: Eto
OD: 4.1 mm



CR-020WHITE

Roller Clamp for 4.1 mm
OD Tubing

Material: ABS
Color: White
Sterilization: Eto



CR-021

Roller Clamp Assembly for
4.1 mm OD Tubing

Material: ABS
Color: White
Sterilization: Eto
OD: 4.1 mm



CR-022

Roller Clamp for 4.1 mm
OD Tubing

Material: ABS
Color: Blue, White
Sterilization: Eto



CR-022W

Roller Clamp for 4.1 mm
OD Tubing

Material: ABS
Color: White
Sterilization: Eto



CR-022WB

Roller Clamp for 4.1 mm
OD Tubing

Material: ABS
Color: Blue, White
Sterilization: Eto



CR-024

NEW

Roller Clamp

Material: PS
Color: White
Sterilization: Eto



CR-117

Roller Clamp for 5.5 mm
OD Tubing

Material: PP
Color: Blue, White
Sterilization: Eto, E-Beam,
Gamma



CR-118

Roller Clamp for 5.5 mm
OD Tubing

Material: PP
Color: Red, White
Sterilization: Eto, E-Beam,
Gamma



CS-001

Slide Clamp for 6.8 mm OD Tubing

Material: PE
Color: White
Sterilization: Eto, E-Beam, Gamma



CS-002

Slide Clamp for 4.1 mm OD Tubing

Material: HDPE
Color: Red
Sterilization: Eto, Gamma



CS-002N

Slide Clamp for 5.5 to 6.8 mm OD Tubing

Material: HDPE
Color: Natural
Sterilization: Eto, E-Beam, Gamma



CS-003

Slide Clamp for 6.8 mm OD Tubing

Material: PP
Color: White
Sterilization: Eto, Gamma, Steam



CS-003B

Slide Clamp for 6.8 mm OD Tubing

Material: PP
Color: Blue
Sterilization: Eto, Gamma, Steam



CS-003R

Slide Clamp for 6.8 mm OD Tubing

Material: PP
Color: Red
Sterilization: Eto, Gamma, Steam



CS-004

Slide Clamp for 6.8 mm OD Tubing

Material: PP
Color: White
Sterilization: Eto



CS-005

Slide Clamp for 6.4 mm OD Tubing

Material: HDPE
Color: Red
Sterilization: Eto, Gamma



CS-006

Slide Clamp for 6.0 mm OD Tubing

Material: HDPE
Color: Red
Sterilization: Eto, Steam



CS-007

Open Side Slide Clamp for 3.5 mm OD Tubing

Material: HDPE
Color: White
Sterilization: Eto



CS-009B

Slide Clamp for 2.4 mm OD Tubing

Material: HDPE
Color: Blue
Sterilization: Eto



CS-009DG

Slide Clamp for 2.4 mm OD Tubing

Material: HDPE
Color: Dark Green
Sterilization: Eto



CS-009LG

Slide Clamp for 2.4 mm OD Tubing

Material: HDPE
Color: Light Green
Sterilization: Eto



CS-009P

Slide Clamp for 2.4 mm OD Tubing

Material: HDPE
Color: Purple
Sterilization: Eto



CS-009R

Slide Clamp for 2.4 mm OD Tubing

Material: HDPE
Color: Red
Sterilization: Eto



CS-009W

Slide Clamp for 2.4 mm OD Tubing

Material: HDPE
Color: White
Sterilization: Eto



CS-009Y

Slide Clamp for 2.4 mm OD Tubing

Material: HDPE
Color: Yellow
Sterilization: Eto



CS-013

Open Side Slide Clamp for 1/4" (6.4 mm) OD Tubing

Material: ABS
Color: Blue
Sterilization: Eto, Gamma



CS-014

Slide Clamp for 4.1 mm OD Tubing

Material: PP
Color: Blue
Sterilization: Eto, Gamma, Steam



CS-014R

Slide Clamp for 4.1 mm OD Tubing

Material: PP
Color: Red
Sterilization: Eto, Gamma, Steam



CS-014W

Slide Clamp for 4.1 mm OD Tubing

Material: PP
Color: White
Sterilization: Eto, Gamma, Steam



CS-016

Slide Clamp for 2.5 mm OD Tubing

Material: HDPE
Color: White
Sterilization: Eto



CS-016B

Slide Clamp for 2.5 mm OD Tubing

Material: HDPE
Color: Blue
Sterilization: Eto



CS-016G

Slide Clamp for 2.5 mm OD Tubing

Material: HDPE
Color: Green
Sterilization: Eto



CS-016P

Slide Clamp for 2.5 mm OD Tubing

Material: HDPE
Color: Purple
Sterilization: Eto



CS-016R

Slide Clamp for 2.5 mm OD Tubing

Material: HDPE
Color: Red
Sterilization: Eto



CS-016Y

Slide Clamp for 2.5 mm OD Tubing

Material: HDPE
Color: Yellow
Sterilization: Eto



CS-018

Open Side Slide Clamp for 2.5 mm OD Tubing

Material: HDPE
Color: White
Sterilization: Eto



CS-018B

Open Side Slide Clamp for 2.5 mm OD Tubing

Material: HDPE
Color: Blue
Sterilization: Eto



CS-031

Open Side Slide Clamp for 4.1 mm OD Tubing

Material: HDPE
Color: White
Sterilization: Eto



CS-032

Open Side Slide Clamp for 4.1 mm OD Tubing

Material: HDPE
Color: Blue
Sterilization: Eto



CS-033

Slide Clamp for 4.1 mm OD Tubing

Material: HDPE
Color: Purple
Sterilization: Eto, Gamma



CS-034

Open Side Slide Clamp for 3.6 mm OD Tubing

Material: HIPS
Color: White
Sterilization: Eto, Gamma



CS-035

Slide Clamp for 6.0 mm OD Tubing

Material: HDPE
Color: White
Sterilization: Eto, Steam



CS-036

Slide Clamp for 10.00 mm OD Tubing

Material: HDPE
Color: White
Sterilization: Eto



CS-037

Slide Clamp for 4.1 mm OD Tubing

Material: HDPE
Color: Purple
Sterilization: Eto



CS-037B

Slide Clamp for 4.1 mm
OD Tubing

Material: HDPE
Color: Blue
Sterilization: N/A



CS-038

Slide Clamp for 4.1 mm
OD Tubing

Material: ABS
Color: Purple
Sterilization: Eto, Gamma



CS-039

Slide Clamp for 2.5 to
6.8 mm OD Tubing

Material: HDPE
Color: Natural
Sterilization: Eto, Gamma



CS-040

Slide Clamp for 10.2 mm
OD Tubing

Material: ABS
Color: White
Sterilization: Eto



CS-042

Open Side Slide Clamp for
6.8 mm OD Tubing

Material: Nylon
Color: Blue
Sterilization: Eto, Gamma,
Steam



CS-043

Open Side Slide Clamp for
4.1 mm OD Tubing

Material: PP
Color: Natural
Sterilization: Eto



CS-047

Slide Clamp

Material: HDPE
Color: White
Sterilization: Eto



CS-047G

Slide Clamp

Material: HDPE
Color: Green
Sterilization: Eto



CS-047R

Slide Clamp

Material: HDPE
Color: Red
Sterilization: Eto



CS-047Y

Slide Clamp

Material: HDPE
Color: Yellow
Sterilization: Eto



CS-048

Open Side Slide Clamp for
1.8 to 3.0 mm OD Tubing

Material: HDPE
Color: Natural
Sterilization: Eto



CV-001

Check Valve with
Male Luer Lock

Material: MABS, SAN,
Silicone, THF
Color: Clear
Sterilization: Eto, E-Beam,
Gamma

OD: 3.1 mm

Pressure:
Backflow: 6 bar (Max)
Cracking: ≤ 12 mbar



CV-007

Check Valve with
Male Luer Lock

Material: Acrylic, Silicone
Color: Clear, Green
Sterilization: Eto
ID: 4.1 mm
Pressure: 100 psi



CV-009

Check Valve, Antisiphon,
with Male Luer Lock

Material: Acrylic, Silicone
Color: Blue, Clear, Natural
Sterilization: Eto, E-Beam,
Gamma

Pressure:
Backflow: 45 psi (Max)
Cracking: 3.5 ± 1 psi
Flow Rate: 203 mL/min



CV-013

Check Valve, Dual Flow,
High Pressure

Material: MABS, Silicone
Color: Clear
Sterilization: Eto
ID: 4.1 mm
Pressure: 280 ± 50 mbar
(Opening - Tubing Conn.)
 200 ± 50 mbar
(Opening - Male Luer Lock)
83 bar / 1200 psi (Max)



CV-045

Check Valve, Dual Flow,
Male Luer Slip to
Female Luer

Material: PC, Silicone
Color: Clear
Sterilization: Eto
Pressure: 200 psi



CV-055

Check Valve, Inline

Material: Acrylic, Silicone
Color: Blue, Clear, White
Sterilization: Eto, Gamma
ID: 4.1 mm
Pressure: 45 psi



CV-060

Pressure Relief "T" Valve

Material: PC, PP, Silicone,
Stainless Steel
Color: Black, Clear, Natural
Sterilization: N/A
ID: 4.0 mm
Pressure: 2.5 ± 0.5 psi



CV-070

Check Valve, Inline

Material: Acrylic
Color: Blue, Clear, White
Sterilization: Eto, E-Beam,
Gamma
ID: 4.1 mm
OD: 0.129" - 0.132"
Pressure: 45 psi



CV-071

Check Valve, Large Bore,
High Flow

Material: MABS, Silicone
Color: Blue, Clear, Natural
Sterilization: Eto
ID: 6.5 mm
Pressure:
Backflow: 40 psi
Cracking: < 1 psi



CV-074

Check Valve, Dual Flow

Material: PC
Color: Blue
Sterilization: Eto, Steam
OD: 6.67 mm
Cracking Pressure: 30 psi



CV-075

Check Valve, Dual Flow

Material: PC
Color: Blue, Clear
Sterilization: Eto, Steam
OD: 6.6 mm
Cracking Pressure: 30 psi



CV-081

Check Valve, Luer
Activated

Material: ABS, PC, Silicone
Color: Blue, Clear, Green
Sterilization: Eto, E-Beam,
Gamma
Pressure: 45 psi
Priming Volume: < 0.3 mL



CV-084

Check Valve, Sealing Tip,
Male to Female Luer Lock

Material: MABS, Silicone,
THF, TPE

Color: Blue, Clear

Sterilization: Eto, Gamma

Pressure:

Opening: 200 ± 50
mbar

Max: 21 bar



CV-086

Check Valve

Material: Polyester, PP,
Stainless Steel

Color: White

Sterilization: Eto, Gamma,
Steam



CV-088

Check Valve with Female
Luer Lock, Needleless

Material: PP, PVC, Silicone

Color: Black, Clear

Sterilization: Eto, E-Beam,
Gamma



CV-094

Check Valve, 1/8" Barb

Material: ABS, PS, Silicone

Color: Natural, Purple

Sterilization: Eto, Gamma

OD: 1/8"

Pressure: ≤ 8 mbar
(Opening)



CV-095

Check Valve, Dual,
High Pressure

Material: PC, Silicone

Color: Clear

Sterilization: Eto, Gamma

Pressure: 300 psi (Max)



CV-096

NEW

Check Valve, Female Luer
to Luer Slip

Material: PC, Silicone

Color: Clear, Tinted

Sterilization: Eto

ID: 2.3 mm | OD: 3.5 mm

Pressure:

Opening: 6 Kpa

Backflow: > 300 Kpa

Priming Volume: 0.25 mL



CV-097

NEW

Check Valve, Dual,
High Pressure

Material: PC, Silicone

Color: Tinted

Sterilization: Eto, Gamma

Pressure: 300 psi (20 bar)



CV-101

Check Valve, Reverse Flow,
with Female Luer Lock

Material: Acrylic, Silicone

Color: Blue, Clear

Sterilization: Eto

ID: 4.1 mm

Pressure:

Backflow: 45 psi (Max)

Cracking: 3.5 ± 1 psi

Flow Rate: 203 mL/min @
5 psi



CV-104

NEW

Check Valve, Duckbill

Material: PC, Silicone

Color: Tinted, White

Sterilization: Eto, Gamma



CV-105

NEW

Check Valve, Duckbill,
with Luer Slips

Material: PC, Silicone

Color: Tinted, White

Sterilization: Eto, Gamma



CV-106

NEW

Check Valve, Male Luer
Lock to Luer Slip

Material: PC, Silicone

Color: Tinted, White

Sterilization: Eto, Gamma



CV-107

NEW

Check Valve, Female Luer
Lock to Luer Slip

Material: PC, Silicone

Color: Tinted, White

Sterilization: Eto, Gamma



CV-108

NEW

Check Valve, Duckbill

Material: PC, Silicone
Color: Tinted, White
Sterilization: Eto, Gamma



CV-122

Check Valve, Dual

Material: COPE, PC,
Silicone
Color: Clear
Sterilization: Eto
Pressure: 200 psi



CV-210

Check Valve, Double, Cap,
Rotating Male Luer Lock

Material: LDPE, MABS,
Silicone, THF
Color: Clear, Natural
Sterilization: Eto, E-Beam,
Gamma

ID: 4.1 mm
Pressure: 400 mbar
Flow Rate: 350 mL/min @
1000 mbar



CV-211

Check Valve, Double, with
Cap and Male Luer Lock

Material: LDPE, MABS,
Silicone, THF
Color: Blue, Clear, Natural
Sterilization: Eto, E-Beam,
Gamma

ID: 4.1 mm
Pressure: 400 mbar
Flow Rate: 350 mL/min @
1000 mbar



CV-237

Check Valve, Female to
Male Luer Lock

Material: PC, Silicone, THF
Color: Clear
Sterilization: Eto, E-Beam,
Gamma

Pressure:
Opening: ≤ 6 mbar
Max: 6 bar



CV-249

Check Valve with
Male Luer Lock

Material: MABS, SAN,
Silicone, THF
Color: Blue, Clear
Sterilization: Eto, E-Beam,
Gamma

ID: 4.1 mm
Pressure:
Opening: 200 ± 50 mbar
Max: 21 bar



CV-250

Check Valve with
Male Luer Lock

Material: MABS, SAN,
Silicone, THF
Color: Blue, Clear
Sterilization: Eto

ID: 4.1 mm
Pressure:
Opening: 200 ± 50
mbar
Max: 21 bar



CV-252

Check Valve, Male to
Female Luer Lock

Material: SAN, Silicone,
THF
Color: Blue, Clear
Sterilization: Eto, E-Beam,
Gamma

Pressure:
Opening: 200 ± 50
mbar
Max: 21 bar



CV-264

Check Valve, Inline

Material: MABS, Silicone,
THF
Color: Blue, Clear
Sterilization: Eto, E-Beam,
Gamma

ID: 4.1 mm
Pressure:
Opening: 200 ± 50
mbar
Max: 21 bar



CV-265

Check Valve

Material: PC, Silicone
Color: Blue, Clear
Sterilization: Eto, Gamma,
Steam

ID: 4.1 mm
Pressure:
Opening: ≤ 6 mbar
Max: 6 bar



CV-267

Check Valve with
Female Luer Lock

Material: MABS, SAN,
Silicone, THF
Color: Blue, Clear
Sterilization: Eto, E-Beam,
Gamma

OD: 3.1 mm
Pressure:
Opening: 200 ± 50 mbar
Max: 21 bar



CV-268

Check Valve with
Female Luer Lock

Material: MABS, SAN,
Silicone, THF
Color: Clear, White
Sterilization: Eto, E-Beam,
Gamma

OD: 3.1 mm
Pressure:
Opening: ≤ 6 mbar
Max: 6 bar



CV-272

Check Valve with
Female Luer Lock

Material: MABS, SAN,
Silicone, THF
Color: Blue, Clear
Sterilization: Eto, E-Beam,
Gamma

ID: 4.1 mm
Pressure:
Opening: 200±50 mbar
Max: 21 bar



CV-278

NEW

Check Valve, Female to
Male Luer Lock

Material: SAN, Silicone,
THF
Color: Tinted
Sterilization: Eto, Gamma
Pressure:

Opening: ≤ 6 mbar
Max: 8 bar



CV-308

Check Valve with
Female Luer Lock

Material: PC, Silicone, THF
Color: Clear, Green
Sterilization: Eto, Gamma
OD: 2.5 mm

Pressure:
Opening: 200 ± 50
mbar
Max: 21 bar



CV-337

Check Valve with
Female Luer Lock

Material: PC, Silicone, THF
Color: Clear, Purple
Sterilization: Eto, Gamma
ID: 4.1 mm

Pressure:
Opening: ≤ 6 mbar
Max: 8 bar



CV-274

Check Valve

Material: MABS, Silicone,
THF
Color: Clear, White
Sterilization: Eto, E-Beam,
Gamma

ID: 4.1 mm
Pressure:
Opening: ≤ 6 mbar
Max: 6 bar



CV-279

Check Valve, Male to
Female

Material: SAN, Silicone,
THF
Color: White
Sterilization: Eto, Gamma
Pressure:

Opening: ≤ 6 mbar
Max: 6 bar



CV-308B

Check Valve with
Female Luer Lock

Material: PC, Silicone, THF
Color: Blue, Clear
Sterilization: Eto, E-Beam,
Gamma

OD: 2.5 mm
Pressure:
Opening: 200 ± 50
mbar
Max: 21 bar



CV-358

Check Valve with
Female Luer Lock

Material: MABS, SAN,
Silicone, THF
Color: Clear, White
Sterilization: Eto, Gamma
ID: 2.5 mm

Pressure:
Opening: ≤ 6 mbar
Max: 6 bar



CV-276

Check Valve with
Male Luer Lock

Material: MABS, SAN,
Silicone, THF
Color: Clear
Sterilization: Eto, E-Beam,
Gamma

OD: 3.1 mm
Pressure:
Opening: ≤ 6 mbar
Max: 6 bar



CV-280

Check Valve, Inline

Material: MABS, Silicone,
THF
Color: Clear, White
Sterilization: Eto, E-Beam,
Gamma

ID: 4.1 mm
Pressure:
Opening: ≤ 6 mbar
Max: 6 bar



CV-335

Check Valve with
Female Luer Lock

Material: PC, Silicone, THF
Color: Clear, White
Sterilization: Eto, E-Beam,
Gamma

ID: 3.0 mm
Pressure:
Opening: 200 ± 50
mbar
Max: 21 bar



CV-370

Check Valve with
Female Luer Lock

Material: PC, Silicone, THF
Color: Clear, White
Sterilization: Eto, Gamma
ID: 2.5 mm

Pressure:
Opening: ≤ 6 mbar
Max: 6 bar



CV-397

Check Valve, Large Bore,
High Flow

Material: MABS, Silicone,
THF

Color: Clear

Sterilization: Eto, Gamma

ID: 8.7 mm

Pressure:

Opening: ≤ 20 mbar

Max: 1 bar



CV-401

Check Valve with
Male Luer Lock

Material: MABS, SAN,
Silicone, THF, TPE

Color: Blue, Clear, White

Sterilization: Eto, E-Beam,
Gamma

ID: 4.1 mm

Pressure:

Opening: 200 ± 50 mbar

Max: 21 bar



CV-407

Check Valve with
Male Luer Slip

Material: MABS, Silicone,
THF

Color: Clear, White

Sterilization: Eto, Gamma

ID: 3.2 mm

Pressure:

Opening: ≤ 6 mbar

Max: 8 bar



CV-427

Check Valve with Male
and Female Luer Locks

Material: PETG, Silicone

Color: Clear

Sterilization: Eto, E-Beam,
Gamma

Pressure:

Opening: ≤ 12 mbar

Max: 6 bar



CV-429

Check Valve, T, Rotating
Male Luer Lock, 2 Female
Luers

Material: PC, Silicone, THF

Color: Clear, Purple

Sterilization: Eto, E-Beam,
Gamma

Pressure:

Opening: 200 ± 50

mbar

Max: 21 bar



CV-430

Check Valve, T with 1 Male
and 2 Female Luer Locks

Material: PC, Silicone, THF

Color: Clear

Sterilization: Eto, E-Beam,
Gamma

Pressure:

Opening: ≤ 6 mbar

Max: 6 bar



CV-041

Hemostasis Valve with
Dual Ports

Material: EPDM, PC, POM,
PTFE, Silicone

Color: Clear

Sterilization: Eto

ID (French Precision):
9 Fr



CV-043

Hemostasis Valve, Dual
Seal, Click Activation

Material: EPDM, PC, PTFE,
Silicone

Color: Blue, Clear, Red

Sterilization: Eto

ID (French Precision):
9 Fr



CV-048

Hemostasis Valve

Material: PC, Silicone

Color: Blue, Clear

Sterilization: Eto, Gamma

ID (French Precision):
6 Fr



CV-049

Hemostasis Valve,
One Hand Activation

Material: PC, Silicone,
Stainless Steel

Color: Blue, Clear, Tinted

Sterilization: Eto

ID (French Precision):
9 Fr

Pressure: 900 psi



CV-050

Hemostasis Valve

Material: PC, Silicone

Color: Blue, Clear, Red

Sterilization: N/A

ID (French Precision):
4.5 Fr



CV-083

Hemostasis Valve,
Click Top, Push and
Pull Function

Material: PC, Silicone

Color: Clear

Sterilization: Eto, Gamma

ID (French Precision):
7 Fr



BC-018A

Drain Valve, T

Material: HDPE, PVC
Color: Blue, White
Sterilization: Eto
OD: 9.2 mm



BC-019

Drain Valve

Material: ABS, PVC,
Silicone
Color: Blue, White
Sterilization: Eto
ID: 6.0 mm



BC-020

Drain Valve, Large

Material: ABS, PVC
Color: Blue, White
Sterilization: Eto
ID: 8.4 mm



BC-049

Bag Connector

Material: PVC
Color: White
Sterilization: Eto
ID: 8.4 mm | OD: 9.5 mm



BC-075

Drain Valve, Push / Pull

Material: LDPE, PVC
Color: Blue, Natural
Sterilization: N/A
ID: 6.0 mm



AP-013

TBA Side Port Assembly

Material: ABS, HDPE, PC,
PP, PVC
DEHP Free,
Silicone

Color: Blue, Clear, Natural,
White

Sterilization: Eto, Steam
ID: 0.060"



CV-015

Tuohy Borst Valve, Female

Material: ABS, PC, Silicone
Color: Blue, Clear, Natural
Sterilization: Eto, Steam
ID: 1.143 mm



CV-030

Tuohy Borst Valve, Male

Material: ABS, PC, Silicone
Color: Blue, Clear
Sterilization: Eto
ID: 0.032"



CV-034

Tuohy Borst Valve, Male

Material: ABS, PC, Silicone
Color: Blue, Clear
Sterilization: Eto, Steam
ID (French Precision):
9 Fr



CV-035A

Tuohy Borst Valve, Male

Material: ABS, PC, Silicone
Color: Blue, Clear
Sterilization: Eto, Steam
ID: 1.52 mm
French Catheter: 3-9 Fr



CV-037

Tuohy Borst Valve, Male

Material: ABS, PC, Silicone
Color: Blue, Clear
Sterilization: N/A
ID: 0.125"



CV-039

Tuohy Borst Valve, Male

Material: PVC, Silicone
Color: Blue, Natural
Sterilization: N/A
ID: 0.055"



CV-051

Tuohy Borst Valve, Female

Material: ABS, Nylon
Color: Blue
Sterilization: N/A
ID: 0.06"



CV-063

Tuohy Borst Valve,
Large Bore, Side Port

Material: PC, Silicone
Color: Blue, Clear
Sterilization: Eto, Steam
ID (French Precision):
0-13 Fr



CV-064

Tuohy Borst Valve,
Large Bore, Side Port

Material: PC, Silicone
Color: Natural, Red
Sterilization: Eto, Steam
ID (French Precision):
6-17 Fr



CV-092

Tuohy Borst Valve, Male,
Side Port

Material: ABS, PC, Silicone
Color: Blue, Clear
Sterilization: Eto, Steam
ID: 0.81, 3.3 mm



CV-100A

Tuohy Borst Valve, Male,
Rotating Male Luer

Material: PC, PTFE,
Silicone
Color: Blue, Clear, White
Sterilization: Eto, Steam
ID (French Precision):
9 Fr



CV-102

Tuohy Borst Valve, T, Female to Male

Material: COPE, PC, PTFE, Silicone
Color: Blue, Clear
Sterilization: Eto
ID (French Precision): 9 Fr



CV-103

Tuohy Borst Valve with Side Port Tubing Extension

Material: COPE, PC, PTFE, PVC-DEHP Free, Silicone
Color: Blue, Clear, Natural
Sterilization: Eto
ID (French Precision): 9 Fr



CV-439

Tuohy Borst Valve, Large Bore

Material: PC, Silicone
Color: Clear, Red, Natural
Sterilization: Eto
ID (French Precision): 0-17 Fr



CV-439B

Tuohy Borst Valve, Large Bore

Material: PC, Silicone
Color: Blue, Clear
Sterilization: Eto
ID (French Precision): 0-17 Fr



CV-440

Tuohy Borst Valve, Large Bore

Material: PC, Silicone
Color: Natural, Red
Sterilization: Eto, Steam
ID (French Precision): 6-17 Fr



CV-441

Tuohy Borst Valve Rotor

Material: PC
Color: Red
Sterilization: Eto, Steam
ID: 16.9 mm



CV-442

Tuohy Borst Valve Gasket

Material: Silicone
Color: Natural
Sterilization: Eto, Steam
ID (French Precision): 6-17 Fr



CV-443

Tuohy Borst Valve Nut

Material: PC
Color: Red
Sterilization: Eto, Steam
ID: 6.9 mm



CV-490

Tuohy Borst Valve, Barbed, Large Bore

Material: PC, Silicone
Color: Blue, Natural
Sterilization: Eto, Steam
ID (French Precision): 0-13 Fr



CV-517R

Tuohy Borst Valve, Barbed, Large Bore

Material: PC, Silicone
Color: Natural, Red
Sterilization: Eto, Steam
ID (French Precision): 6-17 Fr



CV-613B

Tuohy Borst Valve, Barbed, Large Bore, Female Luer Lock

Material: PC, Silicone
Color: Blue, Natural
Sterilization: Eto, Steam
ID (French Precision): 0-13 Fr



CV-617R

Tuohy Borst Valve, Barbed, Large Bore, Female Luer Lock

Material: PC, Silicone
Color: Natural, Red
Sterilization: Eto, Steam
ID (French Precision): 6-17 Fr



APHC0750	105	CA-001	11	CA-095	14	CB-027	20	CC-049	24	CP-120	109	CS-033	114	CV-268	118
AP-013	123	CA-001A	11	CA-099	14	CB-030	20	CH-001DHPF	25	CP-301	109	CS-034	114	CV-272	119
BA-009	8	CA-003	11	CA-100	14	CB-031	20	CH-003	25	CR-002	110	CS-035	114	CV-274	119
BA-023	8	CA-005	11	CA-103	14	CB-032	20	CH-005	31	CR-002B	110	CS-036	114	CV-276	119
BB-016	31	CA-008	11	CA-104	14	CB-033	20	CH-010	31	CR-002BP	110	CS-037	114	CV-278	119
BB-113	19	CA-009	18	CA-106	14	CB-034	20	CM-008	31	CR-002BW	110	CS-037B	115	CV-279	119
BC-001	8	CA-012	11	CA-107	14	CB-035	20	CM-009	31	CR-002W	110	CS-038	115	CV-280	119
BC-001R	8	CA-014	11	CA-108	14	CB-036	20	CM-010	31	CR-003	110	CS-039	115	CV-308	119
BC-002	8	CA-016DHPF	18	CA-109	14	CB-037	20	CP-002ABS	105	CR-003W	110	CS-040	115	CV-308B	119
BC-003	8	CA-016ST	18	CA-109B	14	CB-038	20	CP-004A	105	CR-004	110	CS-042	115	CV-335	119
BC-004	8	CA-017	11	CA-109N	14	CB-039	20	CP-004AB	105	CR-005	110	CS-043	115	CV-337	119
BC-005	8	CA-021	18	CA-1090	15	CB-040	20	CP-004AY	105	CR-005W	110	CS-047	115	CV-358	119
BC-006	8	CA-021MP	18	CA-112	15	CB-041N	20	CP-004R	105	CR-005WW	110	CS-047G	115	CV-370	119
BC-007	8	CA-022	11	CA-113	15	CB-041PP	20	CP-006A	105	CR-006	110	CS-047R	115	CV-397	120
BC-008	8	CA-023	11	CA-114	15	CB-042N	21	CP-006AB	105	CR-008	110	CS-047Y	115	CV-401	120
BC-011	8	CA-024	11	CA-115	15	CB-042PP	21	CP-007	105	CR-009	110	CS-048	115	CV-407	120
BC-012A	9	CA-024PP	11	CA-116	15	CB-043	21	CP-008	105	CR-009ASE	110	CV-001	116	CV-427	120
BC-016	9	CA-025	11	CA-117	15	CB-044N	21	CP-008B	105	CR-009W	111	CV-007	116	CV-429	120
BC-018A	122	CA-025BLUE	11	CA-118	15	CB-044PVDF	21	CP-011A	105	CR-013	111	CV-009	116	CV-430	120
BC-019	122	CA-025BLUEA	11	CA-119	15	CB-045N	21	CP-011AB	105	CR-013W	111	CV-013	116	CV-439	124
BC-020	122	CA-025G	11	CA-120	15	CB-046	21	CP-014	105	CR-014	111	CV-015	123	CV-439B	124
BC-021A	29	CA-025NG	11	CA-120V	15	CB-047	21	CP-015	105	CR-015	111	CV-030	123	CV-440	124
BC-022	29	CA-025RED	11	CA-121	15	CB-048	21	CP-017	106	CR-017	111	CV-034	123	CV-441	124
BC-023	29	CA-027	12	CA-122	15	CB-049	21	CP-017B	106	CR-017ASE	111	CV-035A	123	CV-442	124
BC-024	29	CA-030	12	CA-123	15	CB-050	21	CP-017R	106	CR-017W	111	CV-037	123	CV-443	124
BC-024A	29	CA-030RED	12	CA-124	15	CB-051	21	CP-020	106	CR-019BB	111	CV-039	123	CV-490	124
BC-025	29	CA-031	12	CA-125	15	CB-052	21	CP-024	106	CR-019BBGW	111	CV-041	121	CV-517R	124
BC-026	29	CA-032	12	CA-127	18	CB-053B	21	CP-025	106	CR-019BBWW	111	CV-043	121	CV-613B	124
BC-028	47	CA-032B	12	CA-128	18	CB-054	21	CP-025P	106	CR-019BLUE	111	CV-045	116	CV-617R	124
BC-029	9	CA-032G	12	CA-129	15	CB-055	21	CP-025W	106	CR-019G	111	CV-048	121	CY-001DEHT	33
BC-033	19	CA-032PP	12	CA-131	15	CB-056	21	CP-036BL	106	CR-019GW	111	CV-049	121	CY-002DEHT	33
BC-034	19	CA-032R	12	CA-131V	16	CB-057	21	CP-036BR	106	CR-019WHEEL	111	CV-050	121	CY-002DHPFG	33
BC-036	47	CA-032RG	12	CA-132	16	CB-058	22	CP-036GR	106	CR-020WHITE	112	CV-051	123	CY-002NGDEHT	33
BC-040DHPF	19	CA-033	12	CA-133	16	CB-059	22	CP-036YE	106	CR-021	112	CV-055	116	CY-003DEHT	33
BC-043DEHT	19	CA-039	12	CA-133B	16	CB-060	22	CP-038	106	CR-022	112	CV-060	116	CY-003NGDEHT	33
BC-048DEHT	25	CA-041	12	CA-134	16	CB-061	22	CP-039	106	CR-022W	112	CV-063	123	CY-004DEHT	33
BC-049	122	CA-042	12	CA-135	16	CB-062	22	CP-040	106	CR-022WB	112	CV-064	123	CY-005A	33
BC-052	25	CA-043	12	CA-136	16	CC-010	23	CP-042	107	CR-024	112	CV-070	116	CY-005DHPFG	33
BC-053	25	CA-045	18	CA-148	16	CC-011	23	CP-043	107	CR-117	112	CV-071	116	CY-006A	33
BC-054DHPF	9	CA-046	18	CA-151	16	CC-012	23	CP-044	107	CR-118	112	CV-074	116	CY-006C	33
BC-055	9	CA-048	12	CA-151B	16	CC-013	23	CP-045	107	CS-001	113	CV-075	116	CY-006DEHT	33
BC-056DEHT	9	CA-049W	12	CA-152	16	CC-014	23	CP-045GR	107	CS-002	113	CV-081	116	CY-007DEHT	33
BC-057	9	CA-050	12	CA-153	16	CC-016	23	CP-045YE	107	CS-002N	113	CV-083	121	CY-008DEHT	33
BC-059	9	CA-052	13	CA-155	45	CC-017	23	CP-046	107	CS-003	113	CV-084	117	CY-009DEHT	33
BC-063	9	CA-052A	13	CA-157	65	CC-018	23	CP-047	107	CS-003B	113	CV-086	117	CY-009NGDEHT	34
BC-064	9	CA-053	13	CA-158	16	CC-019	23	CP-049	107	CS-003R	113	CV-088	117	CY-010DEHT	34
BC-065	19	CA-054	13	CA-160	16	CC-020	23	CP-050	107	CS-004	113	CV-092	123	CY-010DHPF	34
BC-066	19	CA-055	13	CA-161	16	CC-021	23	CP-051	107	CS-005	113	CV-094	117	CY-011DHPF	34
BC-067	9	CA-055ABS	13	CA-162	13	CC-022	23	CP-052	107	CS-006	113	CV-095	117	CY-012	34
BC-068	9	CA-055B	13	CA-164	16	CC-023	23	CP-053	107	CS-007	113	CV-096	117	CY-016	34
BC-068E	10	CA-060	13	CA-165	16	CC-025	23	CP-054	107	CS-009B	113	CV-097	117	CY-017	34
BC-069	19	CA-061	13	CA-167	17	CC-027	23	CP-055	107	CS-009DG	113	CV-100A	123	CY-018	34
BC-070	19	CA-062	13	CA-168	17	CC-028	23	CP-056	108	CS-009LG	113	CV-101	117	CY-019	34
BC-071	19	CA-063	18	CA-169	18	CC-029	23	CP-057	108	CS-009P	113	CV-102	124	CY-019W	34
BC-075	122	CA-065S	18	CA-170	18	CC-030	23	CP-058	108	CS-009R	113	CV-103	124	CY-020DEHT	34
BC-078	10	CA-065SW	18	CA-171	17	CC-031	24	CP-059	108	CS-009W	113	CV-104	117	CY-020DHPF	34
BC-083	29	CA-066	13	CA-177	17	CC-032	24	CP-061	108	CS-009Y	113	CV-105	117	CY-020DHPFG	34
BC-084	29	CA-068	13	CA-200	17	CC-033	24	CP-063	108	CS-013	113	CV-106	117	CY-020GSDEHT	34
BC-085	29	CA-068NV	13	CA-324PP	17	CC-034	24	CP-064	108	CS-014	114	CV-107	117	CY-021DEHT	34
BC-091	10	CA-070	13	CB-001A	19	CC-035	24	CP-067	108	CS-014R	114	CV-108	118	CY-022DHPF	35
BC-111	10	CA-080	13	CB-002	19	CC-036	24	CP-068	108	CS-014W	114	CV-122	118	CY-022DHPFGS	35
BC-111TPE	10	CA-080W	13	CB-003	19	CC-038	24	CP-070	108	CS-016	114	CV-157	65	CY-023DEHPFREE	35
BC-112	10	CA-086A	13	CB-006	19	CC-039	24	CP-071	108	CS-016B	114	CV-210	118	CY-026	35
BC-113	10	CA-087	13	CB-008	19	CC-040	24	CP-088B	108	CS-016G	114	CV-211	118	CY-027DHPF	35
BC-119	10	CA-087W	14	CB-012	19	CC-041	24	CP-088G	108	CS-016P	114	CV-237	118	CY-028	35
BC-131	25	CA-089	14	CB-015	19	CC-042	24	CP-088PPW	108	CS-016R	114	CV-249	118	CY-030	35
BH-011	74	CA-091	14	CB-016	19	CC-043	24	CP-088R	108	CS-016Y	114	CV-250	118	CY-031	35
BH-013	74	CA-093	14	CB-022	20	CC-044	24	CP-088W	109	CS-018	114	CV-252	118	CY-031ABS	35
BH-015	74	CA-094	14	CB-023	20	CC-045	24	CP-106B	109	CS-018B	114	CV-264	118	CY-033	35
BH-016	74	CA-094GS	14	CB-024	20	CC-047	24	CP-106R	109	CS-031	114	CV-265	118	CY-033ABS	35
BH-017	74	CA-094V	14	CB-025	20	CC-048	24	CP-113	109	CS-032	114	CV-267	118	CY-035	35

CY-036	35	DC-127	45	DY-018	26	FC-012WHT	28	FF-226	56	IS-058	61	LB-011	89	LF-027A	88
CY-038	35	DC-128DFG	45	DY-019	26	FC-014GRE	28	FF-301	56	IS-060	61	LB-012C	89	LF-027BDEHT	88
CY-039	35	DC-128DFNG	46	DY-021	26	FC-016ORA	28	FF-303	56	IS-062	61	LB-013A	89	LF-027DEHT	88
CY-040	36	DC-130	44	DY-022	26	FC-018RED	28	FF-304	56	IS-063	61	LB-014	83	LF-027DF	88
CY-041	36	DC-131	39	DY-025A	26	FC-020YEL	28	FH-001	74	IS-065	61	LB-015	83	LF-027DHPF	88
CY-041-3	36	DC-133	44	DY-027	26	FC-022PIN	28	FH-002	74	IS-066	65	LB-018A	89	LF-027RDEHT	88
CY-041PVC	36	DC-134	39	DY-029	26	FC-025	58	FH-009	74	IS-067	63	LB-022	22	LF-028	85
CY-042	36	DC-137	44	DY-030	26	FF-003	51	FR-001	58	IS-068	63	LB-028	89	LF-029A	86
CY-043	36	DC-138	44	DY-031	26	FF-003A	51	FR-002	58	IS-069	61	LB-029	83	LF-029DF	86
CY-044	36	DC-140	39	DY-033DEHT	26	FF-004	51	FR-009	58	IS-071	63	LB-030	22	LF-031	88
CY-045	36	DC-141	46	DY-035DHPF	26	FF-005	51	FR-014	58	IS-079	63	LB-034	89	LF-032	93
CY-046	36	DC-142	40	DY-036DHPF	27	FF-007	51	FR-016	58	IS-080	63	LB-036	83	LF-034	86
CY-047	36	DC-143	40	DY-039DEHT	27	FF-008	51	FR-017	58	IS-081	63	LB-037	83	LF-039	86
CY-057	36	DC-144	40	DY-041	27	FF-010	51	FR-018	58	IS-082	63	LB-038	83	LF-041	87
CY-058	36	DC-145	40	DY-045DHPF	51	FF-030	51	FR-019	58	IS-083	63	LB-039	84	LF-041PVC	87
CY-059	36	DC-146	46	DY-045DEHT	51	FF-034	51	FR-021	58	IS-092	61	LB-040	89	LF-042	87
CY-060	36	DC-147	40	DY-047	51	FF-034C	51	FR-025	58	IS-094	63	LB-041	89	LF-048	86
CY-061	36	DC-152	40	DY-047A	51	FF-037	52	FR-030	58	IS-095	64	LB-042	89	LF-059	85
CY-064	37	DC-154	40	DY-047DHPFG	51	FF-038	52	GS-011	74	IS-096	64	LB-108	84	LF-064	88
CY-065	37	DC-156	40	DY-048	104	FF-051	52	GW-010	74	IS-097	64	LB-117A	83	LF-066	86
CY-077	37	DC-158	40	DY-049	104	FF-052	52	GW-011	74	IS-100	61	LB-118A	83	LF-066B	86
CY-079	37	DC-159	44	DY-050DHPF	104	FF-053	52	GW-012	74	IS-101	61	LB-120A	84	LF-066R	86
CY-092	37	DC-160	40	DY-051	104	FF-054A	52	GW-013	75	IS-104	65	LB-121A	84	LF-071	85
CY-092M	37	DC-161	40	DY-052	51	FF-055	52	GW-014	75	IS-105	65	LB-122A	84	LF-074	86
CY-092PVC	37	DC-161B	40	DY-053	27	FF-057	52	GW-014B	75	IS-106	65	LB-130A	84	LF-079	85
CY-093PVC	37	DC-162	41	DY-054	27	FF-059	52	GW-014C	75	IS-107	65	LB-131A	83	LF-079CL	85
CY-101	37	DC-163	44	DY-055	27	FF-060	52	GW-015	75	IS-108	61	LB-132A	83	LF-087	90
CY-161	37	DC-165	44	DY-056	27	FF-061	52	GW-017	75	IS-109	61	LB-133A	84	LF-088	85
DC-004	43	DC-166	44	DY-057	27	FF-064	52	GW-018	75	IS-110	62	LB-135A	89	LF-091GS	85
DC-006	38	DC-167	41	DY-058	27	FF-066	53	IS-001	59	IS-112	62	LB-135APP	89	LF-092	89
DC-007DHPF	38	DC-168	46	DY-059	27	FF-077	53	IS-004	59	IS-112S	62	LB-142A	89	LF-093	86
DC-008ADF	45	DC-169	41	DY-060DEHT	27	FF-077B	53	IS-005	63	IS-113	62	LB-143A	89	LF-094	85
DC-008DF	45	DC-171	41	DY-061	27	FF-083	53	IS-006	59	IS-118	65	LB-145	83	LF-094GS	85
DC-008NG	45	DC-172	41	DY-062DEHT	51	FF-084	53	IS-007A	63	IS-119	65	LB-146	84	LF-095	87
DC-012	38	DC-174	41	DY-077M	104	FF-085	53	IS-008A	63	IS-120	65	LB-147	83	LF-097	86
DC-013F	38	DC-175	41	DY-110	51	FF-086	53	IS-009A	63	IS-121	66	LB-148	83	LF-102	87
DC-013K	68	DC-176	46	DY-110NG	51	FF-087	53	IS-010	59	IS-122	66	LB-149	89	LF-103	88
DC-015	43	DC-177DEHT	41	DY-143	27	FF-088	53	IS-010C	59	IS-124	66	LB-150	83	LF-106	86
DC-018	38	DC-178	46	EN-007	47	FF-089	53	IS-011	59	IS-125	66	LB-227	83	LF-107	87
DC-020DEHT	43	DC-179	41	EN-0080	47	FF-093	53	IS-011CL	59	IS-126	66	LC-011	92	LF-112	85
DC-021DEHT	43	DC-180	41	EN-008P	47	FF-094	53	IS-014	59	IS-127	66	LC-012	92	LF-113	86
DC-024	38	DC-181	44	EN-0090	47	FF-095	54	IS-014S	59	IS-128	62	LC-013	92	LF-114	86
DC-025	38	DC-186	41	EN-0100R	47	FF-096	54	IS-015	63	IS-129	66	LC-014	92	LF-115	87
DC-026	38	DC-187	41	EN-010P	47	FF-097	54	IS-016	63	IS-130	62	LC-015	92	LF-116	87
DC-031A	38	DC-188	42	EN-010W	47	FF-098	54	IS-018	59	IS-131	62	LC-016	92	LF-117	88
DC-031NG	38	DC-189	46	EN-0110R	47	FF-099	54	IS-019	59	IS-132	66	LC-018	92	LF-118	85
DC-032	43	DC-190	46	EN-011P	47	FF-101	54	IS-022	59	IS-133	64	LF-001	87	LF-119	90
DC-033	43	DC-192	42	EN-0150R	47	FF-102	54	IS-025	59	IS-134	66	LF-002	87	LF-119GS	90
DC-034	38	DC-193	42	EN-015P	48	FF-103	54	IS-025A	59	IS-136	66	LF-002PVC	87	LF-120B	90
DC-039NDGS	38	DC-195	44	EN-0160R	48	FF-104	54	IS-026	65	IS-137	66	LF-002PVCNG	87	LF-122	87
DC-040DEHT	39	DC-196	42	EN-016P	48	FF-105	54	IS-028	65	IS-138	64	LF-003	86	LF-123	90
DC-046DFG	39	DC-197-60	44	EN-018	48	FF-106	54	IS-030	59	IS-139	64	LF-004	88	LF-125	88
DC-051DEHT	39	DC-198	46	EN-018B	48	FF-107	54	IS-032	60	IS-141	64	LF-005	86	LF-135	85
DC-054	43	DC-228G	42	EN-019P	48	FF-108	55	IS-033	60	IS-143	64	LF-007	85	LF-155	88
DC-054G	43	DY-001	25	EN-019PA	48	FF-110	55	IS-036DEHT	60	IS-146	64	LF-008	88	LF-178	85
DC-062	45	DY-002	25	EN-020P	48	FF-111	55	IS-039	60	IS-147	64	LF-010	88	LL-002	90
DC-063	45	DY-002R	25	EN-020PA	48	FF-112	55	IS-039G	60	IS-148	64	LF-011	88	LL-002PC	90
DC-064	39	DY-004ADEHT	25	EN-020PT	48	FF-116	55	IS-041DEHT	60	IS-149	62	LF-012	85	LL-003	90
DC-068DEHT	45	DY-006	25	EN-0210	48	FF-120	55	IS-045	60	IS-150B	66	LF-013	85	LL-006	90
DC-069	43	DY-007BDEHT	25	EN-021P	48	FF-122	55	IS-046	60	IS-150G	67	LF-014	90	LL-008	90
DC-076	45	DY-007BG	25	EN-023	49	FF-123	55	IS-047	60	IS-1500R	67	LF-016	87	LL-009	90
DC-079DEHPFREE	39	DY-007DEHT	25	EN-024	49	FF-124	55	IS-048	60	IS-150R	67	LF-016DHPFG	87	LL-014	90
DC-087DEHT	39	DY-007G	26	EN-026	49	FF-125	55	IS-049	60	IS-150W	67	LF-018	85	LL-015	90
DC-095	43	DY-007RDEHT	26	EN-029	49	FF-127	55	IS-052	60	IS-151	62	LF-021DEHT	85	LL-018	90
DC-097GDF20	43	DY-007RG	26	EN-033	49	FF-128	55	IS-052A	60	IS-154C	62	LF-022DEHT	85	LL-019	91
DC-100	45	DY-008DEHT	26	EN-035	49	FF-129	56	IS-054	60	IS-211	62	LF-023DEHT	86	LL-026	91
DC-115	39	DY-011DEHT	51	EN-043	49	FF-130	56	IS-054B	60	IS-245	64	LF-023DHPF	86	LM-002	91
DC-120	39	DY-012	51	FC-006LGR	28	FF-131	56	IS-054R	61	LB-006	84	LF-024	87	LM-003	91
DC-122	39	DY-013ADHPF	51	FC-008LBL	28	FF-141	56	IS-055	61	LB-006PC	84	LF-024DHPF	87	LM-004	91
DC-123	45	DY-014DEHT	51	FC-010	58	FF-178	56	IS-056	61	LB-008	22	LF-026DEHT	88	LM-005	91
DC-126	43	DY-015	51	FC-010BLK	28	FF-224	56	IS-057	61	LB-010	89	LF-026RNG	88	LM-010	79

LM-011	93	LM-110	80	NA-038	96	SP-057	70	SYG-005LSM-LOR	7	TC-143	104
LM-014	79	LM-111	79	NH-014	94	SP-058	72	SYG-005NR-LOR	7	TT-002A	102
LM-014-CAP-A	79	LM-114	78	NH-015	94	SP-059	72	SYG-010LL	6	TT-006A	102
LM-014-CAP-B	79	LM-114GR	78	NH-016	94	SP-060	70	SYG-010LL-LOR	7	TT-007A	102
LM-015	78	LM-114GY	78	NH-017	94	SP-061	70	SYG-010LS	6	TT-008A	102
LM-016	79	LM-115	79	NH-018	94	SP-079	70	SYG-010LS-LOR	7	TT-010	32
LM-017	76	LM-116	91	NH-019	94	SP-087	70	SYG-010LSM	6	TT-038	32
LM-022	76	LM-117	77	NH-020	94	SP-128	73	SYG-010LSM-LOR	7	TT-041	102
LM-024	82	LM-118	79	NH-021	94	SP-140B	57	SYG-020LL	6	TT-042	102
LM-025	81	LM-119	78	NH-022	94	SP-230	70	SYG-020LS	6	TT-046	102
LM-026	79	LM-120	77	PO-004A	31	SR-001	100	SYG-020LSM	6	TT-048	102
LM-027	78	LM-121	77	PO-008A	31	SR-004	100	SYG-030LL	6	TT-059	102
LM-027A	78	LM-122	78	RF-003	56	SR-005	100	SYG-030LL-LOR	7	TT-061DEHT	102
LM-028	77	LM-123	81	RW-009	56	SR-024	100	SYG-030LS	6	TT-062DEHT	102
LM-028A	77	LM-124	78	RW-010	57	SR-025	100	SYG-030LS-LOR	7	TT-063DEHT	102
LM-029	76	LM-125	80	RW-011	71	SR-026	100	SYG-030LSM	6	TT-064	22
LM-029A	76	LM-126	82	RW-012	71	SR-027	100	SYG-030LSM-LOR	7	TT-070	102
LM-030	82	LM-128	82	RW-014	71	ST-006	97	SYG-050LL	6	TT-072	102
LM-030B	82	LM-129	80	RW-017	57	ST-006LL	97	SYG-050LL-LOR	7	TT-073	102
LM-030R	82	LM-129NC	80	RW-018	57	ST-011NC	97	SYG-050LS	6	TT-074	102
LM-031ABS	77	LM-129PVC	81	RW-019	57	ST-033	97	SYG-050LSM	6	TT-075	102
LM-031CAPDEHT	77	LM-131	78	RW-020	57	ST-034	97	SYG-100LL	6	TT-076	102
LM-031DEHT	77	LM-134	79	RW-021	57	ST-035	97	SYG-100LS	6	TT-077	102
LM-032	79	LM-137	79	RW-025	57	ST-036	97	SYG-100LSM	6	TT-080	102
LM-032A	79	LM-138	81	RW-026	68	ST-037	97	SYG-150LL	6	TT-081	103
LM-033	78	LM-139	81	RW-027	68	ST-037NC	97	SYG-150LS	6	TT-088	103
LM-034	79	LM-140	81	RW-122	57	ST-039	97	SYG-150LSM	6	TT-089	103
LM-034B	80	LM-140B	81	SC-009	75	ST-059	97	SYG-200LL	6	TT-099	103
LM-034BODY	80	LM-141	76	SC-011	75	ST-064	97	SYG-200LS	6	TT-100	103
LM-034-11	80	LM-141B	76	SC-015	75	ST-066	98	SYG-200LSM	6	TT-101	103
LM-035	77	LM-142	81	SP-002	68	ST-067	98	SYG-250LL	6	TT-102	103
LM-036	82	LM-144	81	SP-003	71	ST-070	98	SYG-250LS	6	TT-103	103
LM-037	81	LM-144B	81	SP-006	71	ST-071	98	SYG-250LSM	6	VA-01313	73
LM-037B	81	LM-145	76	SP-009	68	ST-072	98	TA-008	32	VA-01320	73
LM-038	82	LM-145B	76	SP-011	68	ST-077B	98	TA-009	32	VA-013F	73
LM-039	81	LM-146	76	SP-012	68	ST-077R	98	TA-011	32	VA-013IS	73
LM-041	79	LM-147	78	SP-013	68	ST-077Y	98	TA-011ABS	32	VA-014	73
LM-042HUB	92	LM-148B	76	SP-014	68	ST-078	98	TA-012	32	VA-015	73
LM-050	81	LM-150	81	SP-014CC	68	ST-079	98	TA-013	32	VA-02020	73
LM-052	76	LM-151	77	SP-014CLEAR	68	ST-081	98	TA-014	32	VA-020F	73
LM-054	76	LM-152	77	SP-015	68	ST-082	98	TA-014PVC	32	VA-020IS	73
LM-055	76	LM-154	81	SP-015A	69	ST-083	99	TA-015	32	VAC160P	101
LM-059	80	LM-174	91	SP-015WC	69	ST-087	99	TA-016	32		
LM-061	91	LS-001A	93	SP-016	71	ST-088	99	TA-018	32		
LM-063	81	LS-001AC	93	SP-017	71	ST-154	99	TA-019	32		
LM-064	78	LS-002	93	SP-017A	71	ST-162	99	TA-020	32		
LM-066	82	LS-003	93	SP-020	69	SY-026	101	TA-021	32		
LM-067	82	LS-004	93	SP-021	69	SY-027	101	TA-024	32		
LM-068	80	LS-005DEHT	93	SP-023	69	SY-032	101	TC-005	30		
LM-074	91	LS-007	93	SP-024	71	SY-034	101	TC-006	104		
LM-074R	91	NA-001	95	SP-024C	71	SY-035	101	TC-008DHPF	30		
LM-081	91	NA-002	95	SP-025	69	SY-037	101	TC-009	31		
LM-082	91	NA-003	95	SP-029	69	SY-045	101	TC-015	32		
LM-082A	91	NA-004	95	SP-031	71	SY-046	101	TC-016	30		
LM-082B	91	NA-005	95	SP-032	71	SY-047	101	TC-030A	30		
LM-084	84	NA-010	95	SP-033	69	SY-063	101	TC-031A	30		
LM-088GS	76	NA-011	95	SP-033C	69	SY-960	101	TC-032A	30		
LM-089	77	NA-020	95	SP-034	69	SYG-001LL	6	TC-034B	30		
LM-089GS	77	NA-021	95	SP-035	69	SYG-001LS	6	TC-035A	30		
LM-090	78	NA-023	95	SP-038	72	SYG-001LSM	6	TC-037A	30		
LM-091	80	NA-025	95	SP-039	73	SYG-002LL	6	TC-037G	30		
LM-093	79	NA-026	95	SP-040	69	SYG-002LS	6	TC-039A	30		
LM-097	82	NA-027	96	SP-042	44	SYG-002LSM	6	TC-040A	30		
LM-098	80	NA-028	96	SP-044	72	SYG-003LL	6	TC-043	22		
LM-099	80	NA-029	96	SP-045	72	SYG-003LS	6	TC-044	31		
LM-099ABS	80	NA-030	96	SP-047	70	SYG-003LSM	6	TC-045	31		
LM-102	76	NA-031	96	SP-048	72	SYG-005LL	6	TC-050	30		
LM-104	77	NA-032	96	SP-049	72	SYG-005LL-LOR	7	TC-053DHPF	30		
LM-105	80	NA-034	96	SP-050	70	SYG-005LS	6	TC-057	30		
LM-106	78	NA-036	96	SP-052	70	SYG-005LS-LOR	7	TC-057R	30		
LM-108	77	NA-037	96	SP-056	72	SYG-005LSM	6	TC-117	104		

ABS	Acrylonitrile Butadiene Styrene	PES	Polyether Sulfone
BIIR	Bromobutyl Rubber	PET	Polyethylene Terephthalate
COPE	Copolyester	PETG	Polyethylene Terephthalate Glycol
DEHT	Dioctyl Terephthalate	PIP	Polyisoprene
ETO	Ethylene Oxide Sterilized	POM	Polyoxymethylene (Acetal)
EPDM	Ethylene Propylene Diene Monomer	PP	Polypropylene
ePTFE	Expanded Polytetrafluoroethylene	PS	Polystyrene
EVA	Ethylene Vinyl Acetate	PSU	Polysulfone
HDPE	High Density Polyethylene	PTFE	Polytetrafluoroethylene
HIPS	High Impact Polystyrene	PUR	Polyurethane
ID	Internal Diameter	PVA	Polyvinyl Acetate
IR	Isoprene Rubber	PVC	Polyvinyl Chloride
LDPE	Low Density Polyethylene	PVDF	Polyvinylidene Fluoride
MABS	Methyl Methacrylate ABS	SAN	Styrene Acrylonitrile
OD	Outer Diameter	SBC	Styrene-Butadiene Copolymer
OP	Operating Pressure	SEBS	Styrene Ethyl Butylene Styrene
PA	Polyamide	SS	Stainless Steel
PBT	Polybutylene Terephthalate	THF	Tetrahydrofuran
PC	Polycarbonate	TPE	Thermoplastic Elastomer
PE	Polyethylene		

FRACTION-DECIMAL CONVERSION CHART

INCHES		MILLIMETERS	INCHES		MILLIMETERS
1/64	.015625	.3969	33/64	.515625	13.096
1/32	.03125	.7938	17/32	.53125	13.493
3/64	.046875	1.1906	35/64	.546875	13.890
1/16	.0625	1.5875	9/16	.5625	14.287
5/64	.078125	1.9844	37/64	.578125	14.684
3/32	.09375	2.3813	19/32	.59375	15.081
7/64	.109375	2.7781	39/64	.609375	15.478
1/8	.125	3.1750	5/8	.625	15.875
9/64	.140625	3.5719	41/64	.640625	16.271
5/32	.15625	3.9688	21/32	.65625	16.668
11/64	.171875	4.3656	43/64	.671875	17.065
3/16	.1875	4.7625	11/16	.6875	17.462
13/64	.203125	5.1594	45/64	.703125	17.859
7/32	.21875	5.5563	23/32	.71875	18.256
15/64	.234375	5.9531	47/64	.734375	18.653
1/4	.250	6.3500	3/4	.750	19.050
17/64	.265625	6.7469	49/64	.765625	19.447
9/32	.28125	7.1438	25/32	.78125	19.843
19/64	.296875	7.5406	51/64	.796875	20.240
5/16	.3125	7.9375	13/16	.8125	20.6375
21/64	.328125	8.3344	53/64	.828125	21.0345
11/32	.34375	8.7313	27/32	.84375	21.431
23/64	.359375	9.1282	55/64	.859375	21.8282
3/8	.375	9.5250	7/8	.875	22.2251
25/64	.390625	9.9219	57/64	.890625	22.6220
13/32	.40625	10.3188	29/32	.90625	23.0188
27/64	.421875	10.7157	59/64	.921875	23.4157
7/16	.4375	11.1125	15/16	.9375	23.8126
29/64	.453125	11.5094	61/64	.953125	24.2095
15/32	.46875	11.9063	31/32	.96875	24.6063
31/64	.484375	12.3032	63/64	.984375	25.0032
1/2	.500	12.7001	1	1.000	25.4001

LENGTH		VOLUME	
1 millimeter (mm)	0.0393701 inches 0.00328084 feet 3 Fr (French Catheter)	1 milliliter (mL)	0.033814 fl oz 0.00105669 quarts 0.000264172 gallons 20 drops
1 centimeter (cm)	0.393701 inches 0.0328084 feet	1 liter (L)	33.814 fl oz 1.05669 quarts 0.264172 gallons 20000 drops
1 meter (m)	39.3701 inches 3.28084 feet	1 fluid ounce (fl oz)	29.5735 mL 0.0295735 L 591.471 drops
1 inch (in)	25.4 mm 2.54 cm	1 quart (qt)	946.353 mL 0.946353 L 18927.1 drops
1 foot (ft)	304.8 mm 30.48 cm 0.3048 m	1 gallon (gal)	3785.41 mL 3.78541 L 75708.2 drops
1 yard (yd)	914.4 mm 91.44 cm 0.9144 m		
1 French Catheter (Fr)	1/3 mm		
WEIGHT		TEMPERATURE	
1 gram (g)	0.035274 oz 0.00220462 lbs	0° Celsius (°C)	32 °F 273.15 K $^{\circ}\text{C} = (^{\circ}\text{F} - 32) / (9/5)$ $^{\circ}\text{C} = \text{K} - 273.15$
1 kilogram (kg)	35.274 oz 2.20462 lbs	0° Fahrenheit (°F)	-17.7778 °C 255.372 K $^{\circ}\text{F} = ^{\circ}\text{C} \times 9/5 + 32$ $^{\circ}\text{F} = \text{K} \times 9/5 - 459.67$
1 ounce (oz)	28.3495 g 0.0283495 kg		
1 pound (lb)	453.592 g 0.453592 kg		
PRESSURE		0° Kelvin (K)	-273.15 °C -459.67 °F $\text{K} = ^{\circ}\text{C} + 273.15$ $\text{K} = (^{\circ}\text{F} + 459.67) \times 5/9$
1 millibar (mbar)	0.0145038 psi		
1 bar	14.5038 psi		
1 pound per sq. in. (psi)	68.9476 mbar 0.0689476 bar		

